



sartorius stedim
biotech

Vivapure®
miniprep and maxiprep
Protein Purification Spin Columns



turning science into solutions



Antibody and His-Tagged Protein Purification Kits



Rapid Purification with High Yields

Vivapure® miniprep and maxiprep Kits are spin column based kits for fast and effective purification of His-tagged proteins and antibodies.

Spin columns have the advantage of speed over gravity drip columns and batch protocols.

With the patented FlowGo regulator the sample residence time is extended to assure adequate sample binding to the resin. Due to this, Sartorius Vivapure® miniprep and maxiprep spin column kits combine the merits of spin columns and gravity drip columns resulting in rapid purification with up to 95% protein recovery and purity.



All spin columns can conveniently be used in a centrifuge. For processing larger sample volumes, e. g. from diluted cell culture supernatants, the Vivapure® maxiprep spin columns can additionally be run (see picture top of page 4) with a peristaltic pump collar (VS-PPCSC).

The Vivapure® miniprep and maxiprep columns come in a variety of different configurations for your convenience. They can be purchased as completely ready-to-use kits with buffers and ultrafiltration devices but just as well as stand alone spin columns in small packs or large packs for frequent users.

Vivapure miniprep and maxiprep Protein A & G Spin Columns



Protein A & G for Antibody Purification

Protein A & G miniprep	Centrifuge
Sample size	0.65 ml
Typical Binding Capacity	1 mg IgG/column
Number of re-uses	3



Protein A & G maxiprep	Centrifuge ¹
Sample size	20 ml
Typical Binding Capacity	20 mg IgG/column
Number of re-uses	5

¹ Use the peristaltic pump accessory (VS-PPCSC) for larger volumes


Applications

- Use Vivapure Protein A or G miniprep or maxiprep spin columns to purify antibodies from serum, ascites and tissue culture supernatant derived from static cultures or bioreactors.
- Purification of a wide range of monoclonal and polyclonal antibodies
- Rapid clean up of immuno-histochemical, immunofluorescence experiments of Western blots
- Removal of abundant IgG from human plasma to simplify samples
- Titration of antibody yield during cell culture


Set up Vivapure[®] spin columns

The plug and play concept for the Vivapure[®] miniprep spin columns


Loading the plug into the spin column



Place the plug into the spin column with the recessed end uppermost.



Push the plug **fully** into the tapered end of the spin column using the plug insertion tool.



Spin
Re-use

It is now ready for pre-equilibration with binding buffer followed by centrifugation.

Vivapure miniprep and maxiprep MC Spin Columns



Coming soon!

IMAC for His-tagged Protein Purification:

Protein MC miniprep Kits	Centrifuge
Sample size	0.65 ml
Typical Binding Capacity	1 mg His-tagged protein
Number of re-uses	2

Protein MC maxiprep Kits	Centrifuge ¹
Sample size	20 ml
Typical Binding Capacity	10 mg His-tagged protein
Number of re-uses	2

¹ Use the peristaltic pump accessory (VS-PPCSC) for larger volumes


Applications

Ready-to-use, robust and reproducible kits for purifying His-tagged proteins from bacteria, insect, mammalian and yeast cells under native or denaturing conditions in Molecular Biology, Biochemistry or Structural Biology laboratories.


- Easy screening for soluble expression of His-tagged proteins
- Run Vivapure[®] maxiprep MC columns in centrifugal mode for volumes < 20 ml or with peristaltic pump for > 20 ml
- Purification of recombinant proteins for use as substrates for enzyme assays, structural studies or for raising antibodies
- Titration of His-tagged protein yield during cell culture
- Remove free His-tag after cleavage from purified recombinant protein
- Multiple expression profiling to generate critical cell culture performance data and for predicting optimal harvest times

The plug and play concept for the Vivapure[®] maxiprep spin columns



Loading the plug into the spin column




Place the plug into the spin column with the recessed end uppermost.



Push the plug **fully** into the tapered end of the spin column using the plug insertion tool.



It is now ready for pre-equilibration with binding buffer followed by centrifugation.



After use, the plug is removed using the plug insertion tool.

Working with samples > 20 ml



Accessory

For working with sample volumes larger than 20 ml, e.g. diluted cell culture supernatants, a sealing cap and peristaltic pump collar (VS-PPCSC) for Vivapure® maxiprep columns offer a fast and easy to use alternative to multiple centrifugation steps. The sample is pumped into the maxiprep spin column through a tube attached to the sealing cap with a connector. To ensure the usage of high sample loading flow rates, the peristaltic pump collar securely locks the sealing cap to the column barrel.

Required Equipment

Variable speed peristaltic pump capable of speeds less than 20 rpm. E. g. Masterflex pump (VFP001, 240V) VFP002 115V, Masterflex easy load pump head-size 16, VF A012.

Ordering information

Product name	Pack size	Product no.
Vivapure A Starter Pack	2	VS-ARSTPKA2
Vivapure miniprepA Kit	16 + UF concentrators + buffer	VS-ARAMINIK
Vivapure miniprepA Bulk Pack	48	VS-ARAMINIB
Vivapure maxiprepA Kit	4 + UF concentrators + buffer	VS-ARAMAXIK
Vivapure maxiprepA Bulk Pack	12	VS-ARAMAXIB
Vivapure A Buffer Pack		VS-ARABUFPK
Vivapure G Starter Pack	2	VS-ARSTPKG2
Vivapure miniprepG Kit	16 + UF concentrators + buffer	VS-ARGMINIK
Vivapure miniprepG Bulk Pack	48	VS-ARGMINIB
Vivapure maxiprepG Kit	4 + UF concentrators + buffer	VS-ARGMAXIK
Vivapure maxiprepG Bulk Pack	12	VS-ARGMAXIB
Vivapure G Buffer Pack		VS-ARGBUFPK
Sealing Cap & Peristaltic Pump Collor 1		VS-PPCSC

Binding Affinities of Protein A and Protein G

Antibody	Protein A	Protein G	Antibody	Protein A	Protein G
Human IgG1	****	****	Hamster IgG	*	**
Human IgG2	****	****	Guinea Pig IgG****		**
Human IgG3	×	****	Bovine IgG	**	****
Human IgG4	****	****	Sheep IgG	* /×	**
Human IgGA	**	×	Goat IgG	* /×	**
Human IgGD	**	×	Pig IgG	***	***
Human IgGE	**	×	Chicken IgG	×	*
Human IgGM	**	×			
Mouse IgG1	*	**	Fragment	Protein A	Protein G
Mouse IgG2a	****	****	Human Fab	*	*
Mouse IgG2b	***	***	Human F(ab)2	*	*
Mouse IgG3	**	***	Human scFv	*	×
Mouse IgGM	* /×	×	Human Fc	**	**
Rat IgG1	×	*	Human κ	×	×
Rat IgG2a	×	****	Human λ	×	×
Rat IgG2b	×	**			
Rat IgG2c	*	**			
Rat IgGM	* /×	×			
Rabbit IgG	****	***			

**** = Strong Affinity
 *** = Moderate Affinity
 ** = Weak Affinity
 * = Slight Affinity
 × = No Affinity

Features

- **Rapid purification with high yields**

The Vivapure Protein A & G maxiprep spin columns purify > 20 mg protein from 20 ml of cleared lysate in under 50 min.

The Vivapure Protein A & G miniprep spin columns purify > 1 mg protein from 0.65 ml of cleared lysate within 20 min

- **Cost effective**

Vivapure Protein A & G maxiprep spin columns can be re-used up to 5 times

while Vivapure Protein A & G miniprep spin columns can be re-used up to 3 times

- **Run multiple spin columns in parallel**

Vivapure® miniprep and maxiprep fit in all bench top centrifuges, allowing purification of multiple samples at the same time.

Use up to 16 Vivapure maxiprep spin columns (320 ml) simultaneously

- **Premium agarose resin**

High quality agarose resin guarantees superior purity and yield in a single step

- **Negligible hold up volume**

Due to the centrifugal driving force, spin columns have negligible hold up volume in comparison to gravity columns and batch methods

Sales and Service Contacts

Europe

Germany

Sartorius Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius-stedim.com

Sartorius BBI Systems GmbH
Schwarzenberger Weg 73-79
34212 Melsungen

Phone +49.5661.71.3400
Fax +49.5661.71.3702

www.sartorius-stedim.com

France

Sartorius Stedim Biotech S.A.
Z.I. des Paluds
Avenue de Jouques – BP 1051
13781 Aubagne Cedex France

Phone +33.442.845600
Fax +33.442.845619

Sartorius Biotechnologie S.A.S.
4, rue Emile Baudot
91127 Palaiseau Cedex

Phone +33.1.6919.2100
Fax +33.1.6920.0922

Austria

Sartorius Biotech Ges.m.b.H.
Franzosengraben 12
A-1030 Vienna

Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium

Sartorius Biotech N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde

Phone +32.2.756.06.80
Fax +32.2.756.06.81

Denmark

Sartorius Nordic A/S
Himmelev Bygade 49
4000 Roskilde

Phone +45.70.23.4400
Fax +45.46.30.4030

Italy

Sartorius S.p.A.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands

Sartorius Biotech B.V.
Edisonbaan 24
3439 MN Nieuwegein

Phone +31.30.6025080
Fax +31.30.6025099

Spain

Sartorius, S.A.
C/Isabel Colbrand 10-12,
Planta 4, Oficina 121
Polígono Industrial de Fuencarral
28050 Madrid

Phone +34.91.3586102
Fax +34.91.3588804

Switzerland

Sartorius Biotech GmbH
Lerzenstrasse 21
8953 Dietikon

Phone +41.1.746.50.00
Fax +41.1.746.50.50

U.K.

Sartorius Biotechnology Ltd.
Longmead Business Park
Blenheim Road, Epsom
Surrey KT19 9 QQ

Phone +44.1372.737100
Fax +44.1372.720799

America

USA

Sartorius Biotech Inc.
131 Heartland Blvd.
Edgewood, New York 11717

Toll-Free +1.800.3687178
Fax +1.631.254.4253

Sartorius BBI Systems, Inc.
2800 Baglyos Circle
Bethlehem, PA 18020

Phone +1.610.866.4800
Fax +1.610.866.4890

Asia | Pacific

India

Sartorius Biotech India Private Ltd.
10, 6th Main, 3rd Phase Peenya
KIADB Industrial Area
Bangalore – 560 058

Phone +91.80.2839.1963|0461
Fax +91.80.2839.8262

Japan

Sartorius Biotech Japan K.K.
KY Building, 8-11
Kita Shinagawa 1-chome
Shinagawa-ku
Tokyo 140-0001

Phone +81.3.3740.5407
Fax +81.3.3740.5406

Malaysia

Sartorius (Malaysia) Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Biotech Singapore Pte. Ltd.
10, Science Park Road, The Alpha
#02-25, Singapore Science Park 2
Singapore 117684

Phone +65.6872.3966
Fax +65.6778.2494

Australia

Sartorius Australia Pty. Ltd.
Unit 17/104 Ferntree Gully Road
Waverley Business Park
East Oakleigh, Victoria 3166

Phone +61.3.9590.8800
Fax +61.3.9590.8828