


# BioPAT® Trace and BioPAT® Multi Trace

Online Measurement of  
Multiple Parameters  
From 1 to 4 Vessels



## Benefits


- Real-time online monitoring of bioprocesses parameters
  - Glucose and lactate
  - Methanol or ethanol
- Enable automated parameter feed | bleed control capability
- Single-use sensor and fluidic elements
- Compact, space-saving design
- Room temperature stable biosensors
- Wide linear measuring range

## Product Information

The control of nutrients and metabolites starts by monitoring the parameter concentration and its rate of change. Once known, then an actuator (e.g. feed pump) can affect a change and bring the parameter to the desired set-point. BioPAT® Trace and BioPAT® Multi Trace enables this control to be automated without the need for operator intervention during the bioprocess.

**BioPAT® Trace** analysis platform is designed for simultaneous online monitoring of up to two parameters in any given setup. The biosensor components measure either glucose and lactate or small-molecule alcohols, such as methanol and ethanol. The system can be used for laboratory or industrial cultivations of microorganisms and mammalian cell lines.


**BioPAT® Multi Trace** has all the measurement functionality of the dedicated system yet is capable of sampling from up to four vessels without the need for operator intervention. The frequency, relative sample volume and accuracy of the parameter concentration data allows automated parameter control. This can be achieved by either actuating a feed pump or increasing the perfusion rate to influence parameter concentration and maintain a defined steady state.



**Figure 1:** BioPAT® Trace practical setup with Univessel® SU

## Biosensor

BioPAT® Trace sample analysis works by coupling an enzymatic reaction in combination with amperometric detection. A two-channel enzyme electrode coated with an oxidase constitutes the biosensor. The biosensor's platinum electrodes are coated with glucose oxidase and lactate oxidase and react with their corresponding substrates. Alternatively, the alcohol bio sensor utilizes alcohol oxidase to detect both ethanol and methanol concurrently.



**Figure 2:** image of the back of the BioPAT® Trace biosensor

## Fluidics Tube Set

To perform the sample analysis, a sample must be transported from the vessel to the biosensor. A complete fluidic setup consisting of 1–2 mm inner diameter tubing, peristaltic pump heads and ergonomically labeled connections is needed for every run. Each sampling method – dialysis or filtration – requires the corresponding Tube set to function correctly. The BioPAT® Multi Trace Tube set design has four sample lines which are controlled by automated pinch valves to sample from up to four different vessels. For easy identification, the back of the biosensor (Fig 2) has a unique serial number printed above the electronic contact points.

## Continuous Analysis

The BioPAT® Trace provides continuous parameter analysis with up to one data point per minute during cultivation. It is independent of the type of cultivation (batch, fed-batch, continuous cultivation) and allows the setup of automated parameter control loops. In addition to the online analysis function, scheduled or manually activated parameter calibrations during a batch ensure measurement robustness without stoppage.

## Range

**The linear measuring range of the BioPAT® Trace extends from**

|                 |                        |             |
|-----------------|------------------------|-------------|
| 1 to 40 g/L     | (56 µM to 0.223 mM)    | Glucose     |
| 0.5 to 10 g/L   | (56 µM to 0.102 mM)    | Lactate     |
| 0.01 to 10 g/L  | (0.56 µM to 0.056 mM)  | Low glucose |
| 0.05 to 2.5 g/L | (0.56 µM to 0.026 mM)  | Low lactate |
| 0.5 to 20 g/L   | (6.24 µM to 0.624 mM)  | Methanol    |
| 1 to 40 g/L     | (10.85 µM to 0.868 mM) | Ethanol     |

BioPAT® Trace and Multi Trace require a minimum two-point calibration with known traceable parameter solutions to generate the linear measuring range. These calibration solutions are available from Sartorius Stedim Biotech.

## Frequency

The measurement frequency of the BioPAT® Trace filtration probe is up to 60 analyses per hour whereas; the dialysis setup yields 30 analyses per hour from a single vessel. The BioPAT® Multi Trace must perform additional flushing between vessel sampling and, therefore provides analyses every 6 (filtration)/7(dialysis) minutes, with four vessels connected. The next sequential data point generated for the first vessel is 24(filtration) or 28(dialysis) minutes later.

If the required glucose|lactate concentration measured is below 0.1 g/L in dialysis mode the accumulation time increases, resulting in biosensor data points every six minutes. This requires a BioPAT® Trace firmware update, with an optional change of the dialysis probe membrane to increase the permeability.

The BioPAT® Trace software can take scheduled samples at any given time point thus increasing or decreasing the sampling frequency to tailor it to specific process control requirements.

## Duration

The service life of the biosensor is 14 days or 5000 analyses depending on the application. The ambient temperature of the BioPAT® Trace can lie between 5°C and 35°C due to internal temperature correction. The ambient humidity should not exceed 90%.

Refer to the Glucose | Lactate: Performance and Accuracy application note for more details.

## Accuracy

The deviation from the average measurement value is less than 3% for a measurement of 5 g/L glucose and 2.5 g/L lactate. Dynamic vessel temperature compensation of the biosensor and membrane diffusion is included if BioPAT® Trace is used with BioPAT® MFCS data acquisition.

## Communication Integration

### BioPAT® DCU Analog Connection

BioPAT® Trace has standard analog outputs in either voltage or milliamp signal ranges in order to connect it to the BioPAT® DCU and to show the two measured parameters. The BioPAT® DCU is configured to interoperate the signal over the fixed calibration range and show the two concentration values on the Biostat® visual display.

BioPAT® Multi Trace is not designed for analog connection to obtain the relevant data from all four vessels. Instead, digital connection is required.

### BioPAT® MFCS digital connection

BioPAT® Trace and BioPAT® Multi Trace can be digitally connected by Modbus|OPC® to BioPAT® MFCS 3.0 bioprocess management module with additional software from Sartorius Stedim Biotech. This option allows parameter data to be imputed directly into the supervisory control and data acquisition software for all sampled cultivations.

# Mechanical Integration

## Reliable Sampling

To analyze substrates in cultivation media, it is necessary for the sample to be taken from the vessel, while maintaining the sterile conditions within. The reproducibility and relevance of the sample taken must remain intact. The BioPAT® Trace has four sampling probe options available: filtration, dialysis single-use and by-pass loop. All probe designs ensure aseptic removal from the vessel and transport to the biosensor module.

## Filtration Probe

Cell-free samples are removed from the vessel using the filtration sampling probe, which incorporates a 0.2 µm polypropylene filter. This hydrophilized probe is mounted inside the bioreactor and subsequently sterilized along with the medium. The concentration of the bulk sample is directly correlated to the linear calibration in the biosensor.

- Continuous sampling removes 1 mL from the sample vessel per sample reducing the overall vessel volume. If in interval mode 5 mL is removed for each initial measurement.



## Dialysis Probe

This probe functions by time-dependent diffusion of sample across a 12 kDa permeable membrane. The transport buffer and vessel solution do not directly mix. However, during the BioPAT® Trace sampling cycle the transport buffer flow is paused and the concentration gradient drives small molecules from the bulk into the transport buffer. When flow resumes after a set time period, this accumulated concentration is transported to the biosensor for analysis. This data is then correlated to the linear calibration to yield the bulk concentration.


- This sampling method uses the diffusion of small molecules across the permeable membrane without any loss of vessel volume.



## Single-Use Probe

The single-use probe functions on the same principles as the dialysis probe but is constructed from known traceable plastics and hollow fibers.

- This sampling method uses the diffusion of small molecules across the permeable membrane, without entailing any loss of vessel volume.



## By-pass loop

The by-pass loop functions when sterilization or a technical limitation exists that a probe cannot be installed in the vessel. It functions on the same principles as the dialysis probe with the membrane contact surface installed into a flexible tubing line with a 6 or 9mm hose-barb connection.



## Manual sample analysis

BioPAT® Trace systems enable a three-way stopcock to be conveniently integrated for manual sampling of off-line samples.

## Calibration Solutions and Transport Buffer

BioPAT® Trace and BioPAT® Multi Trace require two calibration solutions and a transport buffer to function, and accurately determine the parameter concentration in a sampled vessel. These liquids are connected on the labeled end points of the Tube set and are automatically pumped within the system loop. The calibration solutions include expiration dates to ensure the accuracy and longevity of the system analysis. Upon request, test certificates can be provided for each calibration solution.

|         |   |      |
|---------|---|------|
| BPT0006 | Transport buffer solution for CC Glc   Lac (20×);         | 0.5L |
| BPT0060 | Transport buffer Glc Lac (5x) for MO;                     | 1.0L |
| BPT0046 | Transport buffer for MeOH EtOH (5×)                       | 1.0L |
| BPT0007 | Calibration solution 5<br>(10g/L Glc and 5g/L Lac);       | 0.5L |
| BPT0008 | Calibration solution 4<br>(4g/L Glc and 2g/L Lac);        | 0.5L |
| BPT0009 | Calibration solution 3<br>(2g/L Glc and 1g/L Lac);        | 0.5L |
| BPT0010 | Calibration solution 2<br>(1g/L Glc and 0.5g/L Lac);      | 0.5L |
| BPT0011 | Calibration solution 1<br>(0.5g/L Glc and 0.25g/L Lac);   | 0.5L |
| BPT0041 | Calibration solution 6<br>(20g/L Glc and 10g/L Lac);      | 0.5L |
| BPT0043 | Calibration solution 7<br>(0.1 g/L Glc and 0.05 g/L Lac); | 0.5L |
| BPT0051 | Calibration solution, 0.2 g/L Methanol                    | 0.5L |
| BPT0052 | Calibration solution, 1.0 g/L Methanol                    | 0.5L |
| BPT0053 | Calibration solution, 5.0 g/L Methanol                    | 0.5L |
| BPT0054 | Calibration solution, 20.0 g/L Methanol                   | 0.5L |
| BPT0055 | Calibration solution, 40.0 g/L Ethanol                    | 0.5L |
| BPT0056 | Calibration solution, 4.0 g/L Ethanol                     | 0.5L |
| BPT0057 | Calibration solution, 2.0 g/L Ethanol                     | 0.5L |
| BPT0058 | Calibration solution, 0.5 g/L Ethanol                     | 0.5L |
| BPT0044 | Concentrated (200×) cleaning   disinfection solution      | 1.0L |

# Technical Data

## Specifications

### BioPAT® Trace and BioPAT® Multi Trace

|                                      |  |  |                        |
|--------------------------------------|--|--|------------------------|
| Measuring principle                  | Enzymatic amperometric   |  |                        |
| Linear measuring range               | Glucose  | 1 to 40 g/L                                    | (56 µM to 0.223 mM)    |
|                                      | Lactate  | 0.5 to 10 g/L                                  | (56 µM to 0.102 mM)    |
|                                      | low Glucose  | 0.01 to 10 g/L                                 | (0.56 µM to 0.056 mM)  |
|                                      | low Lactate  | 0.05 to 2.5 g/L                                | (0.56 µM to 0.026 mM)  |
|                                      | Methanol   | 0.5 to 20 g/L                                  | (6.24 µM to 0.624 mM)  |
|                                      | Ethanol  | 1 to 40 g/L                                    | (10.85 µM to 0.868 mM) |
| Measurement deviation                | Glucose   Lactate  | <1.5 % 1 to 20 g/L   <1.5 % 1 to 10 g/L F.S.   |                        |
|                                      | Low Glucose   Lactate  | <2.5 % 0.5 to 1 g/L   <2.5 % 1 to 0.5 g/L F.S. |                        |
|                                      | Methanol   | <2.0 % 0.5 to 20 g/L                           |                        |
|                                      | Ethanol  | <2.0 % 1.0 to 40 g/L                           |                        |
| pH range of the medium               | 4.8 – 9.2  |  |                        |
| Service life of the enzyme electrode | Glucose   Lactate  | 30 days or 5,000 analyses                      |                        |
|                                      | Methanol or Ethanol  | 15 days or 5,000 analyses                      |                        |
| Operating temperature                | 15 °C to 35 °C (59 °F to 95 °F)  |  |                        |
| Operating humidity                   | 10 to 90 %   |  |                        |
| System storage conditions            | 5 °C to 60 °C (32 °F to 140 °F);<br>5 % to 75 % RH<br>(non-condensing) |  |                        |
| Storage temperature                  | Glucose   Lactate biosensor  | 5 °C to 25 °C (41 °F to 77 °F)                 |                        |
|                                      | Alcohol biosensor  | 3 °C to 8 °C (37.4 °F to 46.4 °F)              |                        |

### BioPAT® Multi Trace

|   |  |                                |
|---|--|--------------------------------|
| Data frequency<br>(with 4 vessel connected) | Filtration                                       | Up to 2 measurements/hr/vessel |
|   | Dialysis   | Up to 2 measurements/hr/vessel |
| Weight                                      | 2.0 kg   |                                |
| Dimensions (W × H × D)                      | 120 mm × 170 mm × 200 mm<br>(4.7" × 6.7" × 7.9") |                                |
| Serial output                               | 9-pin female connector to<br>external pump       |                                |
| USB interface                               | 1 × USB port (software<br>updates only)          |                                |

## BioPAT® Trace

|                        |  |  |
|------------------------|--|--|
| Data frequency         | Filtration<br>Dialysis                           | Up to 60 measurements/hr<br>Up to 30 measurements/hr |
| Weight                 | 1.8 kg   |  |
| Dimensions (W × H × D) | 120 mm × 170 mm × 200 mm<br>(4.7" × 6.7" × 7.9") |  |
| Analog output options  | 0 to 10 V   0 to 20 mA   4 to 20 mA              |  |
| Serial output          | 9-pin female connector to external pump          |  |
| USB interface          | 1 × USB ports (software updates only)            |  |

## Connections

|                                 |   |
|---------------------------------|---|
| System communication interfaces | RS232, Ethernet   |
| Communication cable length      | Standard: 2 m (5 m and 10 m available)  |
| Communication cable material    | Non-metallic sheathed cable with<br>finely stranded copper conductor;<br>PVC insulation -40 °C to +70 °C<br>(+40 °F to +158 °F) |
| Power supply                    | 100 to 120   220 to 240 V ~; 50   60 Hz   |

## Liquid | Waste Requirements

|                                  |                                 |
|----------------------------------|---------------------------------|
| Transport buffer consumption     | 1.5 L/day to 2.0 L/day          |
| Calibration solution consumption | 5 mL/day to 145 mL/day          |
| Fluid temperature                | 15 °C to 35 °C (59 °F to 95 °F) |
| Collection of waste liquid       | Separate container required     |

## BioPAT® Trace Filtration probe

|   |   |  |
|---|---|--|
| Mechanical port connection  <br>Insertion depth | 12 mm <sup>1</sup> PG 13.5   130 mm<br>12 mm <sup>1</sup> PG 13.5   230 mm<br>25 mm <sup>1</sup> DN 25   130 mm<br>25 mm <sup>1</sup> DN 25   90 mm | BPT0021<br>BPT0022<br>BPT0023<br>BPT0036 |
| Port adapter   insert depth reduction           | 19/12 mm   49 mm<br>25/12 mm   52 mm  | BB-8848630<br>BB-34165225                |
| Product contact materials                       | Construction material: 1.4404/AISI, 316L stainless steel – mechanically polished<br>Filter: Polypropylene, 0.2 µm                                   |  |
| Pressure range                                  | 0.8 to 5 bar   11.6 to 72.52 psi absolute pressure  |  |
| Port gasket                                     | O-ring, 11.00 mm × 3.00 mm EPDM (FDA   USP Class VI)  |  |
| Vessel preparation requirements                 | filter install – hydrophilization – autoclave   steam sterilization   |  |

<sup>1</sup> steel conduit armored thread

## BioPAT® Trace Dialysis Probe

|                                       |  |   |
|---------------------------------------|--|---|
| Mechanical port connection            | Diameter: 12 mm thread: steel conduit, PG 13.5   |   |
| Insertion depth   Univessel® volume   | 132 mm   0.5L<br>165 mm<br>212 mm   ½L<br>232 mm   2L<br>332 mm   5L<br>362 mm   5L<br>432 mm   10L  | BPT0014<br>BPT0015<br>BPT0016<br>BPT0017<br>BPT0018<br>BPT0019<br>BPT0020 |
| Port adapter   insert depth reduction | 19/12 mm   49 mm<br>25/12 mm   52 mm   | BB-8848630<br>BB-34165225   |
| Product contact materials             | 1.4404 AISI – 316L stainless steel – Mechanically polished   |   |
| Construction material                 | Dialysis membrane: Cellulose acetate, approx. 12 kDa MWCO <sup>2</sup><br>Pressure range: 0.8 bar to 3.4 bar   11.6 psi to 49.91 psi absolute pressure |   |
| Port gasket                           | O-ring, 11.00 mm × 3.00 mm EPDM (FDA   USP Class VI)   |   |
| Vessel preparation requirements       | install membrane – fill buffer – autoclave   steam sterilization   |   |

## BioPAT® Trace Single-use Probe

|                           |  |
|---------------------------|--|
| Connection                | Integrated into Flexsafe® STR® or RM bag sample port   |
| Product contact materials | Flexsafe® RM or STR® Validation guide                  |
| Hollow-fiber membrane     | Polyethersulfone PES, approx. 12 kDa MWCO <sup>2</sup> |
| Preparation requirements  | Ready to use   |

<sup>2</sup> Molecular weight cut-off

## BioPAT® Trace Tube set for Filtration


|                                   |                        |
|-----------------------------------|------------------------|
| BioPAT® Trace                     |                        |
| Inlets                            | 4                      |
| Outlets                           | 1                      |
| BioPAT® Multi Trace               |                        |
| Inlet                             | 7                      |
| Outlets                           | 1                      |
| Tube sets for BioPAT® Trace       |                        |
| Glc   Lac                         | BPT0004                |
| Glc   Lac cellulase stable        | BPT0005                |
| MeOH   EtOH                       | BPT0049 and<br>BPT0048 |
| Tube sets for BioPAT® Multi Trace |                        |
| Glc   Lac                         | BPT0M02                |
| Glc   Lac cellulase stable        | BPT0M03                |
| MeOH   EtOH                       | BPT0M04 and<br>BPT0048 |

## BioPAT® Trace Tube set for Dialysis


|                                   |                        |
|-----------------------------------|------------------------|
| BioPAT® Trace                     |                        |
| Inlets                            | 3                      |
| Sample return loops               | 1                      |
| Outlets                           | 1                      |
| BioPAT® Multi Trace               |                        |
| Inlet                             | 3                      |
| Sample return loops               | 4                      |
| Outlets                           | 1                      |
| Tube sets for BioPAT® Trace       |                        |
| Glc   Lac                         | BPT0003                |
| MeOH   EtOH                       | BPT0050 and<br>BPT0048 |
| Tube sets for BioPAT® Multi Trace |                        |
| Glc   Lac                         | BPT0M01                |
| MeOH   EtOH                       | BPT0M05 and<br>BPT0048 |




**BioPAT® Trace**  
**Tube set for Filtration**




**BioPAT® Trace**  
**Tube set for Dialysis**



BioPAT® Multi Trace  
Tube set for Filtration



BioPAT® Multi Trace  
Tube set for Dialysis



# Certificates

## Calibration certificates

**Acceptance Report**  
**Calibration Solution**  
**4 g/l Glucose**  
**2 g/l Lactate**  
**Part.-No.: 850.305.003**



**TRACE**  
Analytics

Richard-Wagner-Str. 1-2  
38106 Braunschweig  
Germany  
Phone: ++49-531-209008-0  
Fax: ++49-531-209008-39  
email: [info@trace.de](mailto:info@trace.de)  
[www.trace.de](http://www.trace.de)

### Certificate of Conformity:

We confirm the agreement of this acceptance report with the original examination report of arotop food environment GmbH, Mainz, December 19, 2011.

TRACE Analytics GmbH

### Examination Report:

Calibration Solution 4 g/l Glucose / 2 g/l Lactate

|                     |   |
|---------------------|---|
| Project-No.:        | 045369  |
| Lot-No.:            | 2011-40376  |
| Date of receive:    | December 09, 2011   |
| Parameter:          | D-Glucose   |
| Result:             | 3.97 g/l  |
| Standard deviation: | ± 0.014 g/l   |
| Parameter:          | L-Lactic Acid   |
| Result:             | 2.007 g/l   |
| Standard deviation: | ± 0.066 g/l   |
| Method:             | R-Biopharm, Boehringer Mannheim, Enzymatic<br>BioAnalysis / Food Analysis<br>According to §64 LFGB, ASU L31.00-12<br>The detection limit of glucose is 4.0 mg/L<br>The detection limit of Lactate is 5.0 mg/L |

Date of Acceptance

Signature

The original report was performed by

arotop food environment GmbH  
Institut für Geschmacksforschung, Lebensmittel- und Umweltanalytik  
Dekan-Leist-Straße 9  
D-55129 Mainz


Accredited in accordance with DIN EN ISO/IEC 17025.

**Germany**

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49 551 308 0

**USA**

Sartorius Stedim North America Inc.  
565 Johnson Avenue  
Bohemia, NY 11716  
Toll-Free +1 800 368 7178

 For further contacts, visit  
[www.sartorius.com](http://www.sartorius.com)

Specifications subject to change without notice.

© 2021 Sartorius Stedim Biotech GmbH, August-Spindler-Strasse 11, 37079 Goettingen, Germany

Publication No.: SLL2006-e | Order No.: 85032-541-03 | Status: 09 | 21 | 2021