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Contents | Inhalt

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
Table of Contents

Sterisart® NF gamma Septum – Overview	4
Product description	5
Packaging	5
Intended application	6
Specification	6
1. Prepare the sterilization area of an isolator	7
2. Perform these steps before filtration of the samples	8
3. Sampling from prefilled syringes	9
4. After filtration of the samples and the rinse washing solution	10
5. Fill containers with liquid culture media	11
6. Incubate the containers according to the respective regulations or conditions	12
7. Sterile sample taking during incubation	12

Sterisart® NF gamma Septum – Overview

Instructions for Use
in Sterility Testing According to
International Pharmacopeias

Pos.	Designation
1	Pre-installed tube clamp
2	Vent filter
3	Sterisart® container
4	Connector with septum for sterile sample taking
5	Wing nut plug
6	Rubber cap
7	Double-needle metal spike for closed containers, modified
8	Silicon tubing
9	Syringe filter holder, adapter for prefilled syringes
10	Double tubing
11	Additional plastic bag

Scope of the instructions for use

These instructions are applicable to the following products:

16469 GSD

Sterisart® NF gamma Septum, system with adapter for prefilled syringes (pack of 10, double-packed, gamma irradiation sterilized)

Product description (see page 4)

The Sterisart® NF systems listed are ready-to-use assemblies. Each system consists of two clear containers that are connected by flexible double tubing (for transferring equal volumes of the sample solution into the two containers) and a Y-distributor equipped with a holder for prefilled syringes. A double-needle metal spike is connected to the holder by a single tube with clamp.

Each container has a validated 50-mm diameter. 0.45 µm pore size membrane filter sealed firmly at its base to prevent the diffusion of liquid into the filter-container boundary area. A validated sterilizing vent filter and a connector with septum for sterile sample taking are fixed to the top of the container. Graduated marks on each container indicate approximately 50, 75 and 100-ml levels.

The adapter is designed to hold prefilled syringes, the contents of which are drained into the Sterisart® NF tubing system and, at the same time, mixed with rinsing| washing solution that is suctioned through the double-needle metal spike from an additional closed container. Equal volumes of the liquids mixed in this way are transferred through the double tubing system into the two Sterisart® containers.

Prefilled syringes with needles are rinsed with the washing solution and can thereby also be tested for sterility. The single tube attached to the double-needle metal spike can be closed off with a (white) tube clamp.

Additional color-coded clamps (two white, two yellow) pre-mounted on the double tubing allow for convenient closing of the tubing as required. All the accessories needed (two rubber caps for the vent filters, two wing nut plugs for the container outlets) are included in a separate package with each system.

Packaging

Sterisart® NF gamma Septum:

A box consists of a Sterisart® NF system packed with appropriate accessories in a gas-impermeable, clear plastic tray which is sealed on top by a gas-impermeable cover sheet. The print on the cover gives details of the enclosed system (brief description, order number, expiration date, sterilization date and lot number).

The gas-impermeable plastic box enable the outer surfaces to be sterilized by usual isolator procedures without passage of sterilant (e.g., H₂O₂) through to the box and the Sterisart® NF system. A change from orange to red of the sterilization indicator located on the outer standard box as well as on the accessory package inside the individual system shows that the product has been reliably sterilized by gamma-irradiation.

A standard box contains ten individually packaged systems, and has a sterilization indicator. A purple dot inside and on the end of the box additionally identifies type 16469.

Intended application

Sterisart® NF gamma Septum systems have been specifically designed for sterility testing using the membrane filter method described in pharmacopoeias, in combination with a dedicated peristaltic pump. Type 16469 allows sterility testing of prefilled syringes of different sizes (versions with and without needles). Sterisart® NF gamma Septum has been optimized for sterility testing in isolators.

Important notes on safe and correct usage

Sterisart® NF systems are only to be used with double-tubing peristaltic pumps designed for sterility testing with such systems, such as the Sterisart® Universal pumps 16419 and 16420.

All necessary precautions must be taken to avoid secondary contamination. Handle the parts only when wearing sterile gloves.

The needles of the double-needle metal spike and the bevels of the prefilled syringes are very sharp. Take great caution to avoid piercing gloves or injuring hands when handling them.

Check the color indicator and the expiration date on the box before removing a Sterisart® NF system from its packaging.

To avoid media leakage during incubation, securely attach the plastic wing nut plug to the outlet on each Sterisart® NF container by turning the plug clockwise. This will prevent leakage of the culture media during incubation.

Never "glue" the rubber caps onto the vent filters! These caps are intended to serve not only as closures, but also as safety valves. Should exceptionally high pressure build up in the containers during pumping (high pump speed or clogged membrane filters), normally fitted caps will blow off before the burst pressure of other system parts is reached.

Specifications


Max. operating pressure:	3 bar at 20°C
Min. operating temperature:	50°C
Chemical compatibility:	Details available on request

Please follow the instructions for use starting on the next page.


1. Prepare the sterilization area of an isolator

Sterisart® NF gamma Septum:


Open the standard package and the plastic bag to remove the quantity of individual packs required for the number of sterility tests to be performed.

Place the individual packs, still enclosed in the gas-impermeable plastic box, in the sterilization area.

Sterilize the outer surfaces of the box according to the isolator validation specifications, for example, with H₂O₂.


Open the box or, if appropriate, pass the closed Sterisart® NF boxes through the air-lock port on the isolator.


2. Perform these steps before filtration of the samples

Sterisart® NF gamma Septum:


Cut open the tray and take out the Sterisart® NF system.
Position the two Sterisart® NF containers into the pump holder.

Thread the Sterisart® NF double tubing in the pump head, following the instructions given for the pump used.


Tightly seal the two vent filters using the enclosed rubber caps.


3. Sampling from prefilled syringes


Make sure that all tube clamps are open.

Detach the plug from the prefilled syringe holder adapter on the Y-distributor.


Insert the prefilled syringe with or without needle.


Insert the double-needle metal spike as far as necessary into the port of the rinsing/washing solution container.

Switch on the pump to transfer equal volumes of rinsing solution into the two Sterisart® NF containers and wet the membrane in the containers. At the same time, sterile air is drawn in through the double-needle metal spike.

After wetting the membrane, inject the contents of the syringe into the system.


After injecting the contents, do not close the clamp on the single tube until the pump has been switched off.

Changing the prefilled syringes:


Open the clamp on the single tube and then start the pump. Perform sampling from other prefilled syringes according to the same procedure.

For rinsing and filling with culture media, the last syringe to be tested remains on the holder. The tube clamp on the single tube should always remain open.

4. After filtration of the samples and the rinse|washing solution


Seal the outlet of each Sterisart® NF container using the enclosed wing nut plugs. Press down on each plug, turning it clockwise to seal correctly.

Remove the rubber caps from the vent filters.


5. Fill containers with liquid culture media

Culture medium 1 (thioglycollate broth):

Close the white tube clamp at the outlet on the Y-distributor.

The clamp on the single tube remains open.


For filling with a double-needle metal spike:

Insert the double-needle metal spike into the port of the closed medium container.

Switch on the pump.

Switch off the pump when the predefined amount has been transferred.


Culture medium 2 (tryptic soybean broth):

Open the white tube clamp and close the adjacent yellow clamp.


The clamp on the single tube remains open.

Transfer the culture medium as described before.


Seal off the two tubes above the Sterisart® NF containers using the pre-installed tube clamps. Cut off the tubing.

6. Incubate the containers according to the respective regulations or conditions

Aerobic incubation and anaerobic incubation.

7. Sterile sample taking during incubation

Remove the transparent plug from the septated connector. Insert the sterile needle tip in the septum and remove samples of the medium inside the container. Next, remove the needle tip from the septum. The septum guarantees a sterile barrier between the outer and inner side of the container during further incubation.

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Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part of
this publication may be reprinted
or translated in any form or
by any means without the prior
written permission of Sartorius
Stedim Biotech GmbH.
The status of the information,
specifications and illustrations
in this manual is indicated by
the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make changes
to the technology, features,
specifications and design of the
equipment without notice.

Status:
November 2010,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany