



Sartorius Pipetting and Dispensing Products

Simplifying Progress

SARTORIUS



Table of Contents

The Sartorius product portfolio, technical data and specifications are subject to change, due to continuous product development.

Please visit the website for an up-to-date overview of Sartorius products:

[sartorius.com](https://www.sartorius.com)

About Sartorius Pipetting and Dispensing

Introduction	4
Quality	6

How to Choose a Pipette

Electronic or Mechanical Pipette	8
----------------------------------	---

Electronic Pipettes

Picus® 2 Electronic Pipettes	12
------------------------------	----

Mechanical Pipettes

Tacta® Mechanical Pipettes	18
Proline® Plus Mechanical Pipettes	23

Stands and Accessories

Pipette Stands	30
Safe-Cone Filters	31
Adjustment Tool	33

Pipette Tips

Pipette Tips	36
Optifit Tips	38
Safetyspace® Filter Tips	39
Packaging Options	40

Maxi-Volume Liquid Handling

Midi Plus Pipetting Controller	48
Prospenser Plus and Prospenser	50
Biotrate Digital Burette	52

Pipetting Academy

Pipetting Academy	58
Pipetting Recommendations	60

Calibration and Maintenance

Pipette Calibration and Maintenance Services	64
Pipette Decontamination Procedure	65
Autoclaving Instructions	66

Sales and Service Contacts	68
----------------------------	----

Introduction

Three key factors – ergonomics, design and reliability – form the cornerstone of our product development. These and other factors have been combined to produce a perfectly balanced mechanical pipette, the Tacta®, and the Picus® 2, which is an excellent example of an electronic pipette that has all these aspects combined. Tacta® and Picus® 2 are among the lightest pipettes on the market, reducing the risk of repetitive strain injury (RSI). Their high reliability, and ease of use, make them valued instruments for professionals, who strive for high quality results.

Their functional and well-rounded design has been recognised with design awards. All Sartorius pipettes are designed and manufactured in Finland, where our R&D team is constantly seeking solutions to further improve liquid handling instruments to make lab work easier.

Ergonomics

When designing a pipette, we always consider the shape and function of the human hand. As we understand the risks of repetitive pipetting, we emphasize ergonomic design in every product we make. Simply put, this means products that you can use in a comfortable posture with minimum muscle power. Our pipettes and dispensers are designed for both right- and left-handed users. Their operating buttons are located sufficiently close together, within ergonomic reach of the thumb.

Design

We provide products with a timeless and light, yet practical, design, suitable for laboratory settings and pleasing to the eye of the user. The Tacta® and the Picus®, won the Red Dot design award, and the German Design Award. The Picus® was also distinguished with the Fennia Prize Honorary Mention.

Responsible Manufacturing

At our manufacturing sites, we have invested in production technologies that generate less waste. Our pipette tip production facility in Kajaani, Finland uses 100% renewable electricity and has achieved a waste recovery rate of 98%. Additionally, plastic waste from Sartorius pipette tip manufacturing is recycled for use as raw materials for other plastic products.

We are continuously working to reduce our environmental impact. To further drive this, our environmental management system in Finland is certified by an independent third party to meet the ISO 14001 standard. Our facilities in Finland work with a Carbon Neutral Waste Management service to offset their waste's carbon footprint through certified afforestation projects. In 2023, investments were made in permanent carbon sinks to neutralize waste emissions.

Read more about product sustainability at Sartorius:
<https://www.sartorius.com/en/company/sustainability/product-sustainability>



Reliability

For us, reliability has many aspects, the most important being accuracy and precision of results and secured purity.

The core of a pipette lies in its **accuracy and precision**. For this reason, we have used the newest technologies together with in-house innovations, to achieve even more reliable pipetting results. Our electronic brake, piston control system and plate tracker for electronic pipettes are our latest innovations. They increase accuracy, precision and reliability of the device. Another important factor in achieving reliable results is the optimal tip fit, which we can guarantee by designing and producing the tips ourselves, so that they perfectly match our pipettes.

To reduce the risk of contaminating the internal components of our pipettes, we offer special Safe-Cone Filters to be used in our pipettes, as we understand that purity is a key concern in many laboratories. We strive to produce as many autoclavable products as possible, both pipettes and tips. Our pipette tips are manufactured in ISO Class 8 Cleanroom conditions. We test every certified tip lot for DNase, RNase and endotoxins at an external laboratory. We also offer an innovative Safetyspace® Filter Tip range for safer and contamination-free pipetting.



Quality

Sartorius' products are developed and manufactured according to the requirements of the ISO 9001, ISO 13485 and ISO 14001 quality and environmental standards. Tip production also follows the ISO 14644-1 standard, in order to fulfil ISO class 8 cleanroom conditions.

We continuously develop our products and processes in order to meet, and often exceed, the demands of regulatory authorities, environmental bodies, and most importantly, our customers.



Sartorius' products are developed and manufactured according to the requirements of the ISO 9001, ISO 13485 and ISO 14001 quality and environmental standards. Tip production also abides by the ISO 14644-1 standard, in order to fulfil ISO Class 8 Cleanroom conditions. ISO 13485 is a specific standard for medical device quality systems, and supplementing the more generic ISO 9001 standard, which applies to many industries.



Sartorius service centers worldwide with accredited pipette calibration laboratories offer calibrations services with precise technical requirements. Their accreditation status has been granted by their national accreditation bodies. Please check [sartorius.com](https://www.sartorius.com) for more information.



During production every pipette is quality tested. Pipette performance testing is carried out according to ISO 8655 standard and delivered with Quality Control certificate. Sartorius accredited pipette calibration at service follows the ISO 17025 standard.



Sartorius offers a 2-year warranty for all mechanical and electronic pipettes. The low lifetime cost and environmental friendliness of our products, which have long warranty periods, give a high return on investment.



The ergonomic design label indicates products, which Sartorius has designed specifically to reduce the risk of work-related hand, arm and shoulder disorders, such as repetitive strain injury (RSI).



The Optiload tip loading mechanism developed by Sartorius in Tacta®, Proline® Plus, and Picus® 2 pipettes allows tips to be loaded with constant force. This secures optimal tip sealing and minimum tip ejection force.



The Optilock®, volume locking system in Tacta® gives the choice of locking and unlocking the volume the traditional way, using both hands, or specially developed convenient method, using one hand.



The Optiject®, is a unique mechanism in Tacta®, where the soft, levered tip ejection feature lets the tip detach in a smooth, controlled way with little force. Additionally the feature allows for hands-free ejection of the Safe-Cone filter for a truly safe disposal of contaminated filters.



Every lot of Sartorius Single Tray, Refill Pack and FlexiBulk® tips are certified to be free of DNase, RNase, human DNA and endotoxins, for the protection of samples from contamination. These lot specific and other certificates can be downloaded from [sartorius.com](https://www.sartorius.com). Sartorius' tip production is ISO 8 cleanroom classified, which ensures a contamination-free manufacturing environment, and products.



Most Sartorius pipetting and dispensing products are autoclavable. Please see details in the following product specific chapters.

How to Choose a Pipette

Electronic or Mechanical Pipette

Are you looking for a pipette for sterile work, or one you could easily calibrate yourself? Or do you seek a really light and ergonomic solution? Perhaps you need a pipette with a certain pipetting mode to speed up your work? By consulting the tables below, you can choose the instrument that best suits your needs.

Electronic or Mechanical Pipette

Features	Electronic Pipettes	Mechanical Pipettes
Highest Ergonomics	▪	
Fastest pipetting	▪	
User-independent results	▪	
Multiple pipetting modes	▪	
Fully autoclavable		▪
Adjustment by user	▪	▪



Electronic Pipettes

Features	Picus® 2
Most ergonomic	▪
Weight ¹	100g
Length ¹	210 mm
Volume range, single-channels	0.5 µL – 10 mL
Volume range, multi-channels	0.5 µL – 1.2 mL
Language options	5
Pipetting modes	8
Advanced functions	8
Repeated blow-out (advanced function)	▪
Microwell plate tracker	▪
Connectivity	Bluetooth, USB
Memory places (for storing programs)	20
Reminders for calibration & service	▪
Password protection	▪
Electronic tip ejection	▪
Calibration adjustment by user	▪
Calibration adjustment in 1, 2 or 3 points	▪
Hot key for stored programs	▪
Use of pipette while charging	▪
Fully charged in 1 hour	▪
Safe-Cone Filters	▪
Autoclavable lower parts ²	▪
Optiload in multi-channels	▪
Color coding on pipette	▪
Warranty for 2 years	▪

¹ 1,000 µL 1-channel models

² Excluding 1,200 µL multi-channel pipettes

Mechanical Pipettes

Features	Tacta®	Proline® Plus
Most Ergonomic	▪	
Ergonomic Finger Hook	▪	▪
Weight ¹	75g	82g
Length ¹	225 mm	239 mm
Volume range, single-channels	0.3 µL – 10 mL	0.3 µL – 10 mL
Volume range, multi-channels	1 – 300 µL	1 – 300 µL
Fixed-volume models		▪
Pipetting force ¹	12 N	15 N
Optiject soft tip ejection	▪	
Light tip ejection		▪
Optiload spring-loaded tip cones	all models	multi-channels only
User adjustment	▪	▪
Optilock on off volume lock	▪	
Volume locking	▪	click stops
Big, and easy to read display	▪	▪
Safe-cone Filters (models > 10 µL)	▪	▪
Filter ejector	▪	
Color coding on pipette	▪	▪
ID tags	▪	
Fully autoclavable	▪	▪
Multipacks	▪	▪
Pipette holder with pipette	▪	▪
Warranty for 2 years	▪	▪

¹ 1,000 µL 1-channel models



Tips
DNA
ee

10 mL

SARTORIUS
2029-03
REF: LH-799000
LOT: 1234567890

Optifit Tips

SARTORIUS
2027-08
REF: 791000
LOT: 501507590

DNase, RNase, DNA,
endotoxin free

1,000 µL

Electronic Pipettes

Table of Contents

Electronic Pipettes Picus® 2 Electronic Pipettes

12



Picus® 2 Electronic Pipettes

The Most Sophisticated and Ergonomic Pipettes Ever!

The Picus® 2 pipettes ensure reliable and repeatable pipetting results and feature an unbeatable ergonomic design that is kind to your hand. Picus® 2 is suited to anyone working in the lab, from a graduate student to an experienced laboratory professional. It is as intuitive to use as a mechanical pipette, but also offers more advanced options for experienced users.

Picus® 2 saves you time in the lab with a wide selection of pipetting modes and customizable programs to meet your every need. The most frequently used modes with specified settings can be saved on your device for quick and easy access. Safety features like password protection and calibration reminders support compliance requirements in regulated labs. By connecting Picus® 2 to your mobile device, you can use the Sartorius pipetting mobile app to smoothly run sample preparation workflows and adjust the pipette setting automatically, taking your productivity to the next level.

Picus® 2 is available in single-channel models, covering a volume range of 0.5 µL to 10,000 µL, and in multi-channel models from 0.5 µL to 1,200 µL.

Superior Ergonomics

The uniquely lightweight and streamlined design of the Picus® 2 pipette ensures an effortless pipetting experience. The pipette rests lightly in the user's hand due to the comfortable handle and finger hook, and minimal gripping force is needed.

The conveniently located soft-touch operating buttons, and electronic tip ejection, help minimize muscle strain for both left- and right-handed users, further enhancing the ergonomics of the pipette.

Reliable Results

The precise electronic pipetting system of the Picus® 2 pipettes enable accurate and precise pipetting results, independent of the user. Using the unique plate tracker increases reliability in microwell plate pipetting, by guiding the user to pipet into the correct microplate wells.




Features and Benefits



Picus® 2 lightweight pipettes, with a comfortable handle and finger hook, requiring minimal gripping force.



Adjustment wheel for quick single-handed operation and speed control.



The Sartorius Pipetting app is your handy companion for sample preparation workflows and pipette management.

Ergonomic Design for Reduced Risk of Strain Injury

- Extremely compact and light (only 100 g) design maximizes user comfort
- Conveniently located, soft-touch operating button and comfortable electronic tip ejection help minimize muscle strain and reduce the risk of repetitive injury
- Comfortable handle design and finger hook allow the pipette to rest effortlessly in hand

Sartorius Pipetting App

The Sartorius pipetting app is your trusty companion for sample preparation workflows, as well as managing Sartorius pipettes. For Picus® 2 pipette users it is a must-have.

This intuitive and user-friendly app is designed to provide step-by-step guidance for various workflows for Picus® 2 pipettes. Whether you are a seasoned scientist or a novice in the field, this app offers a comprehensive and easy-to-use platform for ensuring the accuracy and reproducibility of your liquid handling workflows. The app features an ever-growing database of configurable protocols that can be adapted according to your specific application.

Innovative Technology for Reliable Results


- Fully electronic functionality provides outstanding accuracy and repeatability of results
- Electronic brake stops piston movement rapidly and precisely, ensuring high precision, especially in serial dispensing
- Repeated blow-out function helps dispense the very last droplets of liquid, ensuring complete recovery

Intuitive User-Interface for Ease of Use

- Natural to use buttons
- Easy menu navigation and unique adjustment wheel make for fast setting of volumes and modes
- Hotkeys for quick access to your preset programs
- Enables ergonomic single-handed operation

Comprehensive Range of Pipetting Modes for Speed of Use

- The main modes combined with the advanced functions enable fast and handy execution of various pipetting tasks
- Multi-dispensing mode replaces the need for a stepper
- Save time in the lab with a wide selection of pipetting modes and customizable programs for every need



Microwell plate tracker in use

Microwell Plate Tracker for Improved Efficiency and Reliability

- Unique built-in tracker for 96 and 384 well plates guides the user to pipette into the correct wells
- Improves work efficiency and reliability of results
- Also available in Sartorius Pipetting app

Optiload in Multi-channel Models for Perfect Tip Sealing

- Allows tip loading with an equal constant force onto every channel
- Enables perfect tip sealing onto every individual tip cone

Safe-Cone Filters to Help Prevent Contamination

- Are available for all models >10 µL

Easy to Clean

- Easy to open for deep cleaning
- The pipette lower parts can be decontaminated by autoclaving (excl. 1,200 µL models)

Calibration Adjustment

- Define and save 5 calibration settings for different liquids, in addition to the factory setting
- The adjustment can be made in 1, 2 or 3 points by the user

Fast Charge that Lasts Several Hours

- Battery is charged in 1 hour
- Continue pipetting while charging with micro USB cable
- Most convenient way is to store your pipette in a charging stand

Safety Features

- Password protection helps to prevent unauthorized changes or use of the pipette
- Calibration and service reminders to support maintenance routines
- Pipette locking feature can be used when pipette is awaiting service or cleaning

Pipetting Modes

Pipetting Modes	Advanced functions (ADV)*							
	Plate Tracker	Mixing Before Aspiration	Mixing After Dispensing	Counter	Excess Volume Adjustment	Auto Dispensing	Repeated Blow-out	Fast Dispensing
Pipetting	■	■	■	■			■	
Multi-Dispensing	■	■		■	■	■		
Multi-Aspiration							■	
Reverse pipetting	■	■		■	■			
Sequential Dispensing		■		■	■	■		
Manual pipetting							■	
Diluting			■				■	
Titration								■

* Advanced functions are used in conjunction with the pipetting mode

Ordering Information

Picus® 2 Electronic Pipettes

Tip Selection Guide

Picus® 2 Order Code	Channels	Volume Range	Increment	Pipettes and Tips	Optifit Tip	Safetyspace® Tip	Safe-Cone Filters	
		(µL)	(µL)	Color-Code	Volume	Volume	Standard	Plus
LH-747021	1	0.5 – 10	0.01	■	10 µL	10 µL	–	–
LH-747041	1	5 – 120	0.10	■	200 µL	120 µL	721008	721018
LH-747061	1	10 – 300	0.20	■	350 µL	300 µL	721007	721017
LH-747081	1	50 – 1,000	1.00	■	1,000 µL	1,000 µL	721006	721016
LH-747101	1	100 – 5,000	5.00	■	5,000 µL	5,000 µL	721005	721015
LH-747111	1	500 – 10,000	10.00	■	10,000 µL	–	721005	721015
LH-747321	8	0.5 – 10	0.01	■	10 µL	10 µL	–	–
LH-747421	12							
LH-747341	8	5 – 120	0.10	■	200 µL	120 µL	721008	721018
LH-747441	12							
LH-747361	8	10 – 300	0.20	■	350 µL	300 µL	721007	721017
LH-747461	12							
LH-747391	8	50 – 1,200	1.00	■	1,200 µL	1,200 µL	721006	721016
LH-747491	12							

For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at sartorius.com





Mechanical Pipettes

Table of Contents

Mechanical Pipettes

Tacta® Mechanical Pipettes

18


Proline® Plus Mechanical Pipettes

23



Tacta® Mechanical Pipettes

Perfectly Balanced



Tacta's® Optilock feature helps to prevent accidental volume changes during pipetting.



Easy adjustment ensures accurate results with a variety of liquids. The handy adjustment key is included in the package.



Thanks to the ergonomically designed finger hook, your Tacta® rests comfortably in your hand without extensive gripping force.

Have you ever considered the overall user experience in pipetting? We have.

Sartorius proudly presents Tacta®, the premium mechanical pipette with superb comfort and reliability. Tacta® makes pipetting effortless and safe, while producing accurate and reliable results time after time.

At Sartorius, we listened to our customers and combined their views with our extensive experience, and state-of-the-art R&D, to create the exceptionally balanced design of Tacta®. Tacta® is designed to feel great in your hand, and it is easy and comfortable to use. All materials used for Tacta® have been carefully selected, and each component is designed to meet the highest standards.

Superior Ergonomics

Tacta® rests lightly in your hand thanks to its ergonomically designed handle and finger hook. Tacta® is effortless to use and exceptionally comfortable to hold. Tacta® eases your workload and protects you from strain, even when you pipette for extended periods of your working day. The unique handle and finger hook design lets the pipette rest lightly in your hand, with no need to grip the handle tightly. Tacta® is effortless to use and exceptionally comfortable to hold. Tacta® eases your workload and protects you from strain, even when you pipette for extended periods of your working day. The unique handle and finger hook design lets the pipette rest lightly in your hand, with no need to grip the handle tightly.

Low Pipetting, Tip Loading and Ejection Forces

The forces needed for a pipetting cycle, from tip attachment to pipetting and finally tip ejection, are exceptionally low with Tacta®, protecting you from repetitive injury. The Optiload feature, with spring-loaded tip cones in both the single and multichannel models, secures tip loading with perfect sealing and minimal force. Consequently, Optiload also ensures light tip ejection to protect your hand. The light tip ejection is further enhanced by Sartorius Optiject, the soft, levered tip ejection feature that lets the tip detach in a smooth, controlled way with little force.

Large, Clear Display

The large, easy-to-read display helps you to see all four digits of the set volume from a distance without straining your eyes. The volume is easy to read even when the pipette is angled – eliminating the need to turn your head into an uncomfortable position.

Adjustment

Tacta® is very easy to adjust, for various types of liquids, using a simple adjustment key. The adjustment scale shows the degree of adjustment. By noting this value for a specific liquid, you can return to that setting any time. Easy calibration adjustment ensures accurate results with a variety of liquids.

Volume Adjustment

Use the Sartorius Optilock feature to prevent accidental volume changes during pipetting. Either hold down the volume lock button with one hand and adjust the volume with the other, then release to lock, or with just one hand, slide the volume lock up, adjust the volume and slide down to lock again. The choice is in your hands.



reddot award 2016
winner



Unique to Tacta®, is its Safe-Cone Filter ejector, which enables the removal of used filters without tweezers or human contact with contaminated filters.

Safe-Cone Filters

The exchangeable Safe-Cone Filters, used in pipette tip cones, act as barriers to reduce the risk of contaminating the internal components of the pipette. The unique Safe-Cone Filter ejector enables the removal of used filters without human contact with contaminated filters. Safe-Cone Filters are available for all Tacta® models greater than 10 µL and offer a cost-effective method to reduce contamination. Filters should be changed regularly, and always after over-aspiration.

Easy Cleaning

Tacta® is the quickest and easiest pipette, on the market, to clean, with only three parts to be disassembled for cleaning and no tools are needed for disassembly. Tacta® can also be steam-sterilized or autoclaved as-is, without disassembly. It has also high UV and chemical resistance.

Features and Benefits

Feel the Comfort

- Comfortable to hold due its low weight and ergonomically designed handle
- Low pipetting and tip ejection forces that reduce the risk of repetitive injury
- Controlled and smooth tip ejection with the new Sartorius Optiject technology
- The unique Sartorius Optilock system provides flexibility for volume adjustment and locking

Achieve Reliable Results

- Reliable results, even over lengthy pipetting periods
- Calibration adjustment to provide accurate results for various liquid types
- 4-digit display for accurate and easy volume setting

Safe from Contamination

- Safe-Cone Filters provide cost-effective contamination prevention
- Easy to clean, with only three parts to disassemble
- Fully autoclavable



Ordering Information

Tacta® Mechanical Pipettes

Tip Selection Guide

Tacta® Order Code	Channels	Volume Range	Increment	Pipettes and Tips	Optifit Tip	Safetyspace® Tip	Safe-Cone Filters	
		(µL)	(µL)	Color-Code	Volume	Volume	Standard	Plus
LH-729010	1	0.3 – 3	0.002	■	10 µL	10 µL	-	-
LH-729020	1	1 – 10	0.01	■	10 µL	10 µL	-	-
LH-729030	1	2 – 20	0.02	■	200 µL	20 µL	721014	-
LH-729050	1	10 – 100	0.10	■	200 µL	120 µL	721008	721018
LH-729060	1	20 – 200	0.20	■	200 µL	200 µL	721007	721017
LH-729070	1	100 – 1,000	1.00	■	1,000 µL	1,000 µL	721006	721016
LH-729080	1	500 – 5,000	5.0	■	5,000 µL	5,000 µL	721005	721015
LH-729090	1	1,000 – 10,000	10.0	■	10,000 µL	-	721005	721015
LH-729120	8	1 – 10	0.01	■	10 µL	10 µL	-	-
LH-729220	12							
LH-729130	8	10 – 100	0.10	■	200 µL	120 µL	721008	721018
LH-729230	12							
LH-729140	8	30 – 300	0.20	■	350 µL	300 µL	721007	721017
LH-729240	12							

For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at [sartorius.com](https://www.sartorius.com)



Tacta® Multipacks Smart and Complete Set of Pipettes and Accessories

Choose from a range of eight Multipacks that include a selection of Tacta® Mechanical Pipettes, a Linear Stand, matching Optifit Tips, and Pipette Holders.



Ordering Information

Tacta® Mechanical Pipettes Multipack

Multipack Order Code	Tacta® Pipettes	Optifit Tips	Accessories
Tacta® Pipette 3 – pack 10 LH – 729670	1-channel <ul style="list-style-type: none"> ▪ 1 – 10 µL ▪ 10 – 100 µL ▪ 100 – 1,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> ▪ 10 µL ▪ 200 µL ▪ 1,000 µL 	Linear Stand
Tacta® Pipette 3 – pack 20 LH – 729671	1-channel <ul style="list-style-type: none"> ▪ 2 – 20 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> ▪ 200 µL ▪ 1,000 µL 	Linear Stand
Tacta® Pipette 4 – pack 20 LH – 729672	1-channel <ul style="list-style-type: none"> ▪ 1 – 10 µL ▪ 2 – 20 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> ▪ 10 µL ▪ 200 µL ▪ 1,000 µL 	Linear Stand
Tacta® Pipette 4 – pack 100 LH – 729673	1-channel <ul style="list-style-type: none"> ▪ 1 – 10 µL ▪ 10 – 100 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> ▪ 10 µL ▪ 200 µL ▪ 1,000 µL 	Linear Stand
Tacta® Pipette 5 – pack 10 LH – 729674	1-channel <ul style="list-style-type: none"> ▪ 1 – 10 µL ▪ 10 – 100 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL ▪ 500 – 5,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> ▪ 10 µL ▪ 200 µL ▪ 1,000 µL Tip Rack (50 tips) <ul style="list-style-type: none"> ▪ 5,000 µL 	Linear Stand
Tacta® Pipette 5 – pack 20 LH – 729675	1-channel <ul style="list-style-type: none"> ▪ 2 – 20 µL ▪ 10 – 100 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL ▪ 500 – 5,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> ▪ 200 µL ▪ 1,000 µL Tip Rack (50 tips) <ul style="list-style-type: none"> ▪ 5,000 µL 	Linear Stand
Tacta® Pipette 3+1 – pack LH – 729676	1-channel <ul style="list-style-type: none"> ▪ 1 – 10 µL ▪ 20 – 200 µL ▪ 100 – 1,000 µL 8-channel <ul style="list-style-type: none"> ▪ 30 – 300 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> ▪ 10 µL ▪ 200 µL ▪ 350 µL ▪ 1,000 µL 	Linear Stand
Tacta® Pipette 3 – pack max LH – 729677	1-channel <ul style="list-style-type: none"> ▪ 100 – 1,000 µL ▪ 500 – 5,000 µL ▪ 1,000 – 10,000 µL 	Tip Racks <ul style="list-style-type: none"> ▪ 1,000 µL (96 tips) ▪ 5,000 µL (50 tips) ▪ 10,000 µL (35 tips) 	Linear Stand



Proline® Plus Mechanical Pipettes

Dependable Durability



Fixed volume Proline® Plus pipette



Proline® Plus has a comfortable handle and an ergonomic finger support for effortless pipetting.



Safe-Cone Filters protect the pipette from contamination, and should be changed regularly.

The Proline® Plus mechanical pipette family is designed to offer comfort and quality for your everyday manual pipetting. It combines durable construction with ease and lightness of use, and is therefore the perfect choice for liquid-handling professionals and students alike. In addition, it has the widest pipette range, including fixed volume pipettes for when volumes need to be ready-set to avoid errors.

Ergonomic Design

Proline® Plus has low pipetting forces, a comfortable handle and an ergonomic finger support for effortless pipetting. Good fit in hand minimizes the grip force needed to hold the pipette, thereby reducing the risk of strain injury.

Loading Tips on Multi-channel Pipettes with Minimum Force

Loading and ejecting tips with multi-channel pipettes requires relatively high forces. Proline® Plus multi-channel pipettes have spring loaded tip cones – the Optiload mechanism. Optiload allows tips to be loaded and ejected with minimum force, which reduces the risk of hand injuries. Moreover, it secures even tip sealing onto every individual tip cone.

Protect the Pipette with Safe-Cone Filters

The replaceable Safe-Cone Filter located inside the tip-cone prevents aerosols and fluids from penetrating the pipette, also in case of over-aspiration. The use of Safe-Cone Filters lengthens the maintenance interval of the pipette. Safe-Cone Filters are available for all Proline® Plus models greater than 10 µL. They must be replaced regularly, and always in case of overaspiration.

Easy Maintenance and Adjustment

No opening tools are needed for cleaning and maintaining Proline® Plus pipettes, and only three parts need to be cleaned. These pipettes are also easy to adjust using the adjustment tool that is provided with the pipette.

Features and Benefits



Proline® Plus pipettes have only three parts that need cleaning and maintenance.

- Low pipetting forces that prevent repetitive injury and improve results in long pipetting series
- Ergonomic finger support minimizes the grip force needed to hold the pipette
- Optiload mechanism in multi-channel models for easy and light tip loading with perfect tip sealing
- Wide variety of adjustable single and multi-channel models as well as fixed single-channel models
- Volume range from 3 μ L (5 μ L for fixed) to 10 mL
- Easy volume setting with click stop mechanism
- Big volume display
- Color-coding of different volumes to ease the selection of matching pipette tips
- Safe-Cone Filters available for models > 10 μ L
- Fully autoclavable without disassembly
- Simple to clean and maintain with only three parts to disassemble
- Materials have high chemical and UV-resistance to secure long life-time of the pipette



Ordering Information

Proline® Plus Variable Volume Mechanical Pipettes

Tip Selection Guide

Proline® Plus Order Code	Channels	Volume Range	Increment	Pipettes and Tips	Optifit Tip	Safetyspace® Tip	Safe-Cone Filters	
		(µL)	(µL)	Color-Code	Volume	Volume	Standard	Plus
728010	1	0.3 – 3	0.002	■	10 µL	10 µL	–	–
728020	1	1 – 10	0.01	■	10 µL	10 µL	–	–
728030	1	2 – 20	0.02	■	200 µL	20 µL	721014	–
728040	1	5 – 50	0.10	■	200 µL	120 µL	721008	721018
728050	1	10 – 100	0.10	■	200 µL	120 µL	721008	721018
728060	1	20 – 200	0.20	■	200 µL	200 µL	721007	721017
728070	1	100 – 1,000	1.00	■	1,000 µL	1,000 µL	721006	721016
728080	1	500 – 5,000	10.0	■	5,000 µL	5,000 µL	721005	721015
728090	1	1,000 – 10,000	20.0	■	10,000 µL	–	721005	721015
728120	8	1 – 10	0.01	■	10 µL	10 µL	–	–
728220	12							
728130	8	10 – 100	0.10	■	200 µL	120 µL	721008	721018
728230	12							
728140	8	30 – 300	0.20	■	350 µL	300 µL	721007	721017
728240	12							

For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at sartorius.com



Ordering Information

Tip Selection Guide

Proline® Plus FIXED Volume, Single Channel

Order Code	Channels	Volume	Pipettes and Tips	Optifit Tip	Safetyspace® Tip	Safe-Cone Filters	
		(µL)	Color-Code	Volume	Volume	Standard	Plus
728515	1	5	■	10 µL	10 µL	-	-
728520	1	10	■	10 µL	10 µL	-	-
728530	1	20	■	200 µL	20 µL	721014	-
728535	1	25	■	200 µL	120 µL	721008	721018
728545	1	50	■	200 µL	120 µL	721008	721018
728550	1	100	■	200 µL	120 µL	721008	721018
728560	1	200	■	200 µL	200 µL	721007	721017
728565	1	250	■	1,000 µL	1,000 µL	721006	721016
728567	1	500	■	1,000 µL	1,000 µL	721006	721016
728570	1	1,000	■	1,000 µL	1,000 µL	721006	721016
728575	1	2,000	■	5,000 µL	5,000 µL	721005	721015
728580	1	5,000	■	5,000 µL	5,000 µL	721005	721015
728590	1	10,000	■	10,000 µL		721005	721015

For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at [sartorius.com](#)



Proline® Plus Fixed Volume Pipettes



Proline® Plus Multipacks – Complete Sets of Pipettes and Accessories

Proline® Plus Multipacks offer are sets of mechanical pipettes in an affordable package, including a Linear Stand and racks of matching tips.

These Multipacks allow the affordable renewal of existing pipettes, or to set up new workstations, with a set of highly ergonomic pipettes with high accuracy and precision.

Ordering Information

Proline® Plus Mechanical Pipette Multipacks

Multipack Order Code	Proline® Plus Pipettes	Sartorius Optifit Tips	Accessories
Proline® Plus Pipette 3-pack 10 LH-728670	1-channel <ul style="list-style-type: none"> 1 – 10 µL 10 – 100 µL 100 – 1,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> 10 µL 200 µL 1,000 µL 	Linear Stand
Proline® Plus Pipette 3-pack 20 LH-728671	1-channel <ul style="list-style-type: none"> 2 – 20 µL 20 – 200 µL 100 – 1,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> 200 µL 1,000 µL 	Linear Stand
Proline® Plus Pipette 4-pack 20 LH-728672	1-channel <ul style="list-style-type: none"> 1 – 10 µL 2 – 20 µL 20 – 200 µL 100 – 1,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> 10 µL 200 µL 1,000 µL 	Linear Stand
Proline® Plus Pipette 4-pack 100 LH-728673	1-channel <ul style="list-style-type: none"> 1 – 10 µL 10 – 100 µL 20 – 200 µL 100 – 1,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> 10 µL 200 µL 1,000 µL 	Linear Stand
Proline® Plus Pipette 5-pack 10 LH-728674	1-channel <ul style="list-style-type: none"> 1 – 10 µL 10 – 100 µL 20 – 200 µL 100 – 1,000 µL 500 – 5,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> 10 µL 200 µL 1,000 µL Tip Rack (50 tips) <ul style="list-style-type: none"> 5,000 µL 	Linear Stand
Proline® Plus Pipette 5-pack 20 LH-728675	1-channel <ul style="list-style-type: none"> 2 – 20 µL 10 – 100 µL 20 – 200 µL 100 – 1,000 µL 500 – 5,000 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> 200 µL 1,000 µL Tip Rack (50 tips) <ul style="list-style-type: none"> 5,000 µL 	Linear Stand
Proline® Plus Pipette 3+1-pack LH-728676	1-channel <ul style="list-style-type: none"> 1 – 10 µL 20 – 200 µL 100 – 1,000 µL 8-channel <ul style="list-style-type: none"> 30 – 300 µL 	Tip Racks (96 tips) <ul style="list-style-type: none"> 10 µL 200 µL 350 µL 1,000 µL 	Linear Stand
Proline® Plus Pipette 3-pack max LH-728677	1-channel <ul style="list-style-type: none"> 100 – 1,000 µL 500 – 5,000 µL 1,000 – 10,000 µL 	Tip Racks <ul style="list-style-type: none"> 1,000 µL (96 tips) 5,000 µL (50 tips) 10,000 µL (35 tips) 	Linear Stand



Stands and Accessories

Table of Contents

Stands and Accessories

Pipette Stands

Safe-Cone Filters

Adjustment Tool

30

31

33



Pipette Stands



Charging Carousel Stand: 730991



Charging Stand: 730981



Linear Stand (non-charging):
725620



Carousel Stand (non-charging):
LH-725630

When the pipette is not in use, it should be stored in an upright position in order to avoid contamination from work surfaces. Sartorius provides stands for all of its pipettes. It is recommended that electronic pipettes be stored and charged on a charging stand whenever they are not in use. In this way, their batteries always remain charged for when work begins.

Compact carousel stands are ideal for saving bench space. There is one for mechanical pipettes, and a charging carousel stand for electronic pipettes.

The Linear Stand is designed for all Sartorius mechanical and electronic pipettes. This stand is also compatible with a wide range of other pipette makes.

The simplest of all are the pipette holders which are attached to the front edge of a shelf. These are suitable for mechanical pipettes.

Ordering Information

Pipette Stands

Order Code	Item
730981	Charging Stand for one electronic pipette*
730991	Charging Carousel for 4 electronic pipettes*
725620	Linear Stand for all Sartorius pipette models
LH-725630	Carousel Stand for 6 mechanical pipettes
LH-727650	Adapter for Mechanical Carousel Stand
LH-727640	Holder for one pipette

* Supplied with a universal charger (EU, UK, US | JPN, AUS, KOR and CHN plugs)



Adapter for Mechanical
Carousel Stand: LH-727650



Pipette Holder for one Pipette:
LH-727640

Safe-Cone Filters



Built-in filter ejector in Tacta®

Why Should You Use Safe-Cone Filters?

These unique and replaceable polyethylene (PE) filters prevent any fluids and liquid vapours from reaching the internal components of the pipette.

- Reduce the risk of contaminating the internal components of our pipettes
- Prolong the pipette's lifetime
- Reduce maintenance intervals
- Are cost-effective compared to filter tips
- We do not recommend to use Safe Cone Filters in combination with filter tips.

When Should You Use them?

The ultimate pipette protectors are available in two types:

Plus Filter

For more demanding applications such as radioactive work, cell culture, bacterial and virological work and molecular biology.

Standard Filter


For general applications. Can be used in same type of work as the Plus filter, but needs to be changed more frequently.

How Often Should You Change?

The interval of changing the filter depends completely on the application and the sample. However, according to studies, the filter is recommended to be changed daily (after 50 to 250 pipettings) and immediately in case of over-aspiration.


How to Change?

To ensure that the user is protected from contamination, tweezers should be used when removing used filters from the pipette tip cone. The Tacta® feature a built-in filter ejectors, so tweezers are unnecessary. In addition, the tip cone should be cleaned with ethanol (70%) prior to the insertion of a new filter.



Tweezers for exchanging filters in pipettes are supplied with all pipettes excluding Tacta®.

Contamination in Pipette Interior



Pipette contamination in pipette interior when pipetting liquid culture of bacteria *Micrococcus Luteus*.

Ordering Information

Safe-Cone Filters

Order Code	Item	Qty / Unit
721008	Standard Ø 2.51 mm PE	50
721007	Standard Ø 3.15 mm PE	50
721006	Standard Ø 5.33 mm PE	50
721005	Standard Ø 6.73 mm PE	50
721014	Standard Ø 1.83 mm Polyolefin	50
721018	Plus Ø 2.51 mm PE	50
721017	Plus Ø 3.15 mm PE	50
721016	Plus Ø 5.33 mm PE	50
721015	Plus Ø 6.73 mm PE	50
721009	Tweezers for Safe-Cone Filters	1

PE = polyethylene

See the pipette ordering information charts for corresponding filters and pipettes.

Adjustment Tool



Adjustment Tool is used for adjusting Proline® Plus pipettes.



The Adjustment Tool for adjusting the Tacta® pipette.

The Adjustment Tool is used to adjust the pipette in situations where the factory calibration is not applicable.

Ordering Information

Adjustment Tool

Order Code	Item	Qty
726203	Adjustment Tool for Proline® Plus	1
LH-727080	Adjustment Tool for Tacta®	1





Pipette Tips

Table of Contents

Pipette Tips

Pipette Tips

36

Optifit Tips

38

Safetyspace® Filter Tips

39

Packaging Options

40



Pipette Tips

The Perfect Match for Your Pipette



Optiload for a tight fit and equal sealing on every channel's tip cone

Sartorius pipette tips ensure the performance of Sartorius pipettes and repeatability of your results like no other tip can. They have been designed and manufactured to meet the highest quality and purity standards and to be the inert protectors of your samples. Moreover, correctly fitting tips protect the pipette's tip cone from wear and tear.

Sartorius tip packages are designed to make the daily work of lab professionals easier. Our offering covers a variety of functional tip package options with various purity ratings.

The high purity and consistent quality of Sartorius tips provide your valuable samples with the ultimate protection from contamination. We adhere to strict quality standards and control procedures – from raw material to automated manufacturing and packaging.

Sartorius tips correspond with the color-coding of Sartorius pipettes, to allow easy matching of corresponding volumes.



Premium Quality and Purity

Manufacturing the tips in our own production facility allows us to maintain the highest quality and purity standards, by selecting the best plastic materials and controlling the manufacturing process from beginning to end.

Our quality management system follows not only ISO 9001 and ISO 14001, but also ISO 13485. Tip production also abides by the ISO 14644-1 standard, in order to fulfil ISO Class 8 cleanroom conditions for secured tip purity.

Contamination Free Tips

To avoid contamination through human contact, we have automated the entire tip manufacturing process. Pure virgin polypropylene (PP) plastic is automatically fed from silos into molding machines. Molding machines and robots located in isolated clean cells, load the tips automatically into tip trays and packaging. HEPA filters and higher air pressure are applied for purity within the cell. Sartorius Single Tray tip racks, Single Refill Packs and FlexiBulk® packs are individually and automatically packed in air-tight plastic, in order to rule out any danger of contamination.

Additionally, our highly experienced and trained personnel are equipped with specially designed coveralls, masks, hair nets and gloves, in order to further diminish risks of contamination.

An independent laboratory checks every tip lot of Single Tray and Refill Pack for DNase, RNase, human DNA and endotoxins. Lot-specific purity certificates can be downloaded from sartorius.com.

SARTORIUS	
CERTIFICATE OF RELEASE	
Product Description:	
Product name:	Soft Tip 10-1000 µl
Reference:	281000-001
Batch number:	281000-001
Supply date:	2024-10-17
Quantity:	1000000 pieces
Specifications:	
This product is certified for DNase, RNase, endotoxin and human DNA free and is tested by an independent laboratory.	
Examination:	
Sample preparation was carried out according to ICH guidelines for medical devices by Sartorius. The sample was tested for DNase, RNase, endotoxin and human DNA according to the following methods: DNase: DNase activity was determined by a colorimetric assay. RNase: RNase activity was determined by a colorimetric assay. Endotoxin: Endotoxin activity was determined by a colorimetric assay. Human DNA: Human DNA activity was determined by a colorimetric assay.	
Release:	
The lot was released up and down in an optimal volume of pure water setting meeting the requirements of the product specification. The lot was released up and down in an optimal volume of pure water setting meeting the requirements of the product specification. The lot was released up and down in an optimal volume of pure water setting meeting the requirements of the product specification. The lot was released up and down in an optimal volume of pure water setting meeting the requirements of the product specification.	
Human DNA:	
The lot was released up and down in an optimal volume of pure water setting meeting the requirements of the product specification. The lot was released up and down in an optimal volume of pure water setting meeting the requirements of the product specification. The lot was released up and down in an optimal volume of pure water setting meeting the requirements of the product specification. The lot was released up and down in an optimal volume of pure water setting meeting the requirements of the product specification.	
This is to certify that the product described has been manufactured in accordance with ISO 9001:2015. Certificate No. 0111001 and ISO 14001:2015. Certificate No. 0111002. Quality management processes and the product of Sartorius. All release procedures related to the product.	
The expiry date of this product is 5 years from the manufacturing date and is stated on the product label of the packaging.	
This is an electronic document and valid without signature.	
Signature:	Head Office, Sartorius 1
Signature:	Head Office, Sartorius 2
Signature:	Head Office, Sartorius 3
Signature:	Head Office, Sartorius 4
Signature:	Head Office, Sartorius 5
Signature:	Head Office, Sartorius 6
Signature:	Head Office, Sartorius 7
Signature:	Head Office, Sartorius 8
Signature:	Head Office, Sartorius 9
Signature:	Head Office, Sartorius 10
Signature:	Head Office, Sartorius 11
Signature:	Head Office, Sartorius 12
Signature:	Head Office, Sartorius 13
Signature:	Head Office, Sartorius 14
Signature:	Head Office, Sartorius 15
Signature:	Head Office, Sartorius 16
Signature:	Head Office, Sartorius 17
Signature:	Head Office, Sartorius 18
Signature:	Head Office, Sartorius 19
Signature:	Head Office, Sartorius 20

Lot-specific purity certificate



Definitions	
Human DNA	Sartorius purity certified tips are analyzed for the presence of DNA using quantitative PCR and human DNA specific primers. The determined detection limit of the assay is < 1 pg/μL.
DNase	A deoxyribonuclease (DNase) is any enzyme that catalyzes the degradation of DNA. The absence of DNase is tested by using fluorometric assay. The detection level of the assay is 6.25×10^{-5} U/μL when DNase I is used as a standard.
RNase	Ribonuclease (RNase) is an enzyme that catalyzes the degradation of RNA into smaller components and can be generally found from organisms. The absence of RNase is tested by using fluorometric assay. The detection level of the assay is $< 1 \times 10^{-8}$ U/μL.
Endotoxins	Endotoxins are lipopolysaccharides found in the Gram-negative bacteria and can cause several serious health effects in humans and animals. Limulus Amebocyte Lysate (LAL) Gel Clot method is used to detect the presence of endotoxins on the pipette tips. The detection level of the LAL assay is 0.005 IU/mL (EU/mL).*
Sterilization	The destruction of all microbial life, including bacterial endospores. Can be accomplished e.g. using steam, heating, chemicals, or radiation. We use e-beam irradiation.

* Speciality tip sterility assurance level, DNase, RNase, human DNA and endotoxin testing detection limits may differ from values reported here. Please download a batch specific certificate from sartorius.com

Features and Benefits

Best Fit – Highest Possible Accuracy

- Perfect fitting and sealing with Sartorius pipettes secure the highest possible accuracy and precision
- Compatible with Optiload feature in Sartorius Picus® 2, Tacta® and Proline® Plus pipettes enabling ergonomic and light tip attachment and ejection
- Color-coding of tip trays allows easy matching with a corresponding color-coded Sartorius pipette
- Universal design optimized for Sartorius pipettes

Premium Quality and Purity:

- Strict quality standards, ISO 9001 and ISO 13485, are followed from R&D to production and delivery
- Manufactured in ISO 8 classified clean room conditions
- Manufacturing process free of DNase, RNase, human DNA and endotoxins: Single Trays, Refill Packs and FlexiBulk® packs certified pure by lot number
- Pre-sterilized tips are e-beam irradiated
- All tip packages, including individual racks, are lot numbered for full traceability
- The highest quality virgin polypropylene used as raw material

Tip Selection Guide per Application

Tip Type	Optifit Tips			Safetyspace® Filter Tips
Purity	Standard	Free of DNase, RNase, human DNA & endotoxins	Pre-sterilized & free of DNase, RNase, human DNA & endotoxins	Pre-sterilized & free of DNase, RNase, human DNA & endotoxins
Regular pipetting applications	▪			
Applications where prevention of cross-contamination is vital		▪	▪	▪

Optifit Tips

Standard Tips for Various Needs




Sartorius Optifit tips are an excellent choice for various laboratories and pipetting tasks with their wide range of packaging and purity options. The Optifit tips are packed in single tray racks, refill towers, single refill packs, and bulk packages. Optifit tips are available DNase, RNase, human DNA and endotoxin free, as well as e-beam pre-sterilized.

The Single Tray tip racks are ideal for easy tip loading and for contamination-free pipetting. In order to reuse the empty tip racks and to create less waste, you may fill the empty racks with Refill tips, either using the refill tower or purity certified refill pack tips. The FlexiBulk® pack is the choice, if you need a cost-effective, yet purity-certified solution, in bulk tip format.



Available Packaging Options

- Single Tray Rack
- Refill Tower
- Refill Pack
- FlexiBulk® Pack
- Bulk in a Bag



Features and Benefits

- Standard non-filter tips made to Sartorius quality standards
- Perfect fitting and sealing with Sartorius Picus® 2, Tacta® and Proline® Plus pipettes
- Wide tip volume range from 10 µL to 10 mL
- Wide selection of packaging and purity options
- Available as DNase, RNase, human DNA and endotoxin-free
- e-beam pre-sterilized packaging options available
- Full traceability
- Color-coded trays to match with corresponding Sartorius pipettes
- Fully autoclavable at 121 °C, at 1 bar, for 20 minutes

Optifit tips provide excellent tip sealing

Safetyspace® Filter Tips

Protect Valuable Samples



Safetyspace® Filter Tips offer optimal protection for your work from cross contamination and enable you to use the full volume of the tip with any pipetting mode. They feature the unique Safetyspace® air gap that leaves additional space between the sample and the filter that conventional filter tips do not have. This extra space prevents the liquid from touching, and permeating, the filter and thus guarantees the pipetting accuracy.

Safetyspace® Filter tips are made of virgin polypropylene and are fitted with pure polyethylene filter barriers that effectively capture solid and liquid aerosol particles. The filter is made without self-sealing additives to avoid any interference with the sample or results. In addition to protecting your samples, the filter also protects your pipette.

The Safetyspace® Filter Tips are Ideal for:

- Molecular biology
- Microbiology
- Cell culture applications
- Radioactive work

The Extra Space is Particularly Useful in the Following Applications:

- Pipetting foaming liquids such as buffers and proteins
- Pipetting solvents
- Multiple dispensing functions of electronic pipettes
- Reverse pipetting

Available Packaging Options

- Single Tray Rack



Features and Benefits

- Filter minimizes the risk of aerosol contamination
- Safetyspace® air gap between the sample and filter reduces the risk of contaminating the internal components of our pipettes
- Enables the full use of pipettes' volume range with any pipetting mode
- Available in volumes from 10 µL to 5,000 µL
- Certified free of DNase, RNase, human DNA and endotoxins
- e-beam pre-sterilized
- Full traceability
- Color-coded trays indicate the matching color-coded Sartorius pipette

Packaging Options

Racked Tips

Single Tray Rack

- 96 tips in convenient and reusable tray racks (sales unit contains 10 tray racks, total 960 tips)
- Certified free of DNase, RNase, human DNA and endotoxins
- e-beam pre-sterilized option available
- Lot-specific purity certificates can be downloaded from sartorius.com
- Informative rack labelling: volume, product number, lot number improves tip identification and traceability
- Air-tight plastic wrapping around the rack secures purity during transport and storage (wrapping is regular waste)
- Tip trays are color-coded to indicate the matching, color-coded Sartorius pipette
- Covers a large range of tip volumes from 10 μ L to 10 mL
- Fully autoclavable at 121°C, at 1 bar, for 20 minutes
- Tray racks can be easily reloaded with Refill tips
- Racks, trays and tips are 100 % recyclable polypropylene (PP)



Single Tray Racks

Bulk Tips

FlexiBulk®

- Tips made to the Sartorius quality standard in economical packaging
- Packed orderly in compact re-sealable plastic packages (480 or 960 pcs depending on tip volume)
- Covers a large range of tip volumes from 200 μ L to 1,200 μ L
- Lot-specific purity certificates can be downloaded from sartorius.com
- Tips are fully autoclavable at 121°C, at 1 bar, for 20 minutes
- 100 % recyclable tips (PP) and package (PET)



FlexiBulk®

Bulk in a Bag

- Tips made to the Sartorius quality standard in economical packaging
- Packed in re-sealable bags in cardboard boxes (100, 250 or 1,000 pcs depending on tip volume)
- Covers tip volumes 10 μ L, 5 mL and 10 mL
- Tips are fully autoclavable at 121°C, at 1 bar, for 20 minutes
- 100 % recyclable tips and package



Bulk in a Bag



Single Refill Packs

Refill Tips

Single Refill Packs

- Individually packed air-tight tip trays for maximum purity with less packaging material compared to racked tips
- 10, 15 or 20 trays of 96 tip trays, depending on tip volume
- Certified free of DNase, RNase and endotoxins
- e-beam pre-sterilized option available
- Lot-specific purity certificates can be downloaded from [sartorius.com](https://www.sartorius.com)
- Tip trays are compatible with Single Tray racks for convenient use
- Informative rack labelling: volume, product number, lot number improves tip identification and traceability
- Tip trays are colour-coded to indicate the matching, colour-coded Sartorius pipette
- Covers a large range of tip volumes from 10 μ L to 1,200 μ L
- Trays and tips are fully autoclavable at 121°C, at 1 bar, for 20 minutes
- 100 % recyclable trays and tips (PP). Container lid is regular waste.



Refill Towers

Refill Tower

- Space-saving with 10 × 96 tips in one tower
- Tip trays are compatible with Single Tray racks for convenient use
- Trays are color-coded to indicate the matching, color-coded Sartorius pipette
- Covers the most widely used tip sizes: 10 μ L, 200 μ L and 350 μ L
- Trays and tips are fully autoclavable at 121°C, at 1 bar, for 20 minutes
- 100 % recyclable cardboard packaging, and plastic (PP) trays and tips



The ACT Ecolabel Commitment

With the ACT Environmental Impact Factor Label, our Refill Tower product has earned a third-party verification of its sustainable impacts.





The ACT Label (Accountability, Consistency and Transparency) is an eco-label initiative by My Green Lab, a non-profit organization dedicated to advancing sustainability across the life sciences, industrial sciences, and healthcare sectors.

Read more about pipette sustainability here:





<https://www.sartorius.com/en/products/pipetting/pipetting-and-sustainability>

Ordering Information

Optifit Tips

Volume Range	Length	Packaging	Purity Level		Tips/Unit	Order Code
			Free of DNase, RNase, human DNA & endotoxins	Pre-sterilized		
<div><div></div><div>10 µL</div></div>	31.5 mm	Single Tray	■		10×96	790010
		Single Tray	■	■	10×96	790011
		Refill Tower			10×96	790012
		Refill Pack	■	■	20×96	790013
		Bulk in Bag			1,000	790014
<div><div></div><div>200 µL</div></div>	51 mm	Single Tray	■		10×96	790200
		Single Tray	■	■	10×96	790201
		Refill Tower			10×96	790202
		Refill Pack	■	■	15×96	790203
		FlexiBulk®	■		960	LH-B790204
<div><div></div><div>350 µL</div></div>	54 mm	Single Tray	■		10×96	790350
		Single Tray	■	■	10×96	790351
		Refill Tower			10×96	790352
		Refill Pack	■	■	15×96	790353
		FlexiBulk®	■		960	LH-B790354
<div><div></div><div>1,000 µL</div></div>	71.5 mm	Single Tray	■		10×96	791000
		Single Tray	■	■	10×96	791001
		Refill Pack	■		10×96	791002
		Refill Pack	■	■	10×96	791003
		FlexiBulk®	■		480	LH-B791004

For your guidance the tips are shown here in the actual size.

Volume Range	Length	Packaging	Purity Level		Tips/Unit	Order Code
			Free of DNase, RNase, human DNA & endotoxins	Pre-sterilized		
<div><div></div>1,000µL Wide bore tip</div>	68.5 mm	Single Tray	<div><div></div></div>		10×96	791020
		Single Tray	<div><div></div></div>	<div><div></div></div>	10×96	791021
		FlexiBulk®	<div><div></div></div>		480	LH-B791024
						
<div><div></div>1,200µL</div>	71.5 mm	Single Tray	<div><div></div></div>		10×96	791200
		Single Tray	<div><div></div></div>	<div><div></div></div>	10×96	791201
		Refill Pack	<div><div></div></div>		10×96	791202
		Refill Pack	<div><div></div></div>	<div><div></div></div>	10×96	791203
		FlexiBulk®	<div><div></div></div>		480	LH-B791204
						
<div><div></div>5,000µL</div>	150 mm	Single Tray	<div><div></div></div>		50	780304
		Single Tray	<div><div></div></div>	<div><div></div></div>	50	780305
		Bulk in Bag			100	780300
		Bulk in Bag			1,000	780308
						
<div><div></div>10,000µL</div>	164 mm	Single Tray	<div><div></div></div>		4 x 35	LH-799000
		Bulk in Bag			210	LH-799004
						

Note: The ordering information for 10,000 µL tip for Midi Plus can be found on page 49.






Ordering Information

Safetyspace® Filter Tips




Volume Range	Length	Packaging	Purity Level		Tips/Unit	Order Code
			Free of DNase, RNase, human DNA & endotoxins	Pre-sterilized		
<div><div></div>10 µL</div> <div></div>	31.5 mm	Single Tray	<div></div>	<div></div>	10 × 96	790011F
<div><div></div>20 µL</div> <div></div>	51 mm	Single Tray	<div></div>	<div></div>	10 × 96	790021F
<div><div></div>120 µL</div> <div></div>	51 mm	Single Tray	<div></div>	<div></div>	10 × 96	790101F
<div><div></div>200 µL</div> <div></div>	52.5 mm	Single Tray	<div></div>	<div></div>	10 × 96	790201F
<div><div></div>300 µL</div> <div></div>	52.5 mm	Single Tray	<div></div>	<div></div>	10 × 96	790301F
<div><div></div>1,000 µL</div> <div></div>	78 mm	Single Tray	<div></div>	<div></div>	10 × 96	791001F
<div><div></div>1,200 µL</div> <div></div>	90 mm	Single Tray	<div></div>	<div></div>	10 × 96	791211F
<div><div></div>5,000 µL</div> <div></div>	150 mm	Single Tray	<div></div>	<div></div>	50	LH-795001F

For your guidance the tips are shown here in the actual size.
Filter tips are not recommended to be used simultaneously with Safe-Cone Filters.

Extended Standard Tips

Volume Range	Length	Packaging	Purity Level		Tips/Unit	Order Code
			Free of DNase, RNase, human DNA & endotoxins	Pre-sterilized		
 10 µL 	46 mm	Single Tray	■		10 × 96	783210
		Single Tray	■	■	10 × 96	783211
	77.5 mm	Single Tray	■		10 × 96	LH-X780200
		Single Tray	■	■	10 × 96	LH-X780201
	119 mm	Single Tray	■		10 × 96	LH-X791000
		Single Tray	■	■	10 × 96	LH-X791001
	90 mm	Single Tray	■		10 × 96	791210
		Single Tray	■	■	10 × 96	791211
		Refill Pack	■		10 × 96	791212
		Refill Pack	■	■	10 × 96	791213

Extended Filter Tips

Volume Range	Length	Packaging	Purity Level		Tips/Unit	Order Code
			Free of DNase, RNase, human DNA & endotoxins	Pre-sterilized		
	46 mm	Single Tray	■	■	10 × 96	783201
	77.5 mm	Single Tray	■	■	10 × 96	LH-XF780201
	119 mm	Single Tray	■	■	10 × 96	LH-XF791001

For your guidance the tips are shown here in the actual size.

Extended filter tips are not recommended to be used simultaneously with Safe-Cone Filters.

The liquid handling properties of extended tips might differ from standard Optifit and Safetyspace® Filter tips.



Maxi-Volume Liquid Handling

Table of Contents

Maxi-volume Liquid Handling

Midi Plus Pipetting Controller

Prospenser Plus and Prospenser

Biotrate Digital Burette

48

50

52



Midi Plus Pipetting Controller

Excellent Performance and Ergonomics



The Midi Plus is a lightweight electronic cordless pipetting controller, which allows aspiration from bottles and tubes, without the arm and hand elevations required in the case of serological or volumetric pipettes.

It fits all commonly used 1–100 mL glass or plastic pasteur pipettes, but can also be used with Sartorius 5 mL and 10 mL disposable tips. The speed can be fine-tuned by applying varying finger pressure to the operating buttons.

The Midi Plus is ideal, for example, in microbiological work: dispensing into a culture media dish can be performed carefully, drop by drop, without breaking the fine surface of the medium.


- Stepless speed control
- Hydrophobic autoclavable filter protects the device in case of over-aspiration
- Convenient fold-out bench stand supports the unit and pipette when not in use
- Low battery warning

Features

Pipette types	Plastic or glass 1 – 100 mL Pasteur pipettes 5 mL Optifit and 10 mL Midi Plus Sartorius pipette tips
Rechargeable during use	Yes
Speed control	Stepless adjustable control
Gravity dispensing	Yes
Stand	Attached support
Weight	207 g
Low battery indicator	Yes
Autoclavable parts	Nose cone, pipette holder and filter

Ordering Information

Midi Plus & Accessories and Consumables



Order Code	Item	Qty / Unit
710931	Midi Plus Pipetting Controller with Universal Charger ¹	1
LH-7129120	Replacement filter, 0.45 µm, non-sterile, autoclavable	5
LH-7129131	Replacement filter, 0.2 µm, sterile	5
LH-711019	Adapter set (standard), autoclavable, including nose-cone and silicone adapter (included in 710931)	1
LH-711017	Adapter set for 5 mL pipette tip, autoclavable, including nose-cone and silicone adapter	1
LH-711018	Adapter set for 10 mL pipette tip, autoclavable, including nose-cone and silicone adapter	1
780300	Optifit Tip 5 mL (length 150 mm)	100
780308	Optifit Tip 5 mL (length 150 mm)	1,000
780310	Midi Plus Tip 10 mL (length 242 mm)	100

¹ Supplied with a universal charger (incl. EU, UK, US | JPN, AUS, KOR and CHN plugs)



Prospenser Plus and Prospenser

Dispensing Made Easy

Chemical Resistance

Prospenser and Prospenser Plus offer excellent chemical resistance thanks to their high quality materials. The ceramic piston is chemically stable and compatible with an extremely wide range of liquids.

Light to Use

The operating handle is designed to be comfortable to hold while aspirating and light to press while dispensing. The piston moves smoothly, guaranteeing excellent results in everyday use.

Easily Replaceable Valves

Valves can be easily replaced by the user so there is no need to send the unit to be serviced.

Easy to Attach Aspiration Tube

The aspiration tube is easy to attach with a secure threaded connection. This prevents air entering the system and ensures consistent results.

Volume Setting

Prospenser and Prospenser Plus have an incremental volume setting enabling you to set the exact volume every time.

Versatile Accessories

There is a wide range of accessories available including bottle adapters, dispensing heads, and aspiration tubes.

Fully Autoclavable

Prospenser and Prospenser Plus are easy to disassemble for cleaning and are fully autoclavable.



Prospenser Plus – Advanced Features



Media Recirculation System

Prospenser Plus is equipped with a media recirculation system that allows unused media to be recovered from the device and returned to the bottle. This ensures easy priming and minimizes media loss.

Smart 360° Rotating Dispensing Head

The dispensing head can be turned to the required position without needing to turn the bottle. This ensures the bottle label is always visible to the user.

Connector for Drying Tube and Filters

Prospenser Plus has a luer connection at the back of the unit for drying tubes and filters. Drying tubes are recommended for humidity sensitive liquids.



Biotrate Digital Burette for Convenient Titration

Sartorius Biotrate is a premium digital burette equipped with sophisticated functions for easy titration. Biotrate offers high chemical resistance and it is a smart solution for titrating various liquids accurately and safely. Biotrate is easy and convenient to use due to its large and clear electronic display, smooth operating wheels, and 360° rotating dispensing head.

High Chemical Resistance

- Highly resistant parts in the liquid pathway guarantee excellent chemical resistance.
- The softly-moving piston ensures reliable results.

Premium Quality

- High quality parts and materials guarantee highly accurate and precise results as well as long lifetime.

Convenient to Use

- Biotrate is comfortable to use due to soft touch operating wheels and 360° rotating dispensing head.
- Extremely smooth and light to operate.
- Biotrate offers excellent chemical resistance due to its high-quality materials, allowing the use of wide variety of liquids.

Light to Use

- The operating wheels are light and comfortable to use due to the soft touch surface and smoothly moving piston. These guarantee good ergonomics and excellent results in everyday use.

Long Life Battery

- Biotrate is powered with a long life battery so there is no need to plug in or re-charge the unit. The battery is also easy to replace by the user.

Large and Clear Display

- The display indicates clearly the measured volume, and shows whether you are measuring an aspirated or titrated volume.

Easy to Clean

- Biotrate is easy to disassemble for cleaning and autoclaving. The liquid pathway is fully autoclavable.

User Adjustment

- Biotrate can be adjusted for various liquids and it is easy to revert to factory settings. The display also indicates when a custom adjustment is in use.



Ordering Information

Prospenser and Prospenser Plus

Order Code	Item	Volume Range	Increment
LH-723070	Prospenser Plus	0.2 – 1 mL	0.05 mL
LH-723071	Prospenser Plus	0.4 – 2 mL	0.05 mL
LH-723072	Prospenser Plus	1 – 5 mL	0.10 mL
LH-723073	Prospenser Plus	2 – 10 mL	0.25 mL
LH-723074	Prospenser Plus	5 – 30 mL	0.50 mL
LH-723075	Prospenser Plus	10 – 60 mL	1.00 mL

Note: All Prospenser Plus models have A45 thread as a default and they are supplied with A32, A38, and S40 bottle adapters. Bottle not included
For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at sartorius.com

Prospenser

Order Code	Item	Volume Range	Increment
LH-723060	Prospenser	0.2 – 1 mL	0.05 mL
LH-723061	Prospenser	0.4 – 2 mL	0.05 mL
LH-723062	Prospenser	1 – 5 mL	0.10 mL
LH-723063	Prospenser	2 – 10 mL	0.25 mL
LH-723064	Prospenser	5 – 30 mL	0.50 mL
LH-723065	Prospenser	10 – 60 mL	1.00 mL


Note: All Prospenser models have A32 thread as a default and they are supplied with A28, S40, and A45 bottle adapters. Bottle not included
For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at sartorius.com

Biotrate Digital Burette

Order Code	Item	Volume Range	Increment
LH-723080	Biotrate	10 mL	0.01 mL
LH-723081	Biotrate	20 mL	0.01 mL
LH-723082	Biotrate	50 mL	0.01 mL

Note: All Biotrate models have A45 thread as a default and they are supplied with A32, A38, and 40 bottle adapters. Bottle not included
For pipette specifications (incl. Systematic and Random error limits), please see product datasheet, available at sartorius.com

Accessories




UV-Protection Window
Biotrate is supplied with an amber colored replacement window for protecting light sensitive media.



Bottle support, for bottles
75 – 120 mm



Aspiration tube, 310 mm



Dispensing head, 1.5 mm spiral

Biotrate

Order Code	Item	Material	Biotrate
LH-721650	Dispensing head, standard	FEP	▪
LH-721651	Dispensing head, fine tip	FEP	▪
LH-721648	Dispensing head, 1.5 mm spiral	FEP	▪
LH-721652	Dispensing head, 0.4 m adjustable	FEP	▪
LH-721653	Dispensing head, luer lock	FEP PP	▪
LH-721654	Dispensing head, luer lock	FEP PFA	▪
LH-721678	Aspiration tube, 310 mm	FEP	▪
LH-721679	Aspiration tube, 400 mm	FEP	▪
LH-721682	Recirculation tube, 70 mm	FEP	▪
LH-721689	Bottle adapter, A28	ETFE	▪
LH-721684	Bottle adapter, A28	PP	▪
LH-721683	Bottle adapter, A32	PP	▪
LH-721688	Bottle adapter, A32	ETFE	▪
LH-721686	Bottle adapter, A38	PP	▪
LH-721733	Bottle adapter, A38	ETFE	▪
LH-721685	Bottle adapter, A38 430	PP	▪
LH-721732	Bottle adapter, A38 430	ETFE	▪
LH-721687	Bottle adapter, S40	PP	▪
LH-721734	Bottle adapter, S40	ETFE	▪
LH-721735	Bottle adapter, NS29 32	Silicone	▪
LH-721743	Bottle support, for bottles 75 – 120 mm	Silicone	▪
LH-721744	UV-Protection Window		▪

Prospenser & Prospenser Plus

Order Code	Item	Material	Prospenser Plus	Prospenser	Prospenser
			All models	1, 2, 5 & 10 mL	30 & 60 mL
LH-721647	Dispensing head, standard	FEP	■		
LH-721655	Dispensing head, standard	FEP		■	
LH-721656	Dispensing head, standard	FEP			■
LH-721648	Dispensing head, 1.5 mm spiral	FEP	■		
LH-721649	Dispensing head, 3.0 mm spiral	FEP	■		
LH-721657	Dispensing head, spiral	FEP		■	
LH-721658	Dispensing head, spiral	FEP			■
LH-721653	Dispensing head, luer lock	FEP PP	■		
LH-721654	Dispensing head, luer lock	FEP PFA	■		
LH-721659	Dispensing head, luer lock	FEP PP			■
LH-721660	Dispensing head, luer lock	FEP PFA			■
LH-721678	Aspiration tube, 310 mm	FEP	■		■
LH-721680	Aspiration tube, 310 mm	FEP		■	
LH-721679	Aspiration tube, 400 mm	FEP	■		■
LH-721681	Aspiration tube, 400 mm	FEP		■	
LH-721682	Recirculation tube, 70 mm	FEP	■		
LH-721736	Bottle adapter, A25	PP		■	■
LH-721737	Bottle adapter, A28	PP		■	■
LH-721689	Bottle adapter, A28	ETFE	■		
LH-721684	Bottle adapter, A28	PP	■		
LH-721683	Bottle adapter, A32	PP	■		
LH-721688	Bottle adapter, A32	ETFE	■		
LH-721686	Bottle adapter, A38	PP	■		
LH-721733	Bottle adapter, A38	ETFE	■		
LH-721738	Bottle adapter, A38	PP		■	■
LH-721685	Bottle adapter, A38 430	PP	■		
LH-721732	Bottle adapter, A38 430	ETFE	■		
LH-721741	Bottle adapter, A45	ETFE		■	■
LH-721739	Bottle adapter, A45	PP		■	■
LH-721687	Bottle adapter, S40	PP	■		
LH-721734	Bottle adapter, S40	ETFE	■		
LH-721740	Bottle adapter, S40	PP		■	■
LH-721742	Bottle adapter, S40	ETFE		■	■
LH-721735	Bottle adapter, NS29 32	Silicone	■		
LH-721743	Bottle support, for bottles 75 – 120 mm	Silicone	■	■	■



Pipetting Academy

Table of Contents

Pipetting Academy

Pipetting Academy


58

Pipetting Recommendations

60

Pipetting Academy

Training for Better Performance, Ergonomics and Safety



Sartorius Pipetting Academy offers insight and training for pipetting related topics. It's a comprehensive combination of theory and practice that provides you with useful knowledge and tips for your daily work. Thousands of laboratory professionals around the globe have been trained at the Pipetting Academy over the last two decades.

Pipetting Academy Training Modules

We at Sartorius understand that a laboratory has professionals with different backgrounds and experience levels, and therefore we have created several different ready-made modules.

Pick and choose which modules are relevant for you and your laboratory or discuss with your Sartorius representative to tailor make or combine several different modules.

From each module a certificate can be provided upon completion for internal documentation.



Module	Description	Learning Objectives
1. Basics of Pipetting Module	Basics module brings together the essentials of pipetting techniques, ergonomics, and calibration. It is an excellent package for anyone who is new to pipetting or feels a need to polish pipetting practices.	<ul style="list-style-type: none"> ▪ To choose the right equipment ▪ The basics of pipetting techniques ▪ The best practices in lab ergonomics ▪ Daily maintenance and checking of pipettes
2. Ergonomics Module	The ergonomics module is intended for laboratory personnel and occupational and health professionals who are interested in improving the work ergonomics, minimizing ergonomic risk factors and preventing musculoskeletal disorders.	<ul style="list-style-type: none"> ▪ To recognize the ergonomic risk factors in pipetting and lab work ▪ To select the right tools and working postures ▪ To apply the best ergonomics and make your work more enjoyable
3. Pipetting Techniques Module	Pipetting is a precision task influenced by multiple factors such as equipment, pipetting techniques, liquids, and environmental conditions.	<ul style="list-style-type: none"> ▪ To understand the influence of pipetting techniques on pipetting results ▪ To identify the influence of different error sources on pipetting results ▪ How to select pipetting tools and techniques for different liquids and conditions
4. Pipetting in Cell Culture Applications Module	Working in cell culture lab has special requirements in terms of techniques and the purity of equipment.	<ul style="list-style-type: none"> ▪ How to avoid contamination ▪ The best pipetting practices for cell culture work ▪ To recognize the factors affecting your results
5. Pipetting Standards and Maintenance Module	Regular pipette maintenance and calibration ensures reliable results and conformance to regulation.	<ul style="list-style-type: none"> ▪ To define your pipette maintenance program according to ISO 8655 standard ▪ To handle pipettes as precision instruments ▪ To ensure reliable performance of your pipettes
6. Best Practices for Working with Proteins Module	To achieve the best results, protein samples should have a high concentration and maintain a high biological activity.	<ul style="list-style-type: none"> ▪ Best pipetting practices for protein work ▪ How to ensure your sample quality and yield
7. Best Pipetting Practices for Nucleic Acid Work Module	DNA and RNA extraction and qPCR assays are the cornerstones of many research in life science.	<ul style="list-style-type: none"> ▪ To maintain sample quality and purity ▪ To avoid contamination ▪ To recognize the factors affecting your results
8. Best Practices for Cell-based Sample Preparation Module	Have you ever wondered if sample preparation could be the source of variance you see in your cell-based assays?	<ul style="list-style-type: none"> ▪ The aspects that affect your cell-based assays with Incucyte or iQue® screener ▪ How to reduce variability between experiments ▪ How to speed up the workflow and improve your cell health with smart pipetting practices
9. Electronic Pipettes and Digital Workflows Module	Planning to switch to electronic pipettes? Or already done so but you need some tips and tricks about the different functionalities of an electronic pipette and how it can together with digital workflows make pipetting easier but also more accurate?	<ul style="list-style-type: none"> ▪ Functions and features of an electronic pipette ▪ How to get the most out of your electronic pipette ▪ How digital workflows can speed up your pipetting work
10. Qualification Training According to ISO 8655:10 Module	In the latest revision of ISO 8655, part 10 was introduced regarding pipette user training and qualification requirements. This module is designed to support continuous learning and development of pipette users.	<ul style="list-style-type: none"> ▪ (Re-)Qualification according to the the ISO 8655:10 requirements ▪ Certificate upon completion


Pipetting Recommendations



Hold the pipette in a vertical position during aspiration.



Avoid contamination with Safe-Cone Filters.



Tacta® value lock can be used in two ways, by pressing it while changing the volume, or by sliding it up to open it, and back to lock it.

Preparations Before Pipetting

- Use the tip specified by the manufacturer.
- Ensure that the pipette and the tip have been tested according to ISO 8655 and the tip is seated correctly.
- Make sure pipettes have been correctly calibrated.
- Make sure that the volume that you are pipetting is within usable volume range (ISO 8655).
- Check that the pipette, tip and liquid are all at the same temperature.
- When pipetting liquids with temperatures different to the ambient temperature, do not pre-rinse the tip. Change the tip after each pipetting.
- Ensure that any fluid viscosity variations have been accounted for and the correct technique is employed, i.e. reverse pipetting.
- If handling infectious or radioactive agents make sure appropriate shielding and other precautions protect the operator.
- Use Safe-Cone Filter in the tip cone whenever possible.
- Use racked tips when possible for best ergonomics

While Pipetting

- Hold the pipette in a vertical position, while aspirating. Tilting the pipette at an angle causes a volume greater than the set volume of liquid to enter the tip.
- In most cases, pre-rinsing of the tip is recommended, to achieve accurate results. Do not pre-rinse the tip, if the temperature of the liquid is different to the ambient temperature.
- When aspirating fluid, the pipette tip should normally be immersed to a depth of 2–3 mm.
- When using a mechanical pipette, operate the piston with a smooth and consistent thumb action, for repeated results without foaming or bubbles.
- Tilting the pipette and placing the tip against the inside wall, when possible, can help to get the liquid out from the tip.
- Ensure that the pipette blow out action is fully activated.
- Ensure that the volume is still set at the required position. A pipette with a volume locking mechanism is recommended, in order to avoid accidental volume change during pipetting.
- Avoid laying the pipette on its side while there is liquid in the tip. It may seep up into the mechanism.



Charging it while pipetting is possible with Sartorius electronic pipettes.



Load the tip onto the pipette carefully and take advantage of the Optiload tip loading mechanism.



Clean the pipette before sending it to service.

Other Precautions

- Store the pipette on a stand when not in use – see page 46, on pipette stands, for more information. Electronic pipettes should be returned to their charging stands.
- Remember that the pipette and tip combination is one complete system.
- Avoid dropping the pipette or allowing contact with dirt or grease.
- Change the Safe-Cone Filter regularly (recommendation after 50 to 250 pipetting cycles), and in every case of over-aspiration.
- Never strike the tip cone against the tip tray when loading the tip, as this can damage the pipette.
- Avoid exposing the unit to extreme temperature changes, humidity and dust (operating temperature from 15 °C to 40 °C).
- Service the pipette regularly.
- Clean the pipette thoroughly before sending it in for service. Decontaminate the pipette with 70 % ethanol. Notify the service personnel of the purpose for which the instrument has been used. Postal services may refuse to deliver instruments used for hazardous materials. Make sure that a qualified person services the pipette.



Calibration and Maintenance

Table of Contents

Calibration and Maintenance

Pipette Calibration and Maintenance Services	64
Pipette Decontamination Procedure	65
Autoclaving Instructions	66

Pipette Calibration and Maintenance Services

Why is Calibration and Maintenance Needed?

Pipette calibration is a fundamental part of Good Laboratory Practice (GLP) and quality systems and must be considered a vital part of any laboratory regime where precise volumes of fluid need transferring or diluting. Pipette performance is measured as accuracy and precision or how close the dispensed volumes are to the target and how close the results are to each other.

Adopting a regular calibration and maintenance routine for your pipette has the following benefits:

Confidence

Your pipettes are operating correctly with the accuracy and precision you know.

Reliability

With maintenance and calibration and operational qualification you are able to trust the instrument's operational capability.

Efficiency and Effectiveness

With properly working pipettes you can work uninterrupted and be more efficient.

Sartorius Service Centres

It is vitally important to Sartorius that our customers receive world class service and support, from the first phone call to the moment the engineer leaves, with the customer's equipment in perfect working condition.

Sartorius provides a global network of service centers for calibration of all makes and models of pipettes, burettes, bottle top dispensers and other liquid handling instruments. With several decades of experience in pipettes and liquid handling instrument services, the global organization provides world class services tailored to customers' needs, simultaneously fulfilling standards and regulatory body requirements. With the Sartorius concept of "all makes and models liquid handling services" you can be certain that your instruments are handled according to the international ISO 8655 standard defining the use and calibration of pipettes and other liquid handling instruments. Sartorius provides reliability of operation, reduced instrument downtime and confidence, so that your work is according to the strictest requirements.



Pipette Decontamination Procedure

Mechanical Pipettes (Tacta® and Proline® Plus)



1. Unscrew the tip ejector collar counter clockwise and remove it. Remove the Safe-Cone Filter if fitted.



2. Unscrew the tip cone holder counter clockwise and carefully remove it along with the tip cone.



3. The piston is left on the mechanical pipette.
4. Place the tip ejector collar, tip cone holder, tip cone and tip cone cylinder into a beaker containing 70 % ethanol and leave for at least 30 minutes.



5. After performing the procedure described above, remove the components from the beaker and rinse them with distilled water, then dry preferably with warm air, for at least an hour.
6. Re-grease the piston as described in the instruction manual. Replace all components including new filter if fitted.

Electronic Pipettes (Picus® 2)



1. Unscrew the tip ejector collar counter clockwise and remove it. Remove the Safe-Cone Filter if fitted.



2. Unscrew the tip cone holder counter clockwise and carefully remove the tip cone holder, tip cone and spring.



3. Unscrew the piston counter clockwise from the pipette.
4. Place the tip cone, tip cone holder, tip ejector collar, piston and spring into a beaker containing 70 % ethanol and leave for at least 30 minutes.



5. After performing the procedure described above, remove the components from the beaker and rinse them with distilled water, then dry preferably in warm air, for at least an hour.
6. Re-grease the piston as described in the instruction manual. Replace all components including the new filter if fitted.

Note: When opening pipette lower parts, as a matter of routine the parts should be inspected for possible wear. Always check the pipette performance after opening the lower parts.

Autoclaving Instructions



Tacta® and Proline® Plus Mechanical Pipettes

The entire Tacta® and Proline® Plus mechanical pipette can be steam sterilized by autoclaving at 121°C (252°F), at 1 bar (15 p.s.i.), for 20 minutes. The dispensing head of the multi-channel pipettes must be unscrewed 360° counter clockwise before autoclaving.

- Remove the Safe-Cone Filter (if fitted)
- Place the pipette into the sterilisation bag and place it into the autoclave
- After autoclaving the pipette must be cooled down and left to dry overnight before use

It is recommended that you check the performance of the pipette after every autoclaving, and grease the piston | seal of the pipette after every 10th autoclaving.

Tips and Tip Boxes

- Place the bulk tips into the sterilisation bag and the tip tray as such in the autoclave
- Autoclave, at 1 bar, for 20 minutes at a temperature of 121°C
- Cool before use

Picus® 2 Electronic Pipettes Lower Parts

The dispensing head (tip ejector collar, tip cone holder, tip cone, spring and piston) of the single-channel and multi-channel models (except for multi-channel 1,200 µL) can be autoclaved (121°C, at 1 bar, for 20 minutes). These parts can be autoclaved as one unit or separately as individual parts. It is also possible to clean the parts and grease the piston prior to autoclaving.

- Remove the Safe-Cone Filter (if fitted).
- Put the dispensing head into the sterilisation bag and place it into the autoclave
- After autoclaving the parts must be cooled down and left to dry before use

It is recommended that you check the performance of the pipette after every autoclaving, and grease the piston | seal of the pipette after every 10th autoclaving.



Note:

- Excessive heat or length of time may damage the products. Never place the handle part of the Picus® 2 into the autoclave
- The lower ends of multi-channel pipettes are not interchangeable between 8 and 12-channel pipettes
- The cover of the tip tray should be closed during autoclaving



Germany

Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Strasse 20
37079 Göttingen
Phone +49 551 308 0

USA

Sartorius Corporation
3874 Research Park Drive
Ann Arbor, MI 48108
Phone +1 734 769 1600

 For further information, visit
[sartorius.com](https://www.sartorius.com)