


## Operating Instructions

# BioPAT® MFCS 4

Recipe Control Module  
Software



1000035948



# SARTORIUS



# Contents

<b>1 About these Instructions</b>	<b>5</b>
1.1 Accompanying Documents	5
1.2 Symbols Used	5
1.2.1 Warnings in Operation Descriptions	5
1.2.2 Other Symbols	5
<b>2 Operating Concept</b>	<b>6</b>
2.1 Module Description	6
2.2 Configuring the Module	6
2.2.1 Updating the Basic Version	6
2.2.2 Licensing and Activating the Module	6
<b>3 User Interface</b>	<b>7</b>
3.1 Function Panes	7
3.2 "RECIPES" Section in the "ADMINISTRATION" Function Pane	7
3.3 Additional Selection Keys and Selection Fields	8
3.3.1 [ADMINISTRATION] Function Pane	8
3.3.2 [MONITORING] Function Pane	11
3.4 Drag and Drop Function	14
3.4.1 Changing the Sequence of Operations	14
3.4.2 Positioning Phases in the Planned Sequence	15
3.4.2.1 Positioning Phases	16
3.4.2.2 Repositioning Phases	17
3.4.2.3 Adding a Phase	19
3.5 Linking the Phases	20
3.5.1 Creating a Sequential Path	22
3.5.2 Creating an Alternative Path	23
3.5.3 Creating a Parallel Path	24
3.5.4 Creating a "Jump" Link	25
3.5.5 Disconnecting a Link	27
<b>4 Recipes Module</b>	<b>29</b>
4.1 [RECIPES (#/#)] Menu	29
4.2 Recipe Editor	30
4.2.1 [META DATA] Input / Selection Screens	31
4.2.2 Sequence List of Operations	32
4.2.3 Planned Sequence and Library	33
4.2.4 [EDIT PARAMETER PHASE] Input / Selection Screens	35
4.2.4.1 [META DATA] Input / Selection Field	35
4.2.4.2 [PARAMETERS] Input / Selection Field	36
4.2.4.3 [ADD PARAMETER (#/#)] Selection Menu	38
4.2.5 [EDIT OPERATOR PROMPT PHASE] Input / Selection Screen	39
4.2.5.1 [META DATA] Input / Selection Field	39
4.2.5.2 [DECISIONS] Input / Selection Field	40
4.2.5.3 [NOTIFICATION] Input / Selection Field	40
4.2.6 [EDIT PROFILE PHASE] Input / Selection Screen	41
4.2.6.1 [META DATA] Input / Selection Field	42
4.2.6.2 [PROFILE POINTS] Input / Selection Field	43
4.2.6.3 [PROFILE] Input / Selection Field	44
4.2.6.4 [SELECT PARAMETER (#/#)] Selection Menu	49
4.2.7 [EDIT INSTRUMENT PHASE] Input / Selection Screen	50
4.2.7.1 [SELECT PHASE] Input / Selection Screen	51
4.2.8 [EDIT WAIT PHASE] Input / Selection Screen	52

4.3	Formula Editors .....	53
4.3.1	[SELECT VARIABLE] Selection Menu .....	59
4.3.2	[SELECT CALCULATION] Selection Menu .....	60
4.3.3	Functions and Operations .....	61
4.4	Creating a Recipe .....	63
4.4.1	Opening the Recipe Editor .....	64
4.4.2	Entering Meta Data .....	64
4.4.3	Creating Operations .....	65
4.4.4	Configuring Phases .....	65
4.5	Special Recipe Functions .....	80
4.5.1	Configuration of the "Operator Prompt Phase" .....	80
4.5.2	Setting the Process Time Automatically to Zero .....	82
4.5.3	Setting the Alarm Parameters .....	83
4.5.3.1	Disabling / Enabling the Alarming .....	83
4.5.3.2	Modifying Alarm Limits .....	86
4.5.4	Validating the Recipe Syntax .....	88
4.6	Automated Batch Process .....	91
4.6.1	[RECIPES] Menu in the [MONITORING] Function Pane .....	93
4.6.1.1	[OPERATIONS] Sequence List .....	93
4.6.1.2	Planned Sequence .....	95
4.6.1.3	[RECIPES] Preview Window .....	104
4.6.2	Starting an Automated Batch Process .....	106
4.6.3	Disconnecting a Link .....	107
4.6.4	Starting a Phase Manually .....	108
4.6.5	Pausing and Continuing a Phase .....	109
4.6.6	Stopping and Resetting a Phase .....	111
4.6.7	Stopping a Phase and Activating a Transition .....	112
4.6.8	Ending an Automated Batch Process .....	113
4.7	Application Example for Resetting a Calculation via a Parameter Phase .....	114
4.8	Display of a Phase Configuration .....	117
4.8.1	Opening and Closing the Overview .....	118
4.9	Exporting and Importing a Recipe .....	119
4.9.1	Exporting a Recipe .....	120
4.9.2	Importing a Recipe .....	121
4.10	Phase Structures .....	122
4.10.1	Sequential Path .....	122
4.10.2	Alternative Path .....	123
4.10.3	Parallel Paths .....	124
4.10.4	Path with a "Jump" Link .....	125
5	Index .....	126



# 1 About these Instructions

These instructions form part of the Recipe Control Module of the BioPAT® MFCS 4 software.

Fundamental information about the basic version of BioPAT® MFCS 4 can be found in the BioPAT® MFCS 4 Operating Instructions. Fundamental information includes:

- Software use
- Default settings and installation
- User Interface
- ADMINISTRATION, MONITORING and ANALYSIS function panes

Ensure that these instructions are read by anyone working with Recipe Control Module.

- You can order a printout of these instructions at an additional cost. Contact Sartorius in order to do so.

The figures and illustrations from the Recipe Control Module in these instructions were created with English as the set language (for language settings, see BioPAT® MFCS 4 Operating Instructions, “System Settings”).

## 1.1 Accompanying Documents

- In addition to these instructions, observe the following documents:
  - BioPAT® MFCS 4 Operating Instructions
  - BioPAT® MFCS 4 Calculation Module Operating Instructions

## 1.2 Symbols Used

### 1.2.1 Warnings in Operation Descriptions

---

#### NOTICE

Denotes a danger with the risk that property damage may result if it is **not** avoided.

---

### 1.2.2 Other Symbols

- Required action: Describes activities that must be carried out.
- ▷ Result: Describes the result of the activities carried out.
- [ ] Text inside brackets refers to control and display elements
- [ ] Text inside brackets indicates status, warning, and error messages.

## 2 Operating Concept

### 2.1 Module Description

The Recipe Control Module is an optional program module that enhances the functional scope of the basic version of BioPAT® MFCS 4.

Once the Recipe Control Module is activated, the following additional functions can be used:

- Creation of recipes for automating batch processes
- Editing and deleting recipes
- Linking recipes with batch processes

### 2.2 Configuring the Module

#### 2.2.1 Updating the Basic Version

##### Procedure

- To be able to use the advanced functions of the Recipe Control Module: Update the basic version of BioPAT® MFCS 4 (see BioPAT® MFCS 4 Operating Instructions, Chapter “Program Updates”).

#### 2.2.2 Licensing and Activating the Module

Recipe and recipe-controlled batch process created and configured while using the 90-days demo version continue to be saved even after the demo period expires. However, access to the configured recipes and recipe-controlled batch processes is **not** possible after the demo period ends.

The configured recipe and recipe-controlled batch processes can again be accessed once the Recipe Control Module has been licensed and activated.

##### Procedure

- To have permanent use of the advanced functions of the Recipe Control Module: License and activate the Recipe Control Module (see BioPAT® MFCS 4 Operating Instructions, Chapter “Program Licensing”).

## 3 User Interface

Fundamental information about the user interface of the basic program can be found in the BioPAT® MFCS 4 User Manual in the following chapters:

- Program Start
- Start Screen
- Navigation Aids
- Selection Keys

### 3.1 Function Panes


Fundamental information about the function panes of the basic program can be found in the BioPAT® MFCS 4 User Manual in the following chapters:

- ADMINISTRATION Function Pane
- MONITORING Function Pane
- ANALYSIS Function Pane

### 3.2 “RECIPES” Section in the “ADMINISTRATION” Function Pane

After licensing and activating the Recipe Control Module, the [RECIPES] section is displayed in the [ADMINISTRATION] function pane.

- ▶ Click on the [RECIPES] button.
- ▷ The [RECIPES] overview is shown. Existing recipes are listed in the [RECIPES] overview.



### 3.3 Additional Selection Keys and Selection Fields

The additional selection keys and selection fields are used in the [ADMINISTRATION] and [MONITORING] function panes for the automated sequence of a batch process.

#### 3.3.1 [ADMINISTRATION] Function Pane

In order for a batch process to run automatically, operations, phases and conditions for the batch process are modeled in the form of a recipe. The recipe is created in the [ADMINISTRATION] function pane.






Fig. 1: Recipe editor







Pos.	Field	Description
1	META DATA	Meta data input area: Input and selection screens for meta data
2	OPERATIONS	List of operations: Input and selection screens for the operations
3		Planned sequence for the marked operation: Input and selection screens for the phases and conditions
4	LIBRARY	Library: Phase types for development in the planned sequence

The function elements used to create a recipe are displayed in the tables below:



### [OPERATIONS] menu







Function element	Description
	Displays an existing operation.
	Displays an existing and marked operation. The phases of the operation are displayed in the planned sequence and can be edited.

### Library

Function element	Description
	Selects the “Parameter phase” phase type. The phase can be configured with a formula for a one-off calculation (see Chapter “4.3 Formula Editors”, page 53).
	Selects the “Operator prompt phase” phase type.
	Selects the “Instrument phase” phase type.
	Selects the “Profile phase” phase type.
	Selects the “Wait phase” phase type.
	Selects the “End” phase type.

### Planned sequence

Function element	Description
	“Start” phase: The phase is automatically set as the start phase of the operation.
	“Parameter phase”: The phase is a selection element in the library.

Function element	Description
	"Operator prompt phase": The phase is a selection element in the library.
	"Instrument phase": The phase is a selection element in the library.
	"Wait phase": The phase is a selection element in the library.
	"Profile phase": The phase is a selection element in the library.
	"End" phase: The phase is a selection element in the library and ends the operation.
	Transition: A transition links phases to each other.

### 3.3.2 [MONITORING] Function Pane

The [RECIPES] menu presents the sequence of one or more than one started and automated batch process.

[RECIPES] menu: Started batch process with one unit




Fig. 2: Active [RECIPES] menu for the started batch process of one (example)

Pos.	Field	Description
1	OPERATIONS	List of operations of the automated batch process
2		Planned sequence for the marked operation
3	UNIT DISPLAY	[UNIT DISPLAY] preview window
4	TREND	[TREND] preview window
5	ALARMS	[ALARMS] preview window

[RECIPES] menu: Started batch process with one unit group




Fig.3: Active [RECIPES] menu for a started batch process of a unit group (example)

Pos.	Field	Description
1	<div><div>ECB1</div><div><div>U1</div><div>0</div></div></div> <div><div>---</div><div><div>U2</div><div>0</div></div></div> <div><div>---</div><div><div>U3</div><div>0</div></div></div> <div><div>ECB3</div><div><div>U4</div><div>0</div></div></div>	<p>Overview with units of the unit group (example):</p> <p>Recipe-controlled batch process with unit “U1” and recipe “ECB1” is started. The recipe is displayed in the planned sequence.</p> <p>Batch process with unit “U2” (<b>not</b> recipe-controlled) is started. Clicking the button displays <b>no</b> recipe in the planned sequence.</p> <p><b>No</b> batch process is started for unit “U3”. Clicking the button displays <b>no</b> recipe in the planned sequence.</p> <p>Recipe-controlled batch process with unit “U4” and recipe “ECB3” is started. Clicking the button displays the recipe in the planned sequence.</p>
2	OPERATIONS	List of operations of the automated batch process
3		Planned sequence for the marked operation
4	UNIT DISPLAY	[UNIT DISPLAY] preview window



Pos.	Field	Description
5	TREND	[TREND] preview window
6	ALARMS	[ALARMS] preview window

#### [RECIPES] preview window

The [RECIPES], [UNIT DISPLAY], [TREND] and [ALARMS] menus are displayed in one window. Only one menu is ever active (settings can be configured in this menu), while the other menus are displayed as previews in the right-hand pane (for toggling the active menu display, see BioPAT® MFCS 4 User Manual, Chapter “MONITORING Function Pane”).




Fig. 4: [RECIPES] preview window for the started batch process (example)

Pos.	Field	Description
1	RECIPES	[RECIPES] preview window

## 3.4 Drag and Drop Function

Using the drag-and-drop function, the sequence of operations can be adjusted in the [OPERATIONS] menu and the phase positioned in the planned sequence.

### 3.4.1 Changing the Sequence of Operations


Operations in the sequence list are numbered according to the sequence, e.g. 01 Operation A, 02 Operation B, 03 Operation C.

#### Example


Operation C should follow Operation A.


#### Procedure

- Select [03 Operation C] in the sequence list.



- Move [03 Operation C] between the position of [01 Operation A] and [02 Operation B] by holding down the left mouse button.





- ▶ Release the left mouse button.
- ▷ Operation C follows Operation A.
- ▷ Operation B follows Operation C.

### 3.4.2 Positioning Phases in the Planned Sequence

The planned sequence is constructed as a grid plan with fixed positions for positioning the phases. The phases are positioned in the planned sequence.




Fig. 5: Planned sequence with unlinked phases (example)

#### Pos. Description

- | Pos. | Description   |
|------|---|
| 1    | Displays the automatically positioned "Start" phase.    |
| 2    | Displays the positioned "Parameter phase".              |
| 3    | Shows the dragging and positioning of the "Wait phase". |
| 4    | Indicates a free position in the planned sequence grid. |
| 5    | Vertical scrolling in the planned sequence              |
| 6    | Horizontal scrolling in the planned sequence            |

### 3.4.2.1 Positioning Phases


Phases are dragged from the library as templates, dropped into the planned sequence and positioned.

#### Example


The “Wait phase” is added to the operation and follows the “Parameter phase”.

#### Procedure


- Drag the “Wait phase” from the [LIBRARY] into the planned sequence by holding down the left mouse button.



- Place the “Wait phase” in a free position in the planned sequence.



- The “Wait phase” is positioned.



### 3.4.2.2 Repositioning Phases

Phases that are already positioned in the sequence and are linked to each other via a transition can be moved to a new position within the planned sequence and positioned. The link remains unchanged.

#### Example

The “Wait phase” is positioned in a new position in the planned sequence.





Fig. 6: Linked phases

#### Pos. Description


- |   |   |
|---|---|
| 1 | Displays the “Transition” function element, creates the link and defines the condition. |
| 2 | Displays the “Wait phase” in the current position.                                      |

#### Procedure


- Lift the “Wait phase” from its current position by holding down the left mouse button.



- Place the “Wait phase” in a free position in the planned sequence.



- ▷ The “Wait phase” is repositioned.
- ▷ The “Wait phase” is also linked to the “Parameter phase” via the “Transition” function element.



### 3.4.2.3 Adding a Phase

A phase can be added to a planned sequence, in which linked phases are already present, and linked with those phases. If the phase is to be positioned between two phases, then the link between these two phases must be removed.

#### Example

The “Wait phase” is added to the planned sequence between the “Start” phase and the “Parameter phase”.





Fig. 7: Positioning attempt of a phase (example)


#### Pos. Description

- 1 Displays the “Transition” function element. The function element creates the link and defines the condition.
- 2 Indicates the attempt to position the “Wait phase” between the linked “Start” phase and “Parameter phase”.


#### Procedure

- To disconnect a link: In the “Transition” function element, click on the [Remove] button.
- The “Transition” function element is deleted and the link between the “Start” phase and the “Parameter phase” is disconnected.





- Place the "Wait phase" from the library into a free position between the "Start" phase and the "Parameter phase".



- The "Wait phase" is positioned between the "Start" phase and the "Parameter phase".
- The phase can be linked with the "Start" phase and "Parameter phase" (for linking phases, see Chapter 3.5, page 20).

### 3.5 Linking the Phases

An operation is only run when all phases are linked with each other according to the sequence of the batch process. Linked phases form different paths, depending of the type of phase structure (for overview of the phase structures, see Chapter 4.10, page 122). The following paths can be created in the planned sequence:

- Sequential path
- Alternative path
- Parallel path

For sequential and alternative paths, the process for linking two phases generates the "Transition" function element.

The condition between two phases is defined by means of the "Transition" function element. The phase after the "Transition" function element is only started once the condition is met.









Fig. 8: 3 phases in the planned sequence (example)

Pos.	Symbol	Description
1		Indicates the phase:
		Removes* the phase and all connections.
		Opens* the phase configuration. The "Start" and "End" phases <b>cannot</b> be configured.
2		Displays the output of the phase link:
		Links the phase. The input / output can be linked with another phase.
		Links the phase. The output of the phase is selected for linking with an input / output of another phase. The output of the "Start" phase can be linked exclusively with the input of another phase.
3		Displays the "Transition" function element and links phases.
		Removes* the function element. The link and condition between the phases is removed.
		Opens* the "Condition" configuration and determines the condition between the phases with the formula editor.

Pos.	Symbol	Description
4		Displays the input of the phase link:
		The input can be linked with another phase.
		Connection of the input with the input / output of another phase is <b>not</b> possible.


\* Symbol is visible when approached with the mouse pointer.

### 3.5.1 Creating a Sequential Path


The “Start” phase is linked with the “Parameter phase”.


#### Procedure

► In the “Start” phase, click on the [Link] button.



► In the input of the “Parameter phase”, click on the [Link] button.





- ▷ The "Start" phase and "Parameter phase" are linked.
- ▷ The "Transition" function element is set.

### 3.5.2 Creating an Alternative Path

The alternative path is linked with the "Start" phase.





Fig. 9: Path and alternative path in the planned sequence (example)


#### Pos. Description

- |   |  |
|---|--|
| 1 | Indicates the path that has already been configured. The "Parameter phase" is linked with the "Start" phase. |
| 2 | Indicates the alternative path. The "Parameter phase" is <b>not</b> linked with the "Start" phase.           |


#### Procedure

- ▷ In the input of the "Parameter phase" of the alternative path, click on the [Link] button.





► In the “Start” phase, click on the [Link] button.



- The “Start” phase and “Parameter phase” of the alternative path are linked.
- The “Transition” function element is set for the alternative path.
- In the automated batch process, the path for which the transition condition is first met is processed.

### 3.5.3 Creating a Parallel Path




Fig.10: Path and parallel path in the planned sequence (example)


Pos.	Description
1	Indicates the path that has already been configured. The “Parameter phase” is linked with the “Start” phase.
2	Indicates the parallel path. The input of the “Parameter phase” is <b>not</b> linked with the input of the “Parameter phase” that has already been configured.

## Requirements


The phase of the path that has already been configured must be linked with a preceding phase, as a sequential phase.

## Procedure


- In the input of the “Parameter phase” of the path that has already been configured, click on the [Link] button.



- In the input of the “Parameter phase” of the parallel path, click on the [Link] button.



- The inputs of the phases are linked.
- In the automated batch process, both paths are processed when the condition of the common transition is met.



### 3.5.4 Creating a “Jump” Link


The “Jump” link can fulfill the “Jump” and “Loop” functions:

- Loop: A phase can be linked with itself. An additional transition is created. If the condition for the transition is met, the phase is repeated.
- Jump: A phase can be linked with another phase which does not directly follow the first phase. An additional transition is created. If the condition for the transition is met, the linked phase is started. This phase can be located
  - in the same path
  - in an alternative path
  - in a parallel path


This example demonstrates how a “Jump” link is created with the “Loop” function.

## Procedure


- In the output of the “Wait 20 min” phase, click on the [Link] button.



- In the input of the “Wait 20 min” phase, click on the [Link] button.



- ▷ The output and input of the “Wait 20 min” phase are linked.
- ▷ If the condition for the transition (1) is met: the next phase starts.
- ▷ If the condition for the transition (2) is met: The “Wait 20 min” phase is again started.





### 3.5.5 Disconnecting a Link

#### Example 1


The link between the “Start” phase and the “Parameter phase” is disconnected by deleting the “Transition” function element.

#### Procedure

- In the “Transition” function element, click on the [Remove] button.



- To confirm the deletion process: Click on the [YES] button.




- The link between the “Start” phase and the “Parameter phase” is disconnected.
- The “Transition” function element is deleted.

### Example 2


The link between the “Start” phase, the “Parameter phase”, and the “Wait phase” is disconnected by deleting the “Parameter phase”.

#### Procedure


► In the “Parameter phase”, click on the [Remove] button.



► To confirm the deletion process: Click on the [YES] button.



- ▷ The links between the phases are disconnected.
- ▷ The “Parameter phase” is deleted.
- ▷ The “Transition” function elements are deleted.





## 4 Recipes Module

### 4.1 [RECIPES (#/#)] Menu

The [RECIPES] menu lists the existing recipes. The operating elements can be used to create new recipes or delete or edit existing recipes. Recipes can be exported and imported.




Fig. 11: [RECIPES] menu

Pos.	Symbol	Description
1		Sorts the entries by "Name" and "Short Name".
2		Displays the list of created recipes.
3		Shows marked and selected recipe. Double-click to open the editing window for the recipe.
4		Filters the entries in the list by "Name", "Short Name", "Originator", "Medium", "Strain", and "Product".
5		Refreshes the overview.
6		Exports the selected recipe to an XML file.
7		Imports a recipe from an XML file.
8		Deletes the selected recipe.
9		Opens the editing window.
10		Creates a new recipe.

## 4.2 Recipe Editor

The recipe editor can be used to create new recipes and edit stored recipes. The recipe is defined at three levels:

Level	Description
1	Meta data: The name of the recipe and the information on the batch process are stored in the meta data.
2	Operations of the recipe: An operation consists of phases. When all phases of an operation have been processed, the next operation starts.
3	Planned sequence: The planned sequence consists of phases and conditions that are linked to each other by paths. An operation starts with the "Start" phase and ends with the "End" phase.




Fig.12: Recipe editor, example: showing a stored recipe

Pos.	Field	Description
1	META DATA	Level 1: Input and selection screens for meta data
2	OPERATIONS	Level 2: Sequence list of operations
3		Level 3: Planned sequence with phases for an operation
4	LIBRARY	Level 3: Library of phases

4.2.1 [META DATA] Input / Selection Screens

NEW RECIPE

META DATA

GENERAL SETTINGS

Name\*

Short name\*

Originator

DETAILS

Medium

Strain






Product

DESCRIPTION

Description

Fig. 13: [META DATA] input / selection screens for a new recipe

[GENERAL SETTINGS], [DETAILS], and [DESCRIPTION] Input / Selection Screens

Field	Symbol	Description
META DATA	 	Expands and collapses the input / selection screen.
Name*		Name of the recipe (max. 48 characters).
Short name*		Short name of the recipe (max. 4 characters).
Originator		Name of the recipe controller (max. 48 characters):
		Selects a stored name.
Medium		Name of the medium (max. 48 characters):
		Selects a stored medium.
Strain		Name of the strain (max. 48 characters):
		Selects a stored strain.
Description		Description of the automated batch process (max. 255 characters)

\* Required information

4.2.2 Sequence List of Operations










Fig.14: Sequence list of operations (example)

Field	Symbol	Description
OPERATIONS	 	Collapses and expands the sequence list.
01 Operation 1 02 Operation 2 03 Operation 3		Indicates the position and name of the operation. Example: Once "Operation 1" is complete, "Operation 2" starts. Once "Operation 2" is complete, "Operation 3" starts. Once "Operation 3" is complete, the automated sequence for the batch process is finished.
		Removes the operation from the sequence list: The subsequent operation is linked to the preceding operation.
		Shows the [EDIT OPERATION] input screen.
		Creates the first operation and displays the [ADD OPERATION] input screen.
		Creates another operation and displays the [ADD OPERATION] input screen.
Name*		Name of the operation (max. 48 characters)
Mode		Runs the operation exclusively automatically.
Description		Description of the operation (max. 255 characters).

\* Required information

### 4.2.3 Planned Sequence and Library













Fig. 15: Planned sequence and library (example)

Field	Symbol	Description
LIBRARY		Collapses and expands the library.
Parameter phase		Dragging the phase to a free position in the planned sequence adds the phase to the sequence.
Wait phase		Dragging the phase to a free position in the planned sequence adds the phase to the sequence.
Operator prompt phase		Dragging the phase to a free position in the planned sequence adds the phase to the sequence.
Profile phase		Dragging the phase to a free position in the planned sequence adds the phase to the sequence.
Instrument phase		Dragging the phase to a free position in the planned sequence adds the phase to the sequence.
End		Dragging the phase to a free position in the planned sequence adds the phase to the sequence.
		The phase is already set in the planned sequence. An operation in the planned sequence can only have one "End".
Start		Starting point of the operation
Parameter phase		Controls the phase to the set parameters for the selected control module variable.

Field	Symbol	Description
		Controls the phase in "Semi-Automatic" mode to the set parameters for the selected control module variable. The phase starts after confirmation from the user.
Operator prompt phase		Controls the process sequence depending on a user-specific decision (configuration of a maximum of three decisions is possible).
Wait phase		Sets the wait time until the next phase.
Profile phase		Controls the phase to the configured setpoint profile for the selected control module variable.
		Controls the phase in "Semi-Automatic" mode to the configured setpoint profile for the selected control module variable. The phase starts after confirmation from the user.
Instrument phase		Executes a sequence on a device (execution using the default parameters stored on the device).
		Executes a sequence in "Semi-Automatic" mode on a device (execution using the default parameters stored on the device). The phase starts after confirmation from the user.
End		End point of the operation.
		Displays the "Transition" function element: <ul style="list-style-type: none"> <li>- Links a phase (or phases).</li> <li>- Defines the condition for starting the next phase(s) (see Chapter "4.3 Formula Editors", page 53).</li> </ul>

4.2.4 [EDIT PARAMETER PHASE] Input / Selection Screens

The parameters for the “Parameter phase” are selected and set in the input and selection screens.

EDIT PARAMETER PHASE	
<b>META DATA</b>	
Name*	pO2 Setpoint 20%
Mode	Automatic Semi-Automatic
Description	

PARAMETERS	
Name	Value*
pO2.Mode	auto
pO2.Setpoint	20 %sat

SAVE CANCEL

Fig.16: Input / selection screens for the “pO2.Mode” and “pO2.Setpoint” control modules (example)

Field	Description
META DATA	Shows the meta data.
PARAMETERS	Displays the selected control module variables.

4.2.4.1 [META DATA] Input / Selection Field

META DATA	
Name*	
Mode	Automatic Semi-Automatic
Description	

Fig.17: [META DATA] input / selection field

Field	Description
Name*	Name of the phase (max. 48 characters)
Mode	Selection of the phase control mode:
Automatic	Runs the phase automatically.
Semi-Automatic	Runs the phase after confirmation from the user.
Description	Description for the phase (max. 255 characters).

\* Required information

#### 4.2.4.2 [PARAMETERS] Input / Selection Field

Parameters that are sent to the device during the phase are defined in the input and selection field. Depending on the parameter, the corresponding operating modes, such as controller mode and remote mode, and the parameter values are set. The parameters in the list are processed one after the other from top to bottom.

A process parameter (e.g. a setpoint controller) can be assigned a numerical value or a formula. When assigning the process parameters with a formula, a one-off calculation may be executed in an automatic process, e.g. to save a current process value in a variable and reuse it in subsequent calculations. One application example is the calculation of setpoint values for an exponential feeding profile in fed-batch cultivations.























PARAMETERS					
Name		Value*			
STIRR_1.AlarmDeadband [ rpm ]	▼	rpm			
STIRR_1.AlarmHigh [ rpm ]	▼	rpm			
STIRR_1.AlarmHighHigh [ rpm ]	▼	rpm			
STIRR_1.Alarming	▼	Select value ▼			
STIRR_1.AlarmLow [ rpm ]	▼	rpm			
STIRR_1.AlarmLowLow [ rpm ]	▼	rpm			
STIRR_1.Mode	▼	Select mode ▼			
STIRR_1.Output [ % ]	▼	% ⋮			
STIRR_1.Setpoint [ rpm ]	▼	rpm ⋮			
STIRR_1.State	▼	Select state ▼			

Fig.18: [PARAMETERS] input / selection field for the "STIRR\_1" control modules (example)

Field	Symbol	Description
		Shows the [ADD PARAMETER (#/#)] selection menu.
	 	Moves the selected parameter(s) up or down the list.
		Deletes the selected parameter(s) from the list.
Name		Displays the name and parameter of the control module.
	 	Selects the parameter.
Value*		Displays the parameter (depending on the type of control module):



Field	Symbol	Description
Select mode	 	Selects the “auto”, “off”, or “manual” controller mode for the “Mode” parameters.
Select state	 	Selects “local”, “remote”, or “calculated” <sup>1</sup> remote mode for the “State” parameters.
Select value	 	Deactivate (“Off”) or activate (“On”) the alarming.
Value	 	Selects the “Value” function for the “Value”, “Setpoint”, “Deadband”, “Output”, “OutputMax”, “OutputMin”, “Td”, “Ti”, “Xp”, “AlarmDeadband”, “AlarmHigh”, “AlarmHighHigh”, “AlarmLow”, “AlarmLowLow” parameters and enters a numerical value.
Calculation...	 	Selects the “Calculation” function for the “Value”, “Setpoint”, “Deadband”, “Output”, “OutputMax”, “OutputMin”, “Td”, “Ti”, “Xp” parameters and displays the formula editor. The phase can be configured with a formula for a one-off calculation (see Chapter “4.3 Formula Editors”, page 53).
Reset	 	Resets a recursive calculation that uses the “getValue()” function to its start value.

\* Required information

<sup>1</sup> Available if the Calculation Module is licensed and activated

#### 4.2.4.3 [ADD PARAMETER (#/#)] Selection Menu

The selection menu lists the parameters of all units that can be selected for the "Parameter phase". The recipe can only be selected for the batch process if parameters that are configured for the unit are selected.




Fig.19: [ADD PARAMETER (#/#)] selection menu (example)

Pos.	Field	Symbol	Description
1	Parameter		Shows names of the parameters of the units selected in the "Unit" filter.
		▲ ▼	Sorts the entries into alphabetical order (ascending and descending)
	Engineering unit		Displays the engineering unit (if any) of the parameter.
	Unit		Displays the units for which the parameter is configured.
2	FILTER		Filters the entries in the list by name and by the unit for which the control module variable and parameter are configured.
			Activates / deactivates the "Parameter" and "Unit" filters.
			Saves the filter settings as a template.
			Selects a template for the filter.
	OK		Confirms the selection.

4.2.5 [EDIT OPERATOR PROMPT PHASE] Input / Selection Screen

The decision criteria for the “Operator prompt phase” are defined in the input and selection screens.

Fig.20: Input / selection screens for creating an “Operator prompt phase” (example)

Field	Description
META DATA	Shows the meta data.
DECISIONS	Shows the decision criteria (max. 3) that are possible in an “Operator prompt phase”.
NOTIFICATION	Displays the title and a description of the notification message.

4.2.5.1 [META DATA] Input / Selection Field

Fig.21: [META DATA] input / selection field

Field	Description
Name*	Name of the phase (max. 48 characters)
Mode	Runs the phase exclusively automatically.
Description	Description for the phase (max. 255 characters).

\* Required information

4.2.5.2 [DECISIONS] Input / Selection Field

DECISIONS*		
Decision 1*	<input checked="" type="checkbox"/>	Yes
Decision 2*	<input checked="" type="checkbox"/>	No
Decision 3	<input type="checkbox"/>	

Fig.22: [DECISIONS] input / selection field

Field	Description
Decision 1-3*	Selects min. 1 but max. 3 decisions for a phase.  When activating a decision by ticking the corresponding check mark, it is necessary to enter a corresponding decision value in the following field.

\* Required information, if a check mark is ticked (at least one “Decision” is required)

4.2.5.3 [NOTIFICATION] Input / Selection Field

NOTIFICATION	
Title*	Activate pO2 control
Text*	Should the pO2 control be activated?

Fig.23: [NOTIFICATION] input / selection field

Field	Description
Title*	Name of the titles (max. 48 characters)
Text*	Description of the notification message (max. 128 characters)

\* Required information

4.2.6 [EDIT PROFILE PHASE] Input / Selection Screen

The setpoint profile of the selected parameter for the “Profile phase” is configured in the input and selection screens.

EDIT PROFILE PHASE

META DATA

Name\*

Temp Profil 1

Mode

Automatic

Semi-Automatic

Description

PROFILE

Parameter\*

TEMP:Setpoint

Interpolation

Linear

Time unit

s

min

h

d

PROFILE POINTS

Time\*

Value\*

1 min

20 °C

2 min

22 °C

3 min

24 °C

SAVE

CANCEL

Fig. 24: Input / selection screens for the setpoint profile of the temperature controller (example)

Field	Description
META DATA	Shows the meta data.
PROFILE POINTS	Displays the setpoint profile.
PROFILE	Selects the control module variable for the setpoint profile.
	Selects the type of interpolation.
	Selects the engineering unit of time.

4.2.6.1 [META DATA] Input / Selection Field

META DATA

Name\*

Mode

Description

Automatic

Semi-Automatic

Fig.25: [META DATA] input / selection field

Field	Description
Name*	Name of the phase (max. 48 characters)
Mode	Selects the phase control mode:
Automatic	Runs the phase automatically.
Semi-Automatic	Runs the phase after confirmation from the user.
Description	Description for the phase (max. 255 characters)

\* Required information

#### 4.2.6.2 [PROFILE POINTS] Input / Selection Field

The setpoint profile that is sent to the device during the phase is defined in the input and selection field. The setpoint profile is processed once the "Profile phase" starts. The setpoint for the control module is set to the configured value at the defined times in the setpoint profile.

PROFILE POINTS		+	🗑️
Time*	Value*		
1 min	20 °C		
2 min	22 °C		
3 min	24 °C		

Fig. 26: Input / selection field for a setpoint profile with three profile points (example)

Field	Symbol	Description
	+	Adds a value pair (profile point) to the setpoint profile.
Time*		Time at which the configured setpoint is set
		Numerical value of the time at which the configured setpoint is set:
		– Max. value for "s" (seconds) = 315360000
		– Max. value for "min" (minutes) = 5256000
		– Max. value for "h" (hours) = 87600
		– Max. value for "d" (days) = 3650
Value*		Setpoint that is set at the defined time

\* Required information

#### 4.2.6.3 [PROFILE] Input / Selection Field



The following settings are defined in the input and selection field:

- Selection of the control module
- Selection of interpolation
- Selection of the engineering unit of time

The screenshot shows a control interface titled 'PROFILE'. It contains three input fields:

- Parameter\***: A text box containing 'TEMP.Setpoint' with a folder icon on the right.
- Interpolation**: A dropdown menu currently showing 'Linear'.
- Time unit**: A selection area with four buttons: 's', 'min' (which is highlighted in black), 'h', and 'd'.

Fig. 27: [PROFILE] input / selection field for the setpoint controller of the "TEMP" control module (example)

Field	Symbol	Description
Parameter*		Shows the [SELECT PARAMETER (#/#)] selection menu.
Interpolation		Selects the type of interpolation:
Linear		<p>If profile points that are at least 20 seconds apart are defined: Every 10 seconds, an interpolated setpoint is automatically inserted.</p> <p>If profile points that are less than 20 seconds apart are defined: <b>No</b> interpolated setpoint is inserted.</p>
Step		<p>If profile points that are at least 20 seconds apart are defined: Every 10 seconds, an interpolated setpoint is automatically inserted. This setpoint corresponds to the value of the preceding profile point.</p> <p>If profile points that are less than 20 seconds apart are defined: <b>No</b> interpolated setpoint is inserted.</p>
Time unit		Selects the engineering unit of time: "seconds", "minutes", "hours", or "days".

\* Required information



## “Linear” interpolation

This example demonstrates how the “TEMP.Setpoint” controller is controlled when set to “Linear” interpolation.

The screenshot shows the 'EDIT PROFILE PHASE' interface. On the left, under 'META DATA', the 'Name' is 'Profile Phase linear', 'Mode' is 'Automatic', and 'Description' is empty. Under 'PROFILE', the 'Parameter' is 'TEMP.Setpoint', 'Interpolation' is 'Linear', and 'Time unit' is 's'. On the right, the 'PROFILE POINTS' table lists the following data:

Time*	Value*
60 s	34 °C
100 s	36 °C
115 s	38 °C
150 s	40 °C

Fig. 28: Input / selection screen with the “Linear” setting (example)

The time curve and control of the setpoint controller are presented in the following table:

Time (mm:ss)	Setpoint (°C)	Description
00:00		The “Profile phase” is started.
<b>1:00 AM</b>	<b>34</b>	The first setpoint of the setpoint profile is set (user configuration).
1:10 AM	34.5	Between 01:00 and 01:40, interpolated setpoints are automatically set in 10-second intervals.
1:20 AM	35	
1:30 AM	35.5	
<b>1:40 AM</b>	<b>36</b>	The second setpoint of the setpoint profile is set (user configuration).
		Between 01:40 and 01:55, no interpolation occurs.
<b>1:55 AM</b>	<b>38</b>	The third setpoint of the setpoint profile is set (user configuration).
2:05 AM	38.8	Between 1:55 AM and 2:30 AM, interpolated setpoints are automatically set in 10-second intervals.
2:15 AM	39.6	

Time (mm:ss)	Setpoint (°C)	Description
2:25 AM		The setpoint interpolation is omitted because the value to the next profile point is less than 10 seconds.
2:30 AM	40	The fourth setpoint of the setpoint profile is set (user configuration). The "Profile phase" is processed.

The following figure shows the time curve and control of the setpoint controller in a graph:




Fig.29: "Linear" interpolation of the setpoints (example)

Symbol	Description
●	Setpoint at time "X" configured by the user.
○	Interpolated setpoint at time "X" set by the system.
○	The setpoint interpolation is omitted because the value to the next profile point is less than 10 seconds.

### “Step” interpolation

This example demonstrates how the “TEMP.Setpoint” controller is controlled when set to “Step” interpolation.

The screenshot shows the 'EDIT PROFILE PHASE' interface. On the left, under 'META DATA', the 'Name' is 'Profile Phase step', 'Mode' is 'Automatic', and 'Description' is empty. Under 'PROFILE', the 'Parameter' is 'TEMP.Setpoint', 'Interpolation' is set to 'Step', and 'Time unit' is 's'. On the right, the 'PROFILE POINTS' table lists four points: 60 s at 34 °C, 100 s at 36 °C, 115 s at 38 °C, and 150 s at 40 °C.

EDIT PROFILE PHASE											
<b>META DATA</b> Name*: Profile Phase step Mode: Automatic Semi-Automatic Description: 											
<b>PROFILE</b> Parameter*: TEMP.Setpoint Interpolation: Step Time unit: s min h d											
<b>PROFILE POINTS</b> <table border="1"> <thead> <tr> <th>Time*</th> <th>Value*</th> </tr> </thead> <tbody> <tr> <td>60 s</td> <td>34 °C</td> </tr> <tr> <td>100 s</td> <td>36 °C</td> </tr> <tr> <td>115 s</td> <td>38 °C</td> </tr> <tr> <td>150 s</td> <td>40 °C</td> </tr> </tbody> </table>		Time*	Value*	60 s	34 °C	100 s	36 °C	115 s	38 °C	150 s	40 °C
Time*	Value*										
60 s	34 °C										
100 s	36 °C										
115 s	38 °C										
150 s	40 °C										

Fig. 30: Input / selection screen with the “Step” setting (example)

The time curve and control of the setpoint controller are presented in the following table:

Time (mm:ss)	Setpoint (°C)	Description
00:00		The “Profile phase” is started.
<b>1:00 AM</b>	<b>34</b>	The first setpoint of the setpoint profile is set (user configuration).
1:10 AM	34	Between 01:00 and 01:40, interpolated setpoints are automatically set in 10-second intervals. The setpoints correspond to the value at 01:00.
1:20 AM	34	
1:30 AM	34	
<b>1:40 AM</b>	<b>36</b>	The second setpoint of the setpoint profile is set (user configuration).
		Between 01:40 and 01:55, no interpolation occurs.
<b>1:55 AM</b>	<b>38</b>	The third setpoint of the setpoint profile is set (user configuration).
2:05 AM	38	Between 1:55 AM and 2:30 AM, interpolated setpoints are automatically set in 10-second intervals. The setpoints correspond to the value of the last profile point (at 01:55).
2:15 AM	38	
2:25 AM		
		The setpoint interpolation is omitted because the value to the next profile point is less than 10 seconds.

Time (mm:ss)	Setpoint (°C)	Description
2:30 AM	40	The fourth setpoint of the setpoint profile is set (user configuration). The "Profile phase" is processed.

The following figure shows the time curve and control of the setpoint controller in a graph:




Fig. 31: "Step" interpolation of the setpoints (example)

Symbol	Description
●	Setpoint at time "X" configured by the user.
○	Interpolated setpoint at time "X" set by the system.
○	The setpoint interpolation is omitted because the value to the next profile point is less than 10 seconds.

#### 4.2.6.4 [SELECT PARAMETER (#/#)] Selection Menu

The selection menu lists the parameters of all units that can be selected for the "Profile phase". The recipe can only be selected for the batch process if same parameters that are configured for the unit are used.









Fig. 32: [SELECT PARAMETER (#/#)] selection menu

Pos.	Field	Symbol	Description
1	Parameter		Name of parameter
		 	Sorts the entries into alphabetical order (ascending and descending)
	Engineering unit		Engineering unit of the parameter (if available)
	Unit		Units for which the parameter is configured
2	FILTER		Filters the entries in the list by name and by the unit for which the control module variable and parameter are configured.
			Activates / deactivates the "Parameter" and "Unit" filters.
			Saves the "Unit" filter setting as a template.
			Selects a template for the "Unit" filter.
	OK		Confirms the selection.

#### 4.2.7 [EDIT INSTRUMENT PHASE] Input / Selection Screen

The parameters for the “Instrument phase” are selected and set in the input and selection screens.

Fig.33: [EDIT INSTRUMENT PHASE] input / selection screen

Field	Description
META DATA	Entering meta data:
Name*	Name of the phase (max. 48 characters)
Mode	Selects the phase control mode:
Automatic	– Runs the phase automatically.
Semi-Automatic	– Runs the phase after confirmation from the user.
Description	Description of the phase
INSTRUMENT	Selecting the sequence:
Phase*	Selects the control module that the sequence displays on the device.
Force restart	Selects the behavior of the sequence at the start of the “Instrument phase”:
Yes	– Restarts the phase on the device if this phase is already being executed.
No	– Does <b>not</b> restart the sequence on the device if this phase is already being executed.

\* Required information

#### 4.2.7.1 [SELECT PHASE] Input / Selection Screen

The selection menu lists the phases of all units that can be selected for the "Instrument phase". The recipe can only be selected for the batch process if same phases that are configured for the unit are used.




Fig. 34: [SELECT PHASE (#/#)] selection menu (example)

Pos.	Field	Symbol	Description
1	Name		Name of the control module
		▲ ▼	Sorts the entries into alphabetical order (ascending and descending)
	Name on device		Name of the sequence on the DCU
	Unit		Units for which the control module is configured.
2	FILTER		Filters the entries in the list by name and by the unit for which the control module is configured.
		✎	Activates / deactivates the filters "Name", "Name on device", "Unit"
		💾	Saves the filter settings as a template.
		★	Selects a template.
	OK		Confirms the selection.

4.2.8 [EDIT WAIT PHASE] Input / Selection Screen

The “Wait phase” can be used to define a wait time before the next phase is run.

EDIT WAIT PHASE?

META DATA

Name\*

Wait phase

Mode

Automatic

Description

TIME

Time unit

s

min

h

d

Duration\*

1 min

SAVE

CANCEL

Fig.35: [EDIT WAIT PHASE] input / selection screen with “minutes” set as the unit of time (example)

Field	Description
META DATA	Entering meta data:
Name*	Name of the phase (max. 48 characters)
Mode	Runs the phase exclusively automatically.
Description	Description of the phase
TIME	Wait time setting:
Time unit	Selects the engineering unit of time: “seconds”, “minutes”, “hours”, or “days”.
Duration*	Enters the numerical value: <ul style="list-style-type: none"><li>– Max. value for “s” (seconds) = 86313600</li><li>– Max. value for “min” (minutes) = 1438560</li><li>– Max. value for “h” (hours) = 23976</li><li>– Max. value for “d” (days) = 999</li></ul>

\* Required information



### 4.3 Formula Editors

**[EDIT CONDITION] formula editor (“Transition” function element)**

The [EDIT CONDITION] input screen is shown when a new condition is created or edited in the “Transition” function element.




Fig.36: [EDIT CONDITION] formula editor

Pos.	Field	Description
1	META DATA	Shows the meta data.
2	FUNCTIONS	Selection of operators
3	CONDITION	Entering and previewing the formula (condition)
4	OTHERS	Selection of operators
5	NUM BLOCK	Selection of numerical values and operators

**[META DATA] input / selection screen (CONDITION)**



Fig.37: [META DATA] input / selection screen

Symbol	Field	Description
◀ ▶		Minimizes or maximizes the [META DATA] screen.
	Name*	Name of the condition (max. 48 characters)
	Description	Description of the formula (max. 255 characters)

\* Required information

[ADD CALCULATION] formula editor (“Parameter phase” phase element)

The [ADD CALCULATION] input screen is shown when a formula is assigned to a process parameter.

The screenshot shows the [ADD CALCULATION] formula editor interface. It is divided into several sections. At the top, there is a header bar with the title 'ADD CALCULATION' and a question mark icon. Below this, the interface is split into two main panes. The left pane is labeled 'META DATA' and contains fields for 'Name\*', 'Description', 'Result limitation', 'Result min.', and 'Result max.'. The right pane is labeled 'FORMULA' and contains a large text area for entering the formula, with a placeholder 'y = ' and a prompt 'Enter formula here'. Below the main panes, there is a 'FUNCTIONS' section with a grid of buttons for various mathematical functions and operators. The grid is organized into four columns: 'VARIABLE', 'FORMULA', 'OTHERS', and 'NUM BLOCK'. The 'VARIABLE' column includes buttons for 'getValue', 'ln', 'sqrt', 'floor', 'abs', 'mean', and 'ceil'. The 'FORMULA' column includes buttons for 'calcCycle', '<', '<=', '>', '>=', 'OR', and 'AND'. The 'OTHERS' column includes buttons for '(', ')', '<', '<=', '>', '>=', 'OR', 'AND', and 'NOT'. The 'NUM BLOCK' column includes buttons for '7', '8', '9', '4', '5', '6', '1', '2', '3', '0', and '+'. At the bottom right of the interface, there are two buttons: 'SAVE' and 'CLOSE'.

Fig.38: [ADD CALCULATION] formula editor

Pos.	Field	Description
1	META DATA	Entering meta data
2	FUNCTIONS	Selection of operators
3	FORMULA	Entering and previewing the formula*
4	OTHERS	Selection of operators
5	NUM BLOCK	Selection of numerical values and operators

\* Available if the Calculation Module is licensed and activated

[META DATA] input / selection screen (CALCULATION)

META DATA

Name\*

Description

Result limitation

Result min.

Result max.

Fig. 39: [META DATA] input / selection screen

Symbol	Field	Description
<div><div></div><div></div></div>		Minimizes or maximizes the [META DATA] screen.
	Name*	Name of the calculation (max. 48 characters)
	Description	Description of the calculation (max. 255 characters)
<div><div></div><div><div></div></div></div>	Result limitation	Activates or deactivates the limitation of the calculation result.
	Result min.	Minimum admissible value of the calculation result: If the value of the calculation result is lower than the defined value, the defined value will be used as the calculation result instead.
	Result max.	Maximum admissible value of the calculation result: If the value of the calculation result is greater than the defined value, the defined value will be used as the calculation result instead.

\* Required information

### Operators in the [FUNCTIONS] selection field (CONDITION, CALCULATION)

FUNCTIONS			
VARIABLE		FORMULA	
getValue		calcCycle	
ln	sqrt	log	exp
floor	abs	mean	ceil

Fig.40: Operators in the [FUNCTIONS] selection field

FUNCTIONS			
getValue		calcCycle	
ln	sqrt	log	exp
floor	abs	mean	ceil
sin	cos	tan	x <sup>y</sup>

Fig.41: Operators in the advanced [FUNCTIONS] selection field

Pos./Field	Description
1	Scroll bar for displaying the complete range of options for selection
VARIABLE	Selects a control module variable / system resource, such as ProcessTime.
FORMULA	Selects a copy of a created calculation*.
sin, cos, tan	Selects a trigonometric function.
x <sup>y</sup> , exp, sqrt	Selects an exponential function / root function.
ln, log10	Selects a logarithmic function.
floor	Selects the "rounding down" (floor) function.
abs	Selects the "absolute value" function.
mean	Selects the "mean" function (parameter: control module).
ceil	Selects the "rounding up" (ceiling) function.
getValue	Selects the "getValue" function.
calcCycle	Selects the "calcCycle" function.

\* Available if the Calculation Module is licensed and activated

[CONDITION] / [FORMULA] preview and input screen





Fig. 42: Preview and input screen (“CONDITION” example)

Pos.	Symbol	Description
1		Entering the condition: The condition must be entered as a formula.
2		Displays the preview of the created condition.
		Indicates a syntax error in the formula.

Operators in the [OTHERS] selection field (CONDITION, CALCULATION)

OTHERS		
(	;	)
<	=	>
≤	≠	≥
OR	AND	NOT

Fig. 43: Operators in the [OTHERS] selection field

Field	Description
( )	Selects the bracket function.
;	Selects the separator for dividing up lists of parameters.
<, >, ≤, ≥	Selects a relational function.
=, ≠	Selects an equals function.
OR, AND, NOT	Selects a Boolean function.

Operators in the [NUM BLOCK] selection field (CONDITION, CALCULATION)

NUM BLOCK			
7	8	9	/
4	5	6	*
1	2	3	-
←	0	,	+

Fig.44: Operators in the [NUM BLOCK] selection field

Field	Description
0 - 9	Selects numerical values.
+, -, *, /	Selects arithmetic functions.
←	Deletes the last entered character.

### 4.3.1 [SELECT VARIABLE] Selection Menu

The selection menu (formula editor with the selection “VARIABLE”) lists the variables of all units that can be selected for the formula. If the parameter phase contains a formula in which variables are configured that are **not** available in the selected unit, the formula cannot be selected before the start of the recipe batch. The entry in the selection list of the recipes is grayed out (see Chapter “4.6.2 Starting an Automated Batch Process”, page 106).

#### Example

Recipe A contains a parameter phase with a one-off calculation. The formula is configured with Variable B. However, **neither** variable B nor the corresponding control module has been created in unit C with which recipe A is to be executed. Recipe A **cannot** be used for a batch with which Unit C was selected.




Fig. 45: “SELECT VARIABLE” selection menu

Pos.	Field	Symbol	Description
1	Variable		Displays the name of the control module variable for all units.
		▲ ▼	Sorts the entries into alphabetical order (ascending and descending)
	Units		Displays the units for which the control module variable is configured.
2	FILTER		Filters the entries in the list by the name of the control module variable and by the unit for which the control module variable is configured.
			Activates / deactivates the “Variable” and “Unit” filters.
			Saves the filter settings as a filter template.
			Selects a filter template function.
	OK		Confirms the selection.

### 4.3.2 [SELECT CALCULATION] Selection Menu

The selection menu (formula editor with the selection “FORMULA”) lists the variables of all units that can be selected for the formula if the calculation module is licensed and activated. If the parameter phase contains a resulting formula (after entering another formula), in which variables are configured, that are **not** available in the selected unit, the formula **cannot** be selected before the start of the recipe batch. The entry in the selection list of the recipes is grayed out (see Chapter “4.6.2 Starting an Automated Batch Process”, page 106).

#### Example

Recipe A contains a parameter phase with a one-off calculation. The resulting formula is configured with a formula which contains a Variable B. However, **neither** variable B nor the corresponding control module has been created in unit C with which recipe A is to be executed. Recipe A **cannot** be used for a batch with which Unit C was selected.




Fig.46: [SELECT CALCULATION] selection menu

Pos.	Field	Symbol	Description
1	Name		Name of the calculation
		▲ ▼	Sorts the entries into alphabetical order (ascending and descending)
	Short Description		Short description of the calculation
	Description		Description of the calculation
2	FILTER		Filters the entries in the list according to “Name”, “Short Description” and “Description”.
		✎	Activates / deactivates the filters “Name”, “Short Description” and “Description”.
		💾	Saves the filter settings as a filter template.
		★	Selects a filter template for the filters.
	OK		Confirms the selection.



### 4.3.3 Functions and Operations

#### Terminology

The term "true" is used as a placeholder for a formula. "1" must be the result of the formula.

The term "false" is used as a placeholder for a formula. "0" must be the result of the formula.

The following table describes the functions and structure of the operations.

Field	Argument	Description	Example	
cos	(d); angle, unit of measurement: radian	Calculates the cosine of value "d".	cos(0.5)	= 0.878
sin	(d); angle, unit of measurement: radian	Calculates the sine of value "d".	sin(0.5)	= 0.479
tan	(d); angle, unit of measurement: radian	Calculates the tangent of value "d".	tan(0.5)	= 0.546
abs	(d)	Calculates the absolute value of value "d".	abs(-0.5)	= 0.5
ceil	(d)	Calculates the rounded up value of value "d".	ceil(4.5) ceil(-1.5) ceil(5)	= 5 = -1 = 5
log10	(d)	Calculates the logarithm to the base 10 of value "d".	log(10)	= 1
exp	(d)	Calculates Euler's number to the power of exponent "d".	exp(1)	= 2.71828183
floor	(d)	Calculates the rounded down value of value "d".	floor(4.5) floor(-1.5) floor(5)	= 4 = -2 = 5
AND	(b1; b2; ... bn)	Checks whether all arguments apply. Non-Boolean values are interpreted as Boolean; 1 = true, any other value = false. The result can be used in further calculations. Calculations: false = 0, true = 1	AND(true; true; false) AND(1; 2; 1) AND(1; 1) * 5	= false = false = 5
NOT	(b)	Inverts the argument. Non-Boolean values are interpreted as Boolean; 1 = true, any other value = false. The result can be used in further calculations. Calculations: false = 0, true = 1	NOT(true) NOT(1) NOT(1) * 5	= false = false = 0
OR	(b1; b2; ... bn)	Checks whether an argument is true. Non-Boolean values are interpreted as Boolean; 1 = true, any other value = false. The result can be used in further calculations. Calculations: false = 0, true = 1	OR(true; true; false) OR(0; 0; 0) OR(0; 1) * 5	= true = false = 5

Field	Argument	Description	Example	
mean	(d;count)	Calculates the mean of the last “d” process values. The “count” value is the number of the last process values used for the mean calculation.	mean(TEMP. Value;5)	= 33.75 °C
ln	(d)	Calculates the natural logarithm of value “d” to the base of Euler's number.	ln(10)	= 2.30258509
x <sup>y</sup>	(x; y)	Calculates “x” to the power of exponent “y”.	power(2; 4)	= 16
sqrt	(d)	Calculates the square root of value “d”.	sqrt(9)	= 3
getValue	(d1; d2; [d3])	<p>Returns the value of the control module (d1) (recursive method of value assignment).            d1: Variable, the value of which should be taken from a specified control module            d2: Standard value (numeric or as a variable) to be returned if d1 <b>does not</b> provide a value            d3 (optional): determines whether only values with the quality “good”, or whether also values with the quality “bad”, should be taken into account.</p> <p>Possible values for d3:            Any (or not specified)            GoodValues</p>	getValue(TEMP. Value; 20)	= 20
calcCycle	(d)	<p>Returns the calculation cycle of a control module in seconds.</p> <p>Returns the calculation cycle of the control module allocated to the calculation if calcCycle is used without parameters.</p>	calcCycle(TEMP. Value)	= 15

## 4.4 Creating a Recipe

The following sections use an example to demonstrate how to configure a recipe for an automated batch process.

The batch process consists of three operations: "Preparation", "Inoculation", and "Batch". The following table presents the sequence of the batch process:

Operation	Phase	Condition	Description
Preparation	Start		The "Preparation" operation starts.
	Set temperature and pH		The temperature controller for the medium is activated with a setpoint of 37°C. The pH controller is activated with a setpoint of 6.8.
	Wait phase		After a wait time of 30 minutes to allow the values to stabilize, the next phase starts.
	End		The "Preparation" operation ends.
Inoculation	Start		The "Inoculation" operation starts.
	Inoculate		The culture is inoculated.
	End		The "Inoculation" operation ends.
Batch	Start		The "Batch" operation starts.
		$pO_2 < 35\%$	The condition is analyzed: The first time the $pO_2$ value drops below the value of 35 %, the next phase starts.
	Wait 1 min		After a wait time of one minute, the $pO_2$ value is queried again.
		$pO_2 > 35\%$	The condition is analyzed: If the $pO_2$ value is over 35% again, the "Wait 1 min" phase starts again.
		$pO_2 < 35\%$	The condition is analyzed: If the $pO_2$ value is under 35%, the "Activate $pO_2$ control" phase starts.

Operation	Phase	Condition	Description
	Activate pO <sub>2</sub> control		The pO <sub>2</sub> controller is activated with a setpoint of 35%.
	End		The "Batch" operation ends.

### Requirements

The device and the unit for the batch process must be created (see BioPAT® MFCS 4 User Manual).

### Procedure

- To create a recipe, carry out the following steps:

Tasks	Chapter, page
Open the recipe editor	4.4.1, 64
Enter meta data	4.4.2, 64
Create operations	4.4.3, 65
Configure phases	4.4.4, 65

## 4.4.1 Opening the Recipe Editor

### Procedure

- In the [ADMINISTRATION] function pane, click on the [RECIPES] button.
- ▷ The [RECIPES (#/#)] overview is shown.
- Click on the [Add] button.
- ▷ The [NEW RECIPE] editor is shown.

## 4.4.2 Entering Meta Data

### Procedure

- Enter the meta data in the [GENERAL], [DETAILS] and [DESCRIPTION] input / selection screens.
- To save the entries: Click on the [APPLY] button.

META DATA	
<b>GENERAL</b>	
Name *	E.coli Batch 1
Short name *	ECB1
Originator	User 1
<b>DETAILS</b>	
Medium	Defined Salt Medium
Strain	BL21
Product	GFP
<b>DESCRIPTION</b>	
Description	Example Recipe

### 4.4.3 Creating Operations



#### Procedure


- ▶ To create the first "Preparation" operation: In the sequence list, click on the [Add] button.
- ▷ The [ADD OPERATION] input screen is shown.
- ▶ Enter the name "Preparation" and a description for the operation.
- ▶ To create the "Inoculation" operation: In the sequence list, click on the [Add] button.
- ▷ The [ADD OPERATION] input screen is shown.
- ▶ Enter the name "Inoculation" and a description of the operation.
- ▶ To create the "Batch" operation: In the sequence list, click on the [Add] button.
- ▷ The [ADD OPERATION] input screen is shown.
- ▶ Enter the name "Batch" and a description for the operation.
- ▶ To save the entries: Click on the [APPLY] button.

### 4.4.4 Configuring Phases


#### "Preparation" operation

No additional conditions are defined in the transitions in the "Preparation" operation. Once a phase is processed, the next phase starts immediately. The "End" phase ends the "Preparation" operation.

#### Procedure



- ▶ Highlight the "Preparation" operation in the sequence list.
- ▶ Drag the "Parameter phase" phase type from the library and position it in the planned sequence.
- ▶ Drag the "Wait phase" phase type from the library and position it in the planned sequence.
- ▶ Drag the "End" phase type from the library and position it in the planned sequence.



- ▶ To configure the “Parameter phase” phase type: In the phase, click on the [Edit] button or double-click on the phase.
- ▷ The [EDIT PARAMETER PHASE] input / selection screen is shown.

- ▶ Enter the name “Set temperature and pH” for the phase in the input / selection screen.
- ▶ To start the phase automatically: Select “Automatic” mode.

- ▶ To select the parameters: In the [PARAMETERS] selection menu, click on the [Add] button.
- ▷ The [ADD PARAMETER (#/#)] selection menu is shown.
- ▶ Select the following control module variables:
  - “TEMP.State”
  - “TEMP.Mode”
  - “TEMP.Setpoint”
  - “pH.State”
  - “pH.Mode”
  - “pH.Setpoint”
- ▶ To confirm the selection: Click on the [OK] button.
- ▷ The control module variables are listed in the [PARAMETERS] selection menu.

Name	Value*
TEMP.State	Select state
TEMP.Mode	Select mode
TEMP.Setpoint	°C
pH.State	Select state
pH.Mode	Select mode
pH.Setpoint	pH

- ▶ Set the controller modes:
  - ▶ To set the “State” controller mode for the temperature controller: Next to the “TEMP.State” entry, click on the [Expand] button and select “remote” mode.
  - ▶ To set the “Mode” controller mode for the temperature controller: Next to the “TEMP.Mode” entry, click on the [Expand] button and select “auto” mode.
  - ▶ To set the “State” controller mode for the pH controller: Next to the “pH.State” entry, click on the [Expand] button and select “remote” mode.
  - ▶ To set the “Mode” controller mode for the pH controller: Next to the “pH.Mode” entry, click on the [Expand] button and select “auto” mode.

- ▶ Define the setpoints:
  - ▶ To define the setpoint for the temperature controller: Click in the input screen next to the "TEMP.Setpoint" entry and enter temperature value "37".
  - ▶ To define the setpoint for the pH controller: Click in the input screen next to the "pH.Setpoint" entry and enter pH value "6.8".
- ▶ To confirm the entries for the phase: Click on the [SAVE] button.

PARAMETERS		+	▲	▼	🗑️
Name	Value *				
TEMP.State	remote				
TEMP.Mode	auto				
TEMP.Setpoint	37 °C				
pH.State	remote				
pH.Mode	auto				
pH.Setpoint	6,8 pH				

SAVE
CANCEL

Start

+

[P]

Set temperature and...

+

⌚

Wait Phase

✕

+

🏠

End

+

- ▶ To configure the "Wait phase" phase type: In the phase, click on the [Edit] button or double-click on the phase.
- ▶ The [EDIT WAIT PHASE] input / selection screen is shown.

EDIT WAIT PHASE

?

META DATA

Name \*

Wait 30min

Mode

Automatic

Description

TIME

Time unit

s min h d

Duration \*

30 min

SAVE
CANCEL

- ▶ Enter the name "Wait 30 min" for the phase in the input / selection screen.
- ▶ To select "min" as the engineering unit: Click on the [min] button.
- ▶ Enter value "30" for the wait time.
- ▶ Click on the [SAVE] button.

- ▶ Link the following phases to each other:
  - "Start"
  - "Set temperature and pH"
  - "Wait 30 min"
  - "End"
- ▷ The phases are linked to each other via the transitions.
- ▶ Click on the [APPLY] button.
- ▷ The "Preparation" operation is configured.
- ▶ Configure the "Inoculation" operation.




### "Inoculation" operation


No additional conditions are defined in the transitions in the "Inoculation" operation. The medium in the culture vessel is inoculated while the operation is run. Once a phase is processed, the next phase starts. The "End" phase ends the "Inoculation" operation.

### Procedure

- ▶ Highlight the "Inoculation" operation in the sequence list.
- ▶ Drag the "Parameter phase" phase type from the library and position it in the planned sequence.
- ▶ Drag the "End" phase type from the library and position it in the planned sequence.








- ▶ To configure the “Parameter phase” phase type: In the phase, click on the [Edit] button or double-click on the phase.
- ▶ The [EDIT PARAMETER PHASE] input / selection screen is shown.

EDIT PARAMETER PHASE	
META DATA	
Name *	Inoculate
Mode	<input type="radio"/> Automatic <input checked="" type="radio"/> Semi-Automatic
Description	

- ▶ Enter the name “Inoculate” for the phase in the input / selection screen.
- ▶ To start the phase manually after inoculation: Select “Semi-Automatic” mode.
- ▶ Click on the [SAVE] button.
- ▶ Link the following phases to each other:
  - “Start”
  - “Inoculate”
  - “End”
- ▶ The phases are linked to each other via the transitions.
- ▶ Click on the [APPLY] button.
- ▶ Configure the “Batch” operation.



### “Batch” operation



In the “Batch” operation, a condition is defined for the  $pO_2$  value in the transition after the “Start” phase. Once the condition is met, the “Wait 1 min” phase starts immediately. The  $pO_2$  value is queried between the Wait phase and the following Parameter phase. This query takes place via two transitions. The following conditions are defined in the transitions:

- Transition with condition A: If the  $pO_2$  value (oxygen content) is below a defined value after the wait time, the “Activate  $pO_2$  control” parameter phase starts immediately.
- Transition with condition B: If the  $pO_2$  value (oxygen content) is above a defined value after the wait time, the “Wait 1 min” phase starts immediately again. The Wait phase is repeated until condition A is met.

Once the “Activate  $pO_2$  control” parameter phase is processed, the “End” phase ends the “Batch” operation.


### Procedure

- ▶ Highlight the “Batch” operation in the sequence list.
- ▶ Drag the “Wait phase” phase type from the library and position it in the planned sequence.
- ▶ Drag the “Parameter phase” phase type from the library and position it in the planned sequence.
- ▶ Drag the “End” phase type from the library and position it in the planned sequence.



- ▶ To configure the “Wait phase” phase type: In the phase, click on the [Edit] button or double-click on the phase.
- ▷ The [EDIT WAIT PHASE] input / selection screen is shown.

- ▶ Enter the name “Wait 1 min” for the phase in the input / selection screen.
- ▶ To select “min” as the engineering unit: Click on the [min] button.
- ▶ Enter value “1” for the wait time.
- ▶ Click on the [SAVE] button.



- ▶ To configure the “Parameter phase” phase type: In the phase, click on the [Edit] button or double-click on the phase.
- ▷ The [EDIT PARAMETER PHASE] input / selection screen is shown.

EDIT PARAMETER PHASE
 

META DATA
 

Name\*
 

Activate pO<sub>2</sub> control

Mode
 

Automatic
 Semi-Automatic

Description

- ▶ Enter the name “Activate pO<sub>2</sub> control” for the phase in the input / selection screen.
- ▶ To start the phase automatically: Select “Automatic” mode.

PARAMETERS
 

Name
 Value\*
 

pO<sub>2</sub>.State
 Select state

pO<sub>2</sub>.Mode
 Select mode

pO<sub>2</sub>.Setpoint
 %sat

SAVE
 CANCEL

- ▶ To select the parameter: In the [PARAMETERS] selection menu, click on the [Add] button.
- ▷ The [ADD PARAMETER (#/#)] selection menu is shown.
- ▶ Select the following control module variables:
  - “pO<sub>2</sub>.State”
  - “pO<sub>2</sub>.Mode”
  - “pO<sub>2</sub>.Setpoint”
- ▶ To confirm the selection: Click on the [OK] button.
- ▷ The control module variables are listed in the [PARAMETERS] selection menu.
- ▶ Set the controller modes:
  - ▶ To set the “State” controller mode for the pO<sub>2</sub> controller: Next to the “pO<sub>2</sub>.State” entry, click on the [Expand] button and select “remote” mode.
  - ▶ To set the “Mode” controller mode for the pO<sub>2</sub> controller: Next to the “pO<sub>2</sub>.Mode” entry, click on the Activate pO<sub>2</sub> control [Expand] button and select “auto” mode.
- ▶ To define the setpoint: Click in the input screen next to the “pO<sub>2</sub>.Setpoint” entry and enter pO<sub>2</sub> value “35”.

PARAMETERS	
Name	Value*
pO2.State	remote
pO2.Mode	auto
pO2.Setpoint	35 %sat

SAVE CANCEL

- ▶ To confirm the entries for the phase: Click on the [SAVE] button.
- ▶ Link the following phases to each other.
  - "Start"
  - "Wait 1 min"
  - "Activate pO2 control"
  - "End"



- ▶ Create the condition for starting the "Wait 1 min" phase: In order to do this, click on the [Edit] button in the transition.
- ▶ The [EDIT CONDITION] input screen is shown.

EDIT CONDITION	
META DATA	
Name*	pO2 < 35%
Description	

- ▶ Enter the name "pO2.< 35%" in the [META DATA] input screen.

VARIABLE

- ▶ In the [FUNCTIONS] input screen, click on the [VARIABLE] button.
- ▷ The [SELECT VARIABLE] selection menu is shown:

SELECT VARIABLE		FILTER	
Variable ▲	Units	Variable	Units
pH.Mode	U3, U4pH		
pH.Output	U3, U4pH		
pH.Setpoint	U3, U4pH		
pH.Value	U3, U4pH		
pO2.Mode	U3, U5pO, U1		
pO2.Output	U3, U5pO, U1		
pO2.Setpoint	U3, U5pO, U1		
pO2.Value	U3, U5pO, U1		
ProcessTime	U3, U5pO, U4pH, U3tp, U1		
R.Value	Shared		
T(n).Value	Shared		
TEMP.Mode	U3, U3tp		
TEMP.Output	U3, U3tp		
TEMP.Setpoint	U3, U3tp		
TEMP.Value	U3, U3tp		
		OK	CLOSE

- ▶ Select the “pO2.Value” process variable and click on the [OK] button.
- ▷ The process variable is entered in the input field and the condition is displayed in the preview:

CONDITION

pO2 Value

pO2.Value

<


- ▶ In the [OTHERS] input screen, click on the [<] button.
- ▷ The operator “<” is entered in the input field for the condition.
- ▶ Enter the numerical value “35” in the input field.
- ▷ The condition is displayed in the preview:

CONDITION


pO2 Value <35

pO2.Value<35

- Click on the [SAVE] button.
- ▷ The condition is defined. The name of the condition is displayed in the planned sequence next to the transition. When the condition is met, the "Wait 1 min" phase starts immediately.



- Create condition A for starting the "Activate pO2 control" phase: In order to do this, click on the [Edit] button in the transition.
- ▷ The [EDIT CONDITION] input screen is shown.



- Enter the name "pO2.< 35%" in the [META DATA] input screen.
- Enter a description (e.g., "Condition A").

EDIT CONDITION	
META DATA	
Name*	pO2 < 35%
Description	Condition A

VARIABLE

- ▶ In the [FUNCTIONS] input screen, click on the [VARIABLE] button.
- ▷ The [SELECT VARIABLE] selection menu is shown:

SELECT VARIABLE		FILTER	
Variable ▲	Units	Variable	
pH.Mode	U3, U4pH	Units	
pH.Output	U3, U4pH		
pH.Setpoint	U3, U4pH		
pH.Value	U3, U4pH		
pO2.Mode	U3, U5pO, U1		
pO2.Output	U3, U5pO, U1		
pO2.Setpoint	U3, U5pO, U1		
pO2.Value	U3, U5pO, U1		
ProcessTime	U3, U5pO, U4pH, U3tp, U1		
R.Value	Shared		
T(n).Value	Shared		
TEMP.Mode	U3, U3tp		
TEMP.Output	U3, U3tp		
TEMP.Setpoint	U3, U3tp		
TEMP.Value	U3, U3tp		
		OK	CLOSE

- ▶ Select the “pO2.Value” process variable and click on the [OK] button.
- ▷ The process variable is entered in the input field and the condition is displayed in the preview:

CONDITION

pO2 Value

pO2.Value

<


- ▶ In the [OTHERS] input screen, click on the [<] button.
- ▷ The operator “<” is entered in the input field for the condition.
- ▶ Enter the numerical value “35” in the input field.
- ▷ The condition is displayed in the preview:

CONDITION


pO2 Value <35

pO2.Value<35


- Click on the [SAVE] button.
- ▷ Condition A is defined. The name of the condition is displayed in the planned sequence next to the transition. When the condition is met, the "Activate pO2 control" phase starts immediately.



- Create the link to condition B for repeating the "Wait 1 min" phase:
- In the "Wait 1 min" phase output, click on the [Link] symbol.





- In the "Wait 1 min" phase input, click on the [Link] symbol.





- ▶ The “Jump” link is established between the input and output of the “Wait 1 min” phase:



- ▶ Create condition B for starting the “Wait 1 min” phase: In order to do this, click on the [Edit] button in the transition.
- ▶ The [EDIT CONDITION] input screen is shown.

EDIT CONDITION	
META DATA	
Name*	pO2 > 35%
Description	Condition B

- ▶ Enter the name “pO2.Value > 35%” in the [META DATA] input screen.
- ▶ Enter a description (e.g., “Condition B”).

VARIABLE

- In the [FUNCTIONS] input screen, click on the [VARIABLE] button.
- ▷ The [SELECT VARIABLE] selection menu is shown:

SELECT VARIABLE	
Variable ▲	Units
pH.Mode	U3, U4pH
pH.Output	U3, U4pH
pH.Setpoint	U3, U4pH
pH.Value	U3, U4pH
pO2.Mode	U3, U5pO, U1
pO2.Output	U3, U5pO, U1
pO2.Setpoint	U3, U5pO, U1
<b>pO2.Value</b>	<b>U3, U5pO, U1</b>
ProcessTime	U3, U5pO, U4pH, U3tp, U1
R.Value	Shared
T(n).Value	Shared
TEMP.Mode	U3, U3tp
TEMP.Output	U3, U3tp
TEMP.Setpoint	U3, U3tp
TEMP.Value	U3, U3tp

- Select the “pO2.Value” process variable and click on the [OK] button.
- ▷ The process variable is entered in the input field and the condition is displayed in the preview:

CONDITION

pO2.Value

pO2.Value

>

- In the [OTHERS] input screen, click on the [>] button.
- ▷ The operator “>” is entered in the input field for the condition.
- Enter the numerical value “35” in the input field.
- ▷ The condition is displayed in the preview:


CONDITION

pO2.Value > 35

pO2.Value>35

- Click on the [SAVE] button.

- ▷ Condition B is defined. The name of the condition is displayed in the planned sequence next to the transition. When the condition B is met, the “Wait 1 min” phase is repeated until condition A is met.



- ▶ To save the recipe: Click on the [SAVE] button.
- ▷ The recipe editor is closed and the [RECIPES (#/#)] menu is shown. The recipe is listed in the menu.

## 4.5 Special Recipe Functions

### 4.5.1 Configuration of the “Operator Prompt Phase”

A minimum of one and maximum of three decisions must be configured for transition conditions in the “Operator prompt phase” phase type. The decisions must correlate to the planned sequence with the configuration of 1 – 3 transitions. Each transition is followed by a decision path.




Fig.47: Extract of the planned sequence “Operator prompt phase” with 2 transitions (example)

Pos.	Description
1	Configured “Operator prompt phase”
2	Transition that is correlated and configured for the phase with the decision “Yes”
3	Transition that is correlated and configured for the phase with the decision “No”

## Install a condition

### Procedure

- ▶ Click on the [Edit] button for the corresponding transition.
- ▶ The [EDIT CONDITION] window is shown.

FUNCTIONS				OTHERS			NUM BLOCK			
VARIABLE				(	;	)	7	8	9	/
getValue	calcCycle			<	=	>	4	5	6	*
ln	sqrt	log	exp	≤	≠	≥	1	2	3	-
floor	abs	mean	ceil	OR	AND	NOT	←	0	.	+

- ▶ Enter the data in the [META DATA] input screens ("Yes" in the example).
- ▶ Click on the [VARIABLE] button.

Variable	Units
Operator prompt phase.No	None
Operator prompt phase.Yes	None

- ▶ The [SELECT VARIABLE] selection menu is shown.
- ▶ Select "Units" filter.
- ▶ Enter "None" in the input field of the marked "Units" filter.
- ▶ The decision variables configured for the "Operator prompt phase" element are listed.
  - ▶ The decision variables are only available in conditions which are directly connected with the output of an operator prompt phase.
  - ▶ These decision variables displayed are variables that must not be assigned to any units.
  - ▶ The decision variables are created by the user themselves and can **only** be selected in the transition conditions that directly follow the "Operator prompt phase" in the planned sequence.

- Highlight the desired decision variable (“Operator prompt phase. Yes” in this example) and confirm with [OK].

FUNCTIONS				OTHERS			NUM BLOCK			
VARIABLE				(	;	)	7	8	9	/
getValue	calcCycle			<	=	>	4	5	6	*
ln	sqrt	log	exp	≤	≠	≥	1	2	3	-
floor	abs	mean	ceil	OR	AND	NOT	←	0	.	+

- Click on the [SAVE] button.
- ▷ The decision variable for the transition is set to “Yes” and the decision pathway is configured.
- To configure the second decision pathway with the transition “No”: Apply the procedure in the same way.

#### 4.5.2 Setting the Process Time Automatically to Zero

The process time of a started batch can be automatically set to 0.000 h. The “Set process time to zero” function is only executed automatically for the batch that is linked with the corresponding recipe.

The execution is determined in the “Parameter phase” phase type with the “ProcessTime” control module.

The function can be added multiple times within an operation and multiple times within a parameter phase of the same recipe.

The function is executed when the phase transitions to the “Running” status and previous parameters of this phase have been processed.

#### Install “Parameter phase” phase elements

##### Procedure

- Drag the “Parameter phase” phase type from the library and position it in the planned sequence.
- In the phase, click on the [Edit] symbol or double-click on the phase.
- ▷ The [EDIT PARAMETER PHASE] input / selection screen is shown.
- Enter the name for the phase in the input / selection screen.
- Select the mode:
  - To start the phase automatically, select the “Automatic” mode.
  - To start the phase manually, select the “Semi-Automatic” mode.
- In the [PARAMETERS] selection menu, click on the [Add] button.
- ▷ The [ADD PARAMETER (#/#)] selection menu is shown.
- Select the “ProcessTime” control module.
- Confirm the selection with [OK].

PARAMETERS		+ ▲ ▼ ☒	
Name	Value *		
ProcessTime	Set process time to zero		

- ▷ The “ProcessTime” control module is entered in the “PARAMETERS” list.
- ▶ Resumes the phase configuration.
  - ▶ If required: Configure additional parameters for this phase (see Chapter “4.2.6.4 [SELECT PARAMETER (#/#)] Selection Menu”, page 49).
- ▶ Click on the [SAVE] button to complete the setup of the phase.

### 4.5.3 Setting the Alarm Parameters

#### Disabling / enabling

The alarming of a previously correspondingly configured control module can be disabled/enabled via a parameter phase. For example, a recipe-controlled batch runs through various process phases in which the alarming for a specific control module should be temporarily enabled or disabled.

#### Alarm limits

The absolute alarm limits and the deadband value can be changed via a parameter phase of a recipe. For example, a recipe-controlled batch runs through various process phases in which different alarm limits should be used for a specific control module.

In the next batch, the alarm limits that were configured (as default) in the [ADMINISTRATION] function pane are used again.

#### 4.5.3.1 Disabling / Enabling the Alarming

#### Requirements

- Alarming is enabled and configured in the [ADMINISTRATION] function pane for the control module (setting: [Absolute]).
- The alarm parameters are enabled and configured in the [ADMINISTRATION] function pane for the control module.

#### Disable alarming

#### Procedure

- ▶ Drag the “Parameter phase” phase type from the library and position it in the planned sequence.
- ▶ In the phase, click on the [Edit] symbol or double-click on the phase.
- ▷ The [EDIT PARAMETER PHASE] input / selection screen is shown.
- ▶ Enter the name for the phase in the input / selection screen.
- ▶ Select the mode:
  - ▶ To start the phase automatically, select the “Automatic” mode.
  - ▶ To start the phase manually, select the “Semi-Automatic” mode.
- ▶ In the [PARAMETERS] selection menu, click on the [Add] button.
- ▷ The [ADD PARAMETER (#/#)] selection menu is shown.

- Select the “STIRR\_1.Alarming” control module variable (example: STIRR\_1).




Fig. 48: Selection screen with selection of the “STIRR\_1.Alarming” parameter

- Disable the alarming: In order to do this, select [Off] for “Select value”.




Fig. 49: Select [Off] for “Select value”

- Click on the [SAVE] button to complete the setup of the phase.
- ▷ The alarming is disabled when the parameter phase becomes active during the recipe-controlled batch.




Fig. 50: Enabled and disabled alarming in the unit display (example)

Pos. Description

- |   |  |
|---|--|
| 1 | The parameter phase is not yet active: The alarming is enabled.  |
| 2 | The parameter phase has become active: The alarming is disabled. |



## Enable disabled alarming

### Procedure

- ▶ Drag the “Parameter phase” phase type from the library and position it in the planned sequence.
- ▶ In the phase, click on the [Edit] symbol or double-click on the phase.
- ▶ The [EDIT PARAMETER PHASE] input / selection screen is shown.
- ▶ Enter the name for the phase in the input / selection screen.
- ▶ Select the mode:
  - ▶ To start the phase automatically, select the “Automatic” mode.
  - ▶ To start the phase manually, select the “Semi-Automatic” mode.
- ▶ In the [PARAMETERS] selection menu, click on the [Add] button.
- ▶ The [ADD PARAMETER (#/#)] selection menu is shown.
- ▶ Select the “STIRR\_1.Alarming” control module variable (example: STIRR\_1).

ADD PARAMETER (10/180)			FILTER	
Parameter ▲	Unit	Engineering unit	Parameter	Unit
STIRR_1.AlarmDeadband	UB2	rpm	stir	
STIRR_1.AlarmHigh	UB2	rpm	stir	
STIRR_1.AlarmHighHigh	UB2	rpm	Unit	UB2
STIRR_1.Alarming	UB2		RESET ALL	
STIRR_1.AlarmLow	UB2	rpm		

Fig. 51: Selection screen with selection of the “STIRR\_1.Alarming” parameter

- ▶ Enable the alarming: In order to do this, select [On] for “Select value”.

PARAMETERS		+	▲	▼	🗑
Name	Value*				
STIRR_1.Alarming ▼	Select value ▼				
	On				
	Off				

Fig. 52: Select [On] for “Select value”

- ▶ To modify the alarm limits: See Chapter 4.5.3.2.
- ▶ Click on the [SAVE] button to complete the setup of the phase.
- ▷ The disabled alarming is enabled when the parameter phase becomes active during the recipe-controlled batch.




Fig. 53: Disabled and enabled alarming in the unit display (example)

Pos.	Description
1	The parameter phase is not yet active: The alarming is disabled.
2	The parameter phase has become active: The alarming is enabled.

#### 4.5.3.2 Modifying Alarm Limits

##### Requirements

- Alarming is enabled and configured in the [ADMINISTRATION] function pane for the control module (setting: [Absolute]).
- The alarm parameters that should be edited in the parameter phase are enabled and configured in the [ADMINISTRATION] function pane for the control module.

##### Procedure

- ▶ Drag the “Parameter phase” phase type from the library and position it in the planned sequence.
- ▶ In the phase, click on the [Edit] symbol or double-click on the phase.
- ▷ The [EDIT PARAMETER PHASE] input / selection screen is shown.
- ▶ Enter the name for the phase in the input / selection screen.
- ▶ Select the mode:
  - ▶ To start the phase automatically, select the “Automatic” mode.
  - ▶ To start the phase manually, select the “Semi-Automatic” mode.
- ▶ In the [PARAMETERS] selection menu, click on the [Add] button.
- ▷ The [ADD PARAMETER (#/#)] selection menu is shown.

► Select the control module variable(s):

ADD PARAMETER (10/180)			FILTER	
Parameter ▲	Unit	Engineering unit	Parameter	Unit
STIRR_1.AlarmDeadband	UB2	rpm	stir	
STIRR_1.AlarmHigh	UB2	rpm		
STIRR_1.AlarmHighHigh	UB2	rpm		
STIRR_1.Alarming	UB2			
STIRR_1.AlarmLow	UB2	rpm		
STIRR_1.AlarmLowLow	UB2	rpm		
STIRR_1.Mode	UB2			
STIRR_1.Output	UB2	%		
STIRR_1.Setpoint	UB2	rpm		
STIRR_1.State	UB2			
			RESET ALL	

Fig. 54: Selection screen with selection of parameter for the alarm limits

► Modify the alarm limits: Enter the values for the corresponding alarm limits.

PARAMETERS		+	▲	▼	🗑
Name		Value*			
STIRR_1.AlarmDeadband [ rpm ]	▼	rpm			
STIRR_1.AlarmHigh [ rpm ]	▼	rpm			
STIRR_1.AlarmHighHigh [ rpm ]	▼	rpm			
STIRR_1.AlarmLow [ rpm ]	▼	rpm			
STIRR_1.AlarmLowLow [ rpm ]	▼	rpm			

Fig. 55: Entering the alarm limits (example)

- Click on the [SAVE] button to complete the setup of the phase.
- ▷ The modified alarm limit(s) is(are) enabled when the parameter phase becomes active during the recipe-controlled batch.




Fig.56: Modified alarm limits in the unit display (example)

#### Pos. Description

- |   |  |
|---|--|
| 1 | The parameter phase is not yet active: Alarm limits before the modification. |
| 2 | The parameter phase has become active: Alarm limits have been modified.      |

#### 4.5.4 Validating the Recipe Syntax

During creation of a recipe, recipe syntax validation can be activated at any time. When validation is activated, errors in all operations and phases of a recipe are marked with a color and, depending on the type of error, explained in a separate error message. The status of the validation button is the same across all operations in a recipe. Errors in a recipe hinder operation on a unit. Possible errors are:

- Missing connection lines between the phase elements
- Missing END phase
- Connection error in a parallel structure
- Unconfigured phases

The processing of a recipe can be resumed despite there being error messages. It is also possible to save a recipe if the syntax is invalid. An activated validation can be deactivated again at any time.

The following procedure shows the validation function as an example.















Fig. 57: Displaying the error in the validation of the recipe syntax (variants)

Pos.	Symbol	Description
1		Syntax error on an operation level
2		Syntax error in the "Parameter phase"
3		Validation button
		The validation button <b>cannot</b> be activated until at least one operation has been created.
		Validation is deactivated.
		Validation is activated.
4		Tool tip with details of the faulty parallel structure
5		Invalid transition connection

Pos.	Symbol	Description
6		Tool tip with details on the syntax error of the phase element
7	 Wait phase	Syntax error in the "Wait phase"
8	 End is required.	Separate error message (example): <ul style="list-style-type: none"> <li>- "End phase" phase element is missing.</li> <li>- Error in parallel pathways.</li> </ul>

#### Procedure

- ▶ Start the recipe editor.
- ▷ The validation button (1) **cannot** be activated until at least one operation is present.
  
- ▶ Create one operation.
- ▷ The validation button (1) can be activated and a tool tip indicates the deactivated validation.
- ▶ Click on the validation button.
- ▷ Validation is activated.
  - ▷ Syntax errors are highlighted at operation level and within operation in the "Parameter phase".
  - ▷ A missing end phase and errors in parallel paths are indicated in separate error messages.
  - ▷ Invalid transition connections are indicated, e.g. simultaneously running parallel pathways.
- ▶ Position the mouse pointer over the red area of the "Wait phase" element.
- ▷ Details about syntax errors for the phase are displayed in a tool tip.
- ▶ Correct the syntax error.
- ▷ The colored area in the "Wait phase" element changes from red to gray.
- ▶ Position the mouse pointer over the red symbol next to the transition.
- ▷ Details about the missing parallel structure are displayed in a tool tip.
- ▶ Correct the defective parallel structure.
- ▷ The red symbol next to the transition is hidden.



### 4.6 Automated Batch Process

The batch process is started in the “MONITORING” function pane (for starting a batch process, see BioPAT® MFCS 4 User Manual). Before a batch process is started, the “START BATCHES” input and selection screen is displayed. The batch process can be linked to a master recipe in the input and selection screen. This link automatically runs the batch process according to the recipe specifications.

[START BATCHES] input and selection screen

1

START BATCHES

BIOSTAT\_Bplus/Unit 1

Batch Name\*

U1.20170728

Description

Recipe

Load recipe

Date

28.07.2017


Time

08:45:06

START

CANCEL

Fig. 58: [START BATCHES] input and selection screen

Pos.	Field	Symbol	Description
1			Displays the selected unit of the batch process.
	Batch Name*		Displays the name of the batch process.
	Description		Enters a description for the batch process.
	Recipe		Displays the selected recipe. The recipe specifications determine the automated sequence of the batch process.
			Opens the selection list with recipes.
	START		Starts the batch process.
	CANCEL		Undoes the selection of the recipe with the batch process and closes the input and selection screen.

\* Required information

Recipe selection list

The recipe for the automated batch process is selected in the selection list.




Fig.59: Selection list with recipes (example)

Pos.	Field	Description
1		Indicates a recipe that can be selected for the batch process.
2		Indicates a recipe that <b>cannot</b> be selected for the batch process.
3		Filters the entries in the list by <ul style="list-style-type: none"> <li>- Designation</li> <li>- Short name</li> <li>- Originator</li> <li>- Medium</li> <li>- Strain</li> <li>- Product</li> </ul>
	Applicable to unit	Filters the entries in the list by "usability":
	Yes	Displays the recipes that can be selected for the batch process.
	No	Displays the recipes that <b>cannot</b> be selected for the batch process.



### 4.6.1 [RECIPES] Menu in the [MONITORING] Function Pane

A started batch process is displayed in the [RECIPES] menu. The function elements displayed in the "OPERATIONS" list and the planned sequence show the current status and progress of the automated batch process. In the planned sequence, the progress of the automated batch process can or – depending on the type of phase and condition – needs to be intervened in.

#### 4.6.1.1 [OPERATIONS] Sequence List










Fig. 60: Sequence list of operations with three operations (example)

Pos.	Field	Description
1		Collapses and expands the sequence list.
2		Displays the name of the operation.
3		Displays the activity status of the operation:
	Completed	The operation is complete. The phases of the operation have been processed.
	Running	The operation is running. The phases of the operation have <b>not</b> yet been fully processed.
	Idle	The operation has <b>not</b> been started yet.
4		Displays the sequence of operations.

### Display of operations

The way the operation is displayed in the sequence list indicates the status of the activity and the selection status of the operation. The following table lists the different ways of displaying:

Display	Description
	<p>The operation has <b>not</b> been started yet.</p> <p>The operation is <b>not</b> selected.</p> <p>The phases of the operation are <b>not</b> displayed in the planned sequence.</p>
	<p>The operation has <b>not</b> been started yet.</p> <p>The operation is selected.</p> <p>The phases of the operation are displayed in the planned sequence.</p>
	<p>The operation is running.</p> <p>The operation is <b>not</b> selected.</p> <p>The phases of the operation are <b>not</b> displayed in the planned sequence.</p>
	<p>The operation is running.</p> <p>The operation is selected.</p> <p>The phases of the operation are displayed in the planned sequence.</p>
	<p>The operation is complete.</p> <p>The operation is <b>not</b> selected.</p> <p>The phases of the operation are <b>not</b> displayed in the planned sequence.</p>
	<p>The operation is complete.</p> <p>The operation is selected.</p> <p>The phases of the operation are displayed in the planned sequence.</p>

4.6.1.2 Planned Sequence









Fig. 61: Planned sequence for an operation with four phases (example)

Pos.	Field	Description
1		Displays the operation:
		Position of the operation in the sequence.
		Name of the operation.
2		Displays the transition:
		Status of the condition: If the condition is met, the phase is repeated ("Jump").
		Name of the condition.
3		Displays the name of the recipe.
4		Displays the unit of the automated batch process.
5		Displays the transition:
		Status of the condition: If the condition is met, the next phase starts.
		Name of the condition.
6		Displays the name and status of the phase:
	Completed	The phase is processed.
	Running	The phase is running.
	Idle	The phase has <b>not</b> been started yet.
	Held	The phase is paused:
		<div><div></div><div>Pausing of the phase can be ended.</div><div></div><div>The paused phase can be stopped.</div></div>

Pos.	Field	Description
	Stopped	The phase is stopped: <ul style="list-style-type: none"> <li>– The phase can be repeated.</li> <li>– The next transition can be activated.</li> </ul>
	Value	Displays the time until the “Wait phase” is processed (time specified as “dd:hh:mm:ss”).





### Display of conditions







The way the condition is displayed in the planned sequence indicates the status of the condition. The following table lists the different ways of displaying:

Display	Description
	The transition condition is <b>not</b> activated yet.
	The transition condition is being analyzed.
	Displayed during analysis of the condition when the mouse pointer moves over the “Transition” symbol. Cancels the analysis and starts the following phase.
	The transition condition has been analyzed. The next phase is started.
	Displayed when a preceding phase is stopped. Activates the transition that follows the stopped phase.

### Display of the “Wait phase”











The way the “Wait phase” is displayed in the planned sequence indicates the status of the phase. The following table lists the different ways of displaying:

Display	Symbol	Description
		The “Wait phase” has <b>not</b> been started yet.
		The “Wait phase” is started. The time until the phase ends is displayed.
		Pauses the phase.
		Stops the phase.

Display	Symbol	Description
		The "Wait phase" is paused. The time until the phase ends is displayed.
		Continues the phase.
		Stops the phase.
		The "Wait phase" is stopped.
		Restarts the phase.
		The "Wait phase" is processed.

### Display of the "Parameter phase"

The way the "Parameter phase" is displayed in the planned sequence indicates the status of the phase. The following table lists the different ways of displaying:




















Display	Symbol	Description
		The "Parameter phase" in "Automatic" mode has <b>not</b> been started yet.
		The "Parameter phase" in "Semi-Automatic" mode has <b>not</b> been started yet.
		The "Parameter phase" in "Semi-Automatic" mode has <b>not</b> been started yet. The phase must be started manually.
		Starts the phase manually.
		The "Parameter phase" in "Automatic" mode is started.
		Pauses the phase.
		Stops the phase.
		The "Parameter phase" in "Semi-Automatic" mode is started.
		Pauses the phase.
		Stops the phase.

Display	Symbol	Description
<div> <div>[p]</div> <div>Parameter phase</div> <div>Held ▶ ■</div> </div>		The "Parameter phase" in "Automatic" is paused. The phase is <b>not</b> complete.
	▶	Continues the phase.
	■	Stops the phase.
<div> <div>[p]</div> <div>Parameter phase</div> <div>Held ▶ ■</div> </div>		The "Parameter phase" in "Semi-Automatic" mode is paused. The phase is <b>not</b> complete.
	▶	Continues the phase.
	■	Stops the phase.
<div> <div>[p]</div> <div>Parameter phase</div> <div>Stopped ↶</div> </div>		The "Parameter phase" in "Automatic" mode is stopped. The phase is <b>not</b> complete.
	↶	Restarts the phase.
<div> <div>[p]</div> <div>Parameter phase</div> <div>Stopped ↶</div> </div>		The "Parameter phase" in "Semi-Automatic" mode is stopped. The phase is <b>not</b> complete.
	↶	Restarts the phase.
<div> <div>[p]</div> <div>Parameter phase</div> <div>Completed</div> </div>		The "Parameter phase" in "Automatic" mode is processed.
<div> <div>[p]</div> <div>Parameter phase</div> <div>Completed</div> </div>		The "Parameter phase" in "Semi-Automatic" mode is processed.

### Display of the "Profile phase"
















The way the "Profile phase" is displayed in the planned sequence indicates the status of the phase. The following table lists the different ways of displaying:

Display	Symbol	Description
<div> <div>⏏</div> <div>Profile phase</div> <div>Idle</div> </div>		The "Profile phase" in "Automatic" mode has <b>not</b> been started yet.
<div> <div>⏏</div> <div>Profile phase</div> <div>Idle</div> </div>		The "Profile phase" in "Semi-Automatic" mode has <b>not</b> been started yet.
<div> <div>⏏</div> <div>Profile phase</div> <div>Idle ▶</div> </div>		The "Profile phase" in "Semi-Automatic" mode has <b>not</b> been started yet. The phase must be started manually.

Display	Symbol	Description
		Starts the phase manually.
		The "Profile phase" in "Automatic" mode is started.
		Pauses the phase.
		Stops the phase.
		The "Profile phase" in "Semi-Automatic" mode is started.
		Pauses the phase.
		Stops the phase.
		The "Profile phase" in "Automatic" mode is paused. The phase is <b>not</b> complete.
		Continues the phase.
		Stops the phase.
		The "Profile phase" in "Semi-Automatic" mode is paused. The phase is <b>not</b> complete.
		Continues the phase.
		Stops the phase.
		The "Profile phase" in "Automatic" mode is stopped. The phase is <b>not</b> complete.
		Restarts the phase.
		The "Profile phase" in "Semi-Automatic" mode is stopped. The phase is <b>not</b> complete.
		Restarts the phase.
		The "Profile phase" in "Automatic" mode is processed.
		The "Profile phase" in "Semi-Automatic" mode is processed.

Display of the “Operator prompt phase”

The way the “Operator prompt phase” is displayed in the planned sequence indicates the status of the phase. The following table lists the different ways of displaying:

Display	Symbol	Description
<div><div></div><div>Operator prompt phase</div><div>Idle</div></div>		The “Operator prompt phase” has <b>not</b> been started yet.
<div><div></div><div>Operator prompt phase</div><div>Running  </div></div>		The “Operator prompt phase” is started.
		Pauses the phase.
		Stops the phase.
<div><div></div><div>Operator prompt phase</div><div>Held  </div></div>		The “Operator prompt phase” is paused. The decisions in the notification window <b>cannot</b> be selected or confirmed.
		Continues the phase.
		Stops the phase.
<div><div></div><div>Operator prompt phase</div><div>Stopped </div></div>		The “Operator prompt phase” is stopped. The “Operator prompt phase” notification is no longer visible in the notification window.
		Restarts the phase. At the same time, the notification window reopens. The “Operator prompt phase” notification is now visible here with the new phase entry time and corresponding process time.
<div><div></div><div>Operator prompt phase</div><div>Completed</div></div>		The “Operator prompt phase” is processed.




Display

Symbol Description




When the phase has reached “running” status, the notification window with “title” [1] and “text” [2] is shown.



Select the “Operator prompt phase” notification. Now min. 1 to max. 3 selection buttons [3] are visible for the corresponding decision. When the audio signal is activated in the settings (administration), a noise sounds and a loudspeaker symbol [4] is shown in the status line of the window.

Pressing a decision button (“YES” or “NO” can be selected in this case), changes the “Operator prompt phase” to “Completed” status. During the subsequent transition, agreement with the decision variables for the phase is checked. In the case of agreement, the following process is carried out.



In the notification window, the selected decision must be confirmed by clicking “OK”.

The window contains detailed information [5] about:

- the selected decision path,
- the user and
- the current date / time of the process.

The selected decision can be described in more detail by means of an optional comment.



Short name of the unit



Priority (low)


















Date and time of the phase entry | process time





Phase name

## Display of the “Instrument phase”

The way the “Instrument phase” is displayed in the planned sequence indicates the status of the phase. The following table lists the different ways of displaying:








Display	Symbol	Description
		The “Instrument phase” in “Automatic” mode has <b>not</b> been started yet.
		The “Instrument phase” in “Semi-Automatic” mode has <b>not</b> been started yet. The phase must be started manually.
		The “Instrument phase” in “Automatic” mode is started.
		The “Instrument phase” in “Semi-Automatic” mode is started.
		The “Instrument phase” in “Automatic” mode is paused. The phase is <b>not</b> complete.
		The “Instrument phase” in “Semi-Automatic” mode is paused. The phase is <b>not</b> complete.
		The “Instrument phase” in “Automatic” mode is stopped. The sequence on the device is stopped.
		The “Instrument phase” in “Semi-Automatic” mode is stopped. The sequence on the device is stopped.

Display	Symbol	Description
 Instrument phase Completed		The "Instrument phase" in "Automatic" mode is processed. The sequence on the device has been ended.
 Instrument phase Completed		The "Instrument phase" in "Semi-Automatic" mode is processed. The sequence on the device has been ended.

### Display after crash of MFCS recipe service

A phase assumes the "Aborted" status and an alarm is generated when a phase is active (e.g. in the "Running" status) and the MFCS recipe service crashes. The "Aborted" status can occur for any phase type. The "Aborted" status is displayed for the affected phases after restarting the client/service.

The data recording of the corresponding batch continued, but all active recipes at the time of the crash were **not** completed. The following table lists the different ways of displaying:

Display	Description
 Wait phase Aborted	The execution of the "Wait phase" was aborted.
 Parameter phase Aborted	The execution of the "Parameter phase" in "Automatic" mode was aborted.
 Parameter phase Aborted	The execution of the "Parameter phase" in "Semi-Automatic" mode was aborted.
 Profile phase Aborted	The execution of the "Profile phase" in "Automatic" mode was aborted.
 Profile phase Aborted	The execution of the "Profile phase" in "Semi-Automatic" mode was aborted.
 Operator prompt phase Aborted	The execution of the "Operator prompt phase" was aborted.
 Instrument phase Aborted	The execution of the "Instrument phase" in "Automatic" mode was aborted.
 Instrument phase Aborted	The execution of the "Instrument phase" in "Semi-Automatic" mode was aborted.

4.6.1.3 [RECIPES] Preview Window









Fig. 62: [RECIPES] preview window (example)

Pos.	Field	Symbol	Description
1			Displays the status of the current operation:
			The condition is being analyzed.
			Multiple conditions are analyzed at the same time.
			The "Wait phase" is being processed.
			The "Wait phase" is paused. Input is required in order to continue the phase.
			The "Wait phase" is stopped. Input is required in order to restart the phase or activate the next condition.
			The "Parameter phase" is being processed.
			The "Parameter phase" in "Semi-Automatic" mode has <b>not</b> been started yet. The phase must be started manually.
			The "Parameter phase" is stopped or paused. Input is required in order to restart the phase, continue the phase or activate the next condition.
			The "Profile phase" is being processed.
			The "Profile phase" in "Semi-Automatic" mode has <b>not</b> been started yet. The phase must be started manually.

Pos.	Field	Symbol	Description
			The "Profile phase" is stopped or paused. Input is required in order to restart the phase, continue the phase or activate the next condition.
			The "Operator prompt phase" is running. To continue editing the recipe, a decision must be selected and confirmed in the notification window.
			The "Operator prompt phase" is paused or stopped. Input is required in order to restart the phase or activate the next condition.
			The "Instrument phase" is being processed.
			The "Instrument phase" is stopped, paused or aborted. Input is required in order to restart the phase, continue the phase or activate the next condition.
			Two or more phases are processed in parallel.
			At least one of the phases running in parallel is stopped or paused. An entry is required to restart or continue a stopped or paused phase running in parallel.
	Completed		The operations of the automated batch process are complete.
2			Displays the unit of the batch process.
3			Displays the abbreviated name of the recipe.
4	No recipe running		The started batch process of the unit is <b>not</b> recipe-controlled.
5	No batch running		<b>No</b> batch process is started for the unit. The unit was added to the unit group after the batch process for the unit group was started.
6			Displays the name of the operation currently being processed.

## 4.6.2 Starting an Automated Batch Process

The recipe can be edited while the batch process is running. Editing the recipe does not affect the recipe specifications for the ongoing batch process.

### Requirements

- The unit/unit group of the batch process has been configured (for configuration of the unit/unit group, see BioPAT® MFCS 4 User Manual).
- The recipe has been created.

### Procedure

- ▶ To show the “MONITORING” function pane: On the start screen, click on the [MONITORING] button.
- ▷ The [DISPLAY UNIT] or [TREND] menu is shown.
- ▶ In the footer, select the unit for which you would like to start a batch process.
- ▶ Click on the [START] button.
- ▷ The [START BATCHES] input screen is shown.



**START BATCHES**

BIOSTAT\_Bplus/Unit 1

**Batch Name\***

**Description**

**Recipe**

1

- ▶ Edit the name of the batch process and the description.
- ▶ Click on the button (1).
- ▷ The recipe selection list opens.

RECIPES (3/3)					
Name	Short name ▲	Originator	Medium	Strain	Product
Recipe 1	Rp.1	User 1	Typ of medium (1)	Strain (1)	Product (1)
Recipe 2	Rp.2	User 2	Typ of medium (2)	Strain (2)	Product (2)
Recipe 3	Rp.3	User 3	Typ of medium (3)	Strain (3)	Product (3)

- Select the recipe:
  - To filter the recipe selection: Use the filter function.
  - The selection list is refreshed once the filter is set.

RECIPES (3/3)					
Name	Short name ▲	Originator	Medium	Strain	Product
Recipe 1	Rp.1	User 1	Typ of medium (1)	Strain (1)	Product

- Select the recipe.

RECIPES (1/3)					
Name	Short name ▲	Originator	Medium	Strain	Product
Recipe 1	Rp.1	User 1	Typ of medium (1)	Strain (1)	Product

- If the recipe selection does **not** need to be filtered: Select the recipe.

RECIPES (3/3)					
Name	Short name ▲	Originator	Medium	Strain	Product
Recipe 1	Rp.1	User 1	Typ of medium (1)	Strain (1)	Product
Recipe 2	Rp.2	User 2	Typ of medium (2)	Strain (2)	Product
Recipe 3	Rp.3	User 3	Typ of medium (3)	Strain (3)	Product

- To confirm the selection: Click on the [OK] button.
- The recipe specifications are linked to the batch process.
- The automated batch process is started and the process sequence of the unit/unit group is recorded.
- The sequence of the automated batch process is shown in the [RECIPES] menu in the sequence list of the operations and in the planned sequence.
- If semi-automatic phases are defined in the automated batch process: Start the corresponding phase manually (see Chapter 4.6.4).
- To pause and then continue an operation: Pause and continue the operation (see Chapter 4.6.5).

### 4.6.3 Disconnecting a Link

#### Requirements

- A recipe has been selected for the automated batch process.
- The batch process has not been started yet.

#### Procedure

- To cancel the recipe selection: Click on the [CANCEL] button.
- The [START BATCHES] input and selection screen closes.

START BATCHES	
BIOSTAT_Bplus/Unit 1	
Batch Name*	U1.20170728
Description	
Recipe	Recipe 1

### 4.6.4 Starting a Phase Manually


A phase which has been configured in the recipe as a semi-automatic phase (mode: "Semi-Automatic") must be started manually during the recipe-controlled batch process. Once the phase is manually started, the operation is processed further.

**Requirements**


- The phase has been configured as a semi-automatic phase.
- The preceding condition has been met.

**Procedure**


- In the "Parameter phase", click on the button (1).



- To start the phase: In the message, click on [YES].



- ▷ The phase is started and processed.
- ▷ The phase ends when the "Completed" status is displayed.





#### 4.6.5 Pausing and Continuing a Phase

An operation can be paused and continued for the following phase types:

- Parameter phase
- Operator prompt phase
- Profile phase
- Wait phase

A phase of the “Instrument phase” type **cannot** be paused manually. The phase can be paused by the system if another sequence is already being carried out.

The following example demonstrates the process of pausing and continuing a phase during the ongoing “Wait 30 min” phase.



##### Requirements

The preceding condition has been met.

##### Pause a phase

##### Procedure

- In the input of the “Wait 30 min” phase, click on button (1).




- To pause the phase: In the message, click on [YES].
- ▷ The phase is paused.


Continue a phase

Procedure


- In the input of the “Wait 30 min” phase, click on button (1).




- To continue the phase: In the message, click on [YES].



- The phase is continued and processed.



- The phase ends when the “Completed” status is displayed.



### 4.6.6 Stopping and Resetting a Phase

A phase can be stopped and reset for all phase types. The following example demonstrates the process of stopping and resetting a phase during the ongoing “Wait 30 min” phase.



#### Requirements

The preceding condition has been met.

#### Stop a phase

##### Procedure

- In the input of the “Wait 30 min” phase, click on button (1).




- To stop the phase: In the message, click on [YES].
- ▷ The phase will be stopped.


## Reset a phase

### Procedure

- In the input of the “Wait 30 min” phase, click on button (1).



- To reset the phase: In the message, click on [YES].
- ▷ The phase is reset and restarts.



## 4.6.7 Stopping a Phase and Activating a Transition

A phase can be stopped for all phase types and the next transition can be activated. Phase processing is aborted when this occurs.

The following example demonstrates the process of stopping a phase and activating the next transition during the ongoing “Wait 30 min” phase.


### Requirements

The preceding condition has been met.


## Stop a phase

### Procedure

- In the input of the “Wait 30 min” phase, click on button (1).





- To stop the phase: In the message, click on [YES].
- ▷ The phase will be stopped.



## Activate a transition

### Procedure

- ▶ In the input of the “Wait 30 min” phase, click on button (1).



- ▶ To activate the transition: In the message, click on [YES].
- ▷ The transition is activated and the condition is analyzed.

## 4.6.8 Ending an Automated Batch Process

### Procedure

- ▶ In the footer, select the unit/unit group for which you would like to stop recording.
- ▶ Click on the [STOP] button.
- ▷ The “Stop Batch Process” confirmation window is shown:
  - ▷ If all operations of the recipe are complete: The confirmation window on the left is shown.
  - ▶ To confirm the message: Click on the [STOP] button.
  - ▷ The batch process and the recording of the process sequence are stopped.



- ▷ If **not** all operations of the recipe are complete: The confirmation window on the left is shown.
- ▶ To confirm the message: Click on the [STOP] button.
- ▷ The batch process and the recording of the process sequence are stopped.
- ▶ To analyze the process sequence: Switch to the “ANALYSIS” function pane (for process sequence analysis, see BioPAT® MFCS 4 User Manual).

## 4.7 Application Example for Resetting a Calculation via a Parameter Phase

In this application example, the resetting of a calculation via a parameter phase will be displayed using the GetValue method.

Objective: A calculated totalizer should be reset to a defined start value during the course of the recipe-controlled batch process (in this example, the start value "0").

### Create a formula and link with control module

FUNCTIONS				OTHERS			NUM BLOCK			
VARIABLE		FORMULA		(	;	)	7	8	9	/
getValue		calcCycle		<	=	>	4	5	6	*
ln	sqrt	log	exp	≤	≠	≥	1	2	3	-
floor	abs	mean	ceil	OR	AND	NOT	←	0	,	+

Fig. 63: Formula editor, example with calculation to totalize the introduced oxygen (calculated "O2\_Totalizer" system variable)

### Pos. Description

1	Formula: Totalizer_O2 = getValue(O2_Totalizer.Value;0) + O2_Flow.Value × (calcCycle() / 60)
2	getValue (O2_Totalizer.Value;0): Returns the value of the "O2_Totalizer.Value" system variable. The value "0" returns the start value "0" to the "O2_Totalizer.Value" system variable if the control module does <b>not</b> supply a value or if the reset is carried out with a correspondingly configured parameter phase during the recipe-controlled batch.
3	O2_Totalizer.Value: The system variable totalizes the air flow rate of the introduced oxygen (ml).
4	O2_Flow.Value: The process variable represents the current air flow rate (ml/min) of oxygen.
5	(calcCycle()/60): Returns the calculation cycle for the "O2_Totalizer.Value" system variable in minutes.

### Requirements

In the [ADMINISTRATION] function pane, a device is configured that provides process values for the "O2\_Flow" control module. A unit with the "O2\_Flow" and "O2\_Totalizer" control modules is created and the control modules are configured:

- "O2\_Flow" control module:
  - Device: Configured device
  - Type: Process variable

- O2\_Totalizer control module:
  - Device: System device
  - Type: Process variable

#### Procedure

- ▶ Create the formula.
- ▶ Link the "O2\_Totalizer" control module with the formula (see BioPAT® MFCS 4 User Manual, Calculation Module).

#### Create a recipe

Once the batch has started, in order optimally to follow the status of the recipe sequence and changes to the process values before and after resetting the calculation, there is a wait phase before the phase which indicates the resetting. The resetting of the calculation is implemented with a parameter phase.

#### Requirements

- The "O2\_Flow" and "O2\_Totalizer" control modules are created.
- The formula "Totalize\_O2" is created.
- The control module "O2\_Totalizer" is linked with the formula "Totalize\_O2".

#### Procedure

- ▶ Open the recipe editor.
- ▶ Enter the meta data (example of the name of the recipe: "Reset Calculation").
- ▶ Create one operation (example of the name: "OP1") and highlight in the sequence list.
- ▶ Position the "Wait phase" type in the planned sequence and configure (example of waiting time: 2 min).
- ▶ Position the "Parameter phase" type from the library in the planned sequence and configure (example of the name: "Reset calc. O2 Totalizer"):
  - ▶ Select the "O2\_Totalizer.Value" system variable in the [ADD PARAMETER (#/#)] selection menu and confirm the selection with [OK].
  - ▶ The system variable is listed in the [PARAMETERS] selection menu.

EDIT PARAMETER PHASE					
<div> <div>META DATA</div> <div> <div>Name*</div> <div>Reset Calc. O2 Totalizer</div> </div> <div> <div>Mode</div> <div>Automatic</div> <div>Semi-Automatic</div> </div> <div> <div>Description</div> <div></div> </div> </div> <div> <div>PARAMETERS</div> <table> <thead> <tr> <th>Name</th> <th>Value*</th> </tr> </thead> <tbody> <tr> <td>O2_Totalizer.Value</td> <td>ml</td> </tr> </tbody> </table> <div> <div>Value</div> <div>Calculation ...</div> <div>Reset</div> </div> </div>		Name	Value*	O2_Totalizer.Value	ml
Name	Value*				
O2_Totalizer.Value	ml				

SAVE

CANCEL

- ▶ Select the reset. To do this, click on the button (1).
- ▶ Click on the [Reset] entry and confirm the configuration with [SAVE].
- ▶ The parameter phase which indicates the resetting of a calculation is configured.
- ▶ Drag the "End" phase type from the library and position it in the planned sequence.
- ▶ Link the "Start", "Wait 2 min", "Reset Calc. O2 Totalizer" and "End" phases to each other.
- ▶ To save the recipe, click on the [SAVE] button.

## Start the recipe-controlled batch

### Requirements

The "Reset Calculation" recipe has been created.

### Procedure

- ▶ To monitor the process sequence: Display the [UNIT DISPLAY] menu in the [MONITORING] function pane.
- ▶ Select the "URC" unit (Unit-ResetCalc) in the footer.
- ▶ Select "O2\_Flow" and "O2\_Totalizer" control modules for the display in the [UNIT DISPLAY].
- ▶ In the [UNIT DISPLAY] menu, click on the [START] button.
- ▷ The [START BATCHES] input screen is shown.
- ▶ In the selection list, select the "Reset Calculation" recipe and confirm the selection with [OK].
- ▶ Click on the [START BATCHES] input screen, click on the [START] button.
- ▷ The batch has started.

Process sequence: Displaying in the [MONITORING] (RECIPE, UNIT DISPLAY) function pane

Process sequence	RECIPE	UNIT DISPLAY
Process value for "O2_Totalizer" before resetting the calculation		
Process value for "O2_Totalizer" immediately after resetting the calculation		

Fig. 64: Example of calculated value for the "O2\_Totalizer" parameter before and after the reset



Pos.	Description
1	Displays the calculated value for the "O2_Totalizer" parameter in the [UNIT DISPLAY] menu.
2	Displays the not yet fully processed parameter phase (status before reset) in the [RECIPES] menu.
3	Displays the reset value for the "O2_Totalizer" parameter in the [UNIT DISPLAY] menu.
4	Displays the fully processed parameter phase (status after reset) in the [RECIPES] menu.

### 4.8 Display of a Phase Configuration

For a recipe-controlled batch that has been started, the configuration of any phase can be displayed in the [RECIPES] menu. The structure of the window that is shown corresponds to the dialog box for phase configuration in the recipe editor. While the window is open, the other application functions **cannot** be accessed.

The window with the phase configuration

- can be shown for any phase type.
- can be shown irrespective of the status of the phase (e.g. "Idle", "Stopped", "Held" or "Completed").
- only opens in the read-only mode. The phase **cannot** be edited.
- can be closed at any time.




Fig. 65: [RECIPES] menu with an overview of a phase configuration (example)


Pos.	Description
1	[RECIPES] menu in the [MONITORING] function pane
2	Completed phase
3	Overview of the completed phase

### 4.8.1 Opening and Closing the Overview

## Requirements

- A recipe-controlled batch is started.
- The [RECIPES] menu is selected in the [MONITORING] function pane.

## Procedure



- ▶ Select the operation and navigate to the phase for which the phase configuration is to be displayed.
  - ▶ Move the mouse over the phase element.
  - ▶ The [Details] button and the tool tip are shown.
  - ▶ Click on the [Details] button.
- 
- ▶ The window with the phase configuration is shown. Other application functions **cannot** be accessed.

[illegible]

- To close the window, click on the [CLOSE] button.
- After the window is closed, the other application functions can be accessed again.

## 4.9 Exporting and Importing a Recipe

A recipe can be exported and imported. The export / import file contains the complete configuration. The recipe configuration includes the meta data, the operations and phases. The export file (import file) can be saved (uploaded) in (from) a target folder located on a local drive, a removable drive or a network drive. The export / import process is visualized.

### Export

- Individual operations or phases of the recipe **cannot** be excluded from the export.
- The export path configured in the system settings is set by default. The default file name of the export file is created from the recipe name and the export date.


### Import

- The export / import file can be used for importing the recipe if
  - the signature of the export / import file has not changed. (Editing the export / import file after exporting renders the signature invalid and the recipe **cannot** be imported).
  - the export / import file was created with a MFCS program version which corresponds to the installed version or an older version. (The MFCS program must be installed in at least the version with which the export / import file was created.)
- Name/abbreviated name, syntax check and invalid references:
  - The name / abbreviated name of the imported recipe must be changed during the import process if a recipe with the same name / abbreviated name is already in the system.
  - Before the import is complete, the recipe can be edited and the syntax check is activated.
  - The recipe can be imported regardless of invalid references (as long as the recipe contains invalid references, e.g. control modules, the recipe cannot be used for recipe-controlled batches).

## 4.9.1 Exporting a Recipe


### Procedure

- ▶ On the start screen, click on the [ADMINISTRATION] button.
- ▶ In the function pane, click on the [RECIPES] button.
- ▶ Select the recipe for export in the [RECIPES (#/#)] overview.
- ▶ Click on the [Export] button.
- ▷ The [EXPORT RECIPE] window is shown.




### Default file location

- ▶ Click on the [EXPORT] button.
- ▷ The progress of the export is displayed and the recipe is saved as an XML file.
- ▶ To directly access the export file and the complete the export process: Click on the [OPEN FOLDER] button.
- ▶ To complete the export process: Click on the [CLOSE] button.



### Individual file location



- ▶ Click on the [Storage location] button.
- ▷ The window for setting the file location is shown.
- ▶ Select the file location and confirm the selection with [Save].
- ▶ To return to the default file location: Click on the [Reset] button.
- ▶ Click on the [EXPORT] button.
- ▷ The progress of the export is displayed and the recipe is saved as an XML file.
- ▶ To directly access the export file and the complete the export process: Click on the [OPEN FOLDER] button.
- ▶ To complete the export process: Click on the [CLOSE] button.




## 4.9.2 Importing a Recipe

### Procedure

- ▶ On the start screen, click on the [ADMINISTRATION] button.
- ▶ In the function pane, click on the [RECIPES] button.
- ▶ Click on the [Import] button.
- ▶ The [IMPORT RECIPE] window is shown.
- ▶ To select the export / import file: Click on the [Storage location] button.





- ▶ Select the export / import file and confirm the selection.



- ▶ The import progress is displayed.

- ▶ Click on the [IMPORT] button.
- ▶ The export / import file are scanned and verified.



- ▶ The recipe is **not** imported if the signature of the export / import file has changed.
- ▶ To select the export / import file without an altered signature: Click on the [BACK] button and repeat the selection process with an export / import file whose signature has not been changed.
- ▶ The recipe is imported if the signature of the export / import file has **not** been changed.
- ▶ Click on the [CONTINUE] button.
- ▶ The recipe is opened in the editor mode and the phase validation is activated.
  - ▶ Operations / phases which are not correctly configured are in the "OPERATIONS" pane and highlighted in color in the planned sequence. However, the recipe is imported.
  - ▶ Imported recipes which already exist with the same name/abbreviated name in the recipe list must have their names changed.
  - ▶ Change the name/abbreviated name of the imported recipe in the [META DATA] menu in the [GENERAL SETTINGS] pane.
- ▶ To complete the recipe import: Click on the [SAVE] button.
- ▶ The recipe is listed in the recipe menu.

## 4.10 Phase Structures

The corresponding phase structures can be created in the planned sequence for different process sequences. The phase structures are presented in examples in the following chapters.

### 4.10.1 Sequential Path

The phase starts with the “Start” phase. The sequence follows the sequential phase path and ends with the “End” phase.




Fig. 66: Phase with a sequential path (example)

### 4.10.2 Alternative Path

The operation starts with the “Start” phase. The path for which the transition condition is first met is processed. The paths for which the conditions are met later are **not** processed. If multiple conditions are met at the same time, the furthest path on the left is processed. The operation ends with the “End” phase.




Fig. 67: Phase with three alternative paths (example)

### 4.10.3 Parallel Paths

The operation starts with the “Start” phase. As soon as the next condition is met, the phases are started at the same time in the first line of the parallel structure. The next phases are processed in parallel. As soon as all phases in the last line of the parallel structure are complete, the next condition is analyzed. As soon as this condition is met, the operation ends with the “End” phase.




Fig. 68: Phase with three parallel paths (example)



#### 4.10.4 Path with a “Jump” Link

The phase is processed and two transitions follow. If the condition for transition 1 is met, the next phase is started. If the condition for transition 2 is met, the same phase is restarted and processed.




Fig. 69: Phase with a “Jump” link (example)

Pos.	Description
1	Transition 1 starts the next phase once the condition is met.
2	Transition 2 repeats and starts the same phase once the condition is met.

## 5 Index

### A

Alarms	
activate function .....	85
deactivate function .....	83
modify limits .....	86
Application examples	
reset a calculation .....	114
Automate a batch process .....	91

### C

Condition formula editor .....	53
--------------------------------	----

### D

Display in the planned sequence	
crashed MFCS recipe service...	103
instrument phase .....	102
operations .....	94
operator prompt phase .....	100
parameter phase .....	97
profile phase .....	98
transition/condition .....	96
wait phase .....	96
Drag and drop function	
change the sequence	
of operations .....	14
position the phases in the	
planned sequence .....	15

### O

One-off calculation .....	36
Operation	
continue .....	109
create .....	65
pause .....	109
Operators .....	61

### P

Parameter phase and one-off	
calculation .....	59, 60
Phase	
configure .....	65
link .....	20
pause and continue .....	109
reset .....	111
start manually .....	108
stop .....	111
Phase library .....	9
Phase structures	
alternative paths .....	123
parallel paths .....	124
path with link .....	125
sequential path .....	122
Planned sequence and library .....	33

### R

Recipe	
edit .....	30
export .....	119, 120
import .....	119, 121
validate .....	88

### S

Selection menu	
parameters .....	36
variable .....	59
Sequence list of operations .....	32
Set process time to zero .....	82
Start an automated batch	
process .....	106

### T

Transition/condition .....	96
----------------------------	----

### V

Validate syntax .....	88
-----------------------	----

Sartorius Stedim Systems GmbH  
Robert-Bosch-Strasse 5-7  
34302 Guxhagen, Germany

Phone: +49 5665 407 0  
[www.sartorius.com](http://www.sartorius.com)

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote the other gender as well.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:

06 | 2020