Installation Instructions | Installationsanleitung

Sartorius YCP04MS | YCP04-HW | VF 988

Pipette Calibration Kit | Pipettenkalibrierset

Contents

English: page 2

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Equipment Supplied

YCP04MS	3
YCP04-HW	4
VF 988	5

Overview of the Equipment

YCP04MS	3
YCP04-HW	4
VF 988	5

Getting Started


Test Conditions	6
-----------------------	---

Installation

YCP04MS: Installing the Moisture Trap	7
YCP04-HW: Installing the Moisture Trap	8
VF 988: Installing the Moisture Trap	9

Disposal	10
----------------	----

Equipment Supplied YCP04MS

The equipment supplied includes the following:

- 1 Moisture trap with cover
- 1 Balance adapter: for models MSx124 to MSx524
- 1 Balance adapter: for models MSx225 and MSx125
- 1 Centering ring for 21 ml vessel
- 1 Reduction fitting: pan for 21 ml vessel
- 3 Sampling vessels 21 ml
- 1 Carrying case

Overview of the Equipment YCP04MS

- 1 Cover
- 2 Moisture trap
- 3 Sampling vessel 21 ml
- 4 Reduction fitting: Pan for 21 ml vessel
- 5 Balance adapter for models MSx124 to MSx524
- 6 Balance adapter for models MSx225 and MSx125
- 7 Centering ring for 21 ml vessel


Equipment Supplied YCP04-HW

The equipment supplied includes the following:


- 1 Moisture trap with cover
- 1 Centering ring for moisture trap
- 1 Balance adapter
- 1 Centering disk
- 1 Reduction fitting
- 3 21 ml sampling vessels
- 1 Carrying case

Overview of the Equipment YCP04-HW

- 1 Cover
- 2 Moisture trap
- 3 Centering ring
- 4 21 ml sampling vessel
- 5 Reduction fitting: Pan for 21 ml vessel
- 6 Centering disk
- 7 Balance adapter

Equipment Supplied VF 988

The equipment supplied includes the following:

- 1 Moisture trap with cover and plug
- 1 Draft shield
- 1 Balance adapter
- 3 2.5 ml sampling vessels

Overview of the Equipment VF 988

- 1 Plug for cover
- 2 Cover
- 3 Draft shield
- 4 Moisture trap
- 5 2.5 ml sampling vessel
- 6 Balance adapter

Getting Started

Test Conditions

- Place of installation of the balance must not be exposed to drafts or vibrations.
- No exposure to direct sunlight
- Temperature: 20–25°C,
Constancy during measurement: $\pm 0.5^\circ\text{C}$.
Make sure that the liquid in the sampling vessel and the test piece are conditioned to the ambient temperature.
- Test liquid: Use only de-aerated, bi-distilled water at room temperature. The water must be changed every hour.
- Pipettes: All parts of the pipette being tested must be at the temperature in the weighing chamber before testing can proceed. (Store the pipette in the stand for at least 2 hours.)

For more detailed information on preparing and taking measurements, see the standards issued in your country for pipette calibration.

The standard generally applied is ISO 8655.

In Germany, this replaces the DIN 12650 standard.


Installation

YCP04MS:

Installing the Moisture Trap

Remove the following parts from the balance:

- Pan support (only for models with a readability of 0.1 mg, MSx124 to MSx524)
- Weighing pan

Insert the following components into the balance in the order given:

- Centering disk or ring for the moisture trap as appropriate for the particular model
- Balance adapter as appropriate for the particular model
- Reduction fitting for sampling vessel
- Sampling vessel (21 ml)
- Moisture trap
- Fill moisture trap with tap water
- Cover of moisture trap

Caution

To prevent damage due to corrosion inside the weighing chamber, be sure to empty the moisture trap when not in use. Leave the moisture trap filled with water only for the time it takes to prepare for and perform pipette calibration.

YCP04-HW:

Installing the Moisture Trap

Remove the following parts from the balance:

- Weighing pan
- Shield disk


Place the following components on the balance in the order given:


- Centering ring for the moisture trap
- Centering ring (required for use with 21 ml sampling vessel or ME balance with a short-design draft shield chamber)
- System adapter
- Reduction fitting for sampling vessel
- Sampling vessel (21 ml)
- Moisture trap
- Fill moisture trap with tap water
- Cover of moisture trap

The centering ring raises the position of the moisture trap by approximately 20 mm. Thus it is possible to place the cover on the moisture trap even when using a 21-ml sampling vessel. When using an ME balance with a short-design draft shield chamber (to preclude the necessity of having to open the upper half of the draft shield door when taking samples), this centering ring must be installed so that the moisture trap can be reached even when using short pipettes.

Caution

To prevent damage due to corrosion inside the weighing chamber, be sure to empty the moisture trap when not in use. Leave the moisture trap filled with water only for the time it takes to prepare for and perform pipette calibration.

VF 988:

Installing the Moisture Trap

Remove the following parts from the balance:

- Draft shield
- Weighing pan

Insert the following components into the balance in the order given:

- Balance adapter
- Sampling vessel (2.5 ml)
- Moisture trap on the shield disk
- Fill moisture trap with tap water
- Draft shield
- Cover for the draft shield
- Plug for the cover

Caution

To prevent damage due to corrosion inside the weighing chamber, be sure to empty the moisture trap when not in use. Leave the moisture trap filled with water only for the time it takes to prepare for and perform pipette calibration.

Disposal

The packaging provided for protection of Sartorius products during shipping is constructed from environmentally friendly materials and should be recycled after successful installation of the equipment. For information on recycling options, including recycling of old weighing equipment, contact your municipal waste disposal center or local recycling depot.

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Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius.com

Copyright by Sartorius, Goettingen,
Germany.

All rights reserved. No part of this
publication may be reprinted or
translated in any form or by any
means without the prior written
permission of Sartorius.

The status of the information,
specifications and illustrations in this
manual is indicated by the date given
below. Sartorius reserves the right
to make changes to the technology,
features, specifications and design of
the equipment without notice.

Status:
April 2015,
Sartorius Lab Instruments GmbH & Co. KG,
Goettingen, Germany