

# BIOSTAT STR<sup>®</sup> Generation 3 Bioreactors and Flexsafe STR<sup>®</sup> Bags

Engineered for Ultimate  
Upstream Performance

## Benefits

- Optimize capital costs
- Reduce operating expenses
- Enhanced safety and optimized space
- Improve productivity and ease of use
- Flexible configurations and simplified maintenance
- Increase speed to clinic



## Product Information

BIOSTAT STR<sup>®</sup> hardware, software, and consumables comprise a state-of-the-art single use biomanufacturing platform that delivers exceptional performance to meet your process development and commercial manufacturing needs. BIOSTAT STR<sup>®</sup> bioreactors and Flexsafe STR<sup>®</sup> bags range in size from 50L to 2000L and are engineered to maximize the flexibility and scalability of your operations. This next generation of our BIOSTAT STR<sup>®</sup> system is engineered for quality, precision and seamless scalability so that you can accelerate your success with a stable process that delivers a safe, reliable biologic medicine for patients around the world.

# Technical Specifications

## Bag Holder with TCU

	BIOSTAT STR® 50	BIOSTAT STR® 200	BIOSTAT STR® 500	BIOSTAT STR® 1000	BIOSTAT STR® 2000
Material	AISI 304L stainless steel				
Dimensions W × D × H	788 × 1061 × 1952 mm 31 × 47.8 × 77 in.	788 × 1291 × 1978 mm 31 × 50.8 × 78 in.	1032 × 1614 × 2312 mm 40.7 × 63.5 × 91 in.	1225 × 1664 × 2692 mm 48.2 × 65.5 × 106 in.	3308 × 2722 × 3894 mm 130.2 × 107.2 × 153.3 in. (incl.platform)
Footprint	0.84 m <sup>2</sup> 9.04 ft <sup>2</sup>	1.02 m <sup>2</sup> 10.98 ft <sup>2</sup>	1.67 m <sup>2</sup> 17.98 ft <sup>2</sup>	2.04 m <sup>2</sup> 21.96 ft <sup>2</sup>	9 m <sup>2</sup> 96.88 ft <sup>2</sup> (incl.platform)
Weight	315 kg 695 lb	410 kg 904 lb	700 kg 1,543 lb	875 kg 1,929 lb	2,045 kg 4,509 lb (incl.platform)
Packaging dimensions	1450 × 1150 × 2200 mm 57.1 × 45.3 × 82.7 in.	1450 × 1150 × 2200 mm 57.1 × 45.3 × 82.7 in.	1800 × 1200 × 2400 mm 70.9 × 47.2 × 94.5 in.	1800 × 1500 × 2600 mm 70.9 × 59 × 102.4 in.	
Installed on skid	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Double wall	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Electro-polished	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Single front door	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Holder for gas filters	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Viewing window	1	1	1	1	1
Lateral window for sensors and ports	2	2	2	2	2
Top drive motor	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

available      - not available      () optional, needs to be ordered separately



## Control Tower

### Specifications for all BIOSTAT STR® models and sizes

Material	AISI 304 stainless steel
Dimensions (W × D × H)	800 × 850 × 1820 mm   31.5 × 33.5 × 71.7 in.
Footprint	0.68 m <sup>2</sup>   7.32 ft <sup>2</sup>
Weight	240 kg   529 lb
Packaging dimension	1000 × 1050 × 2100 mm   39.4 × 41.3 × 82.7 in.
Single version	<input type="checkbox"/>
Installed on skid	<input type="checkbox"/>
Color touch screen	22"
Safety measurement and shut-off	<input type="checkbox"/>
Different user level log in	( <input type="checkbox"/> )
Audit Trail function	( <input type="checkbox"/> )

### Facility and Utility Requirements

	BIOSTAT STR® 50	BIOSTAT STR® 200	BIOSTAT STR® 500	BIOSTAT STR® 1000	BIOSTAT STR® 2000
<b>Power Supply</b>					
Power   Frequency   Consumption					
208 VAC   60 Hz   23 A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
400 VAC   50 Hz   18 A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Gas Supply</b>					
Gas specification according to ISO 8573-1: dry, free of oil and dust					
Compressed air (bar)	2	2	2	4	4
Gas pressure O <sub>2</sub> , N <sub>2</sub> , CO <sub>2</sub> (bar)	2	2	2	4	4
<b>Water Supply</b>					
Cooling water - supply pressure (bar)	1.5	1.5	1.5	1.5	2
Drain for water	25 lpm	25 lpm	25 lpm	25 lpm	40 lpm
Temperature	min. 8°C min. 46°F	min. 8°C min. 46°F	min. 8°C min. 46°F	min. 8°C min. 46°F	min. 8°C min. 46°F
Degree of hardness	max. 12 dH	max. 12 dH	max. 12 dH	max. 12 dH	max. 12 dH
<b>Environmental Requirements</b>					
Ambient temperature	5 - 40°C 41 - 104°F				
Relative humidity range	< 85%				

available      - not available      () optional, needs to be ordered separately

## Process Control

	BIOSTAT STR® 50	BIOSTAT STR® 200	BIOSTAT STR® 500	BIOSTAT STR® 1000	BIOSTAT STR® 2000
<b>Agitation Module</b>					
Max. stirrer speed*	240 rpm	150 rpm	110 rpm	90 rpm*	70 rpm
<b>Temperature Control Unit</b>					
Type	Heating   cooling	Heating   cooling	Heating   cooling	Heating   cooling	Heating   cooling
<b>Heating kW   cooling HP</b>					
Temperature control, double wall	8°C (46°F) above cooling water up to 40°C (104°F)	8°C (46°F) above cooling water up to 40°C (104°F)	8°C (46°F) above cooling water up to 40°C (104°F)	8°C (46°F) above cooling water up to 40°C (104°F)	8°C (46°F) above cooling water up to 40°C (104°F)
Over-temperature protection	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Connection to pressure rated cooling water up to 6 bar	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Aeration Module</b>	4-gas mix (O <sub>2</sub> , N <sub>2</sub> , CO <sub>2</sub> , air) with 3 outlets				
<b>Mass Flow Controllers (MFC)</b>					
■ Number of MFCs	6	6	6	6	6
■ Accuracy of MFC	± 1% full-scale	± 1% full-scale	± 1% full-scale	± 1% full-scale	± 1% full-scale
For sparger line	Ring sparger, micro- or combi-sparger (O <sub>2</sub> , N <sub>2</sub> , CO <sub>2</sub> , air)				
■ Flow rates (lpm)**	(0.025 – 5)	(0.1 – 20)	(0.25 – 50)	(0.5 – 100)	(1.0 – 200)
For overlay line	(CO <sub>2</sub> , air)				
■ Flow rates (lpm)**	(0.025 – 5)	(0.1 – 20)	(0.25 – 50)	(0.5 – 100)	(1.0 – 200)
Softkey to switch air   O <sub>2</sub> between ring and micro-sparger	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Advanced DO controller	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
max. gassing rate	0.1 vvm for micro sparger				
max. total gassing rate	0.2 vvm for all STR sizes				
<b>Pump Module</b>					
Max. 3 internal pumps installed					
WM314 speed-controlled	1	1	1	1	1
WM314 fixed speed	2	2	2	2	2
<b>External Pumps</b>					
Speed-controlled	(up to 5)	(up to 5)	(up to 5)	(up to 5)	(up to 5)

\* Valid for 2 × 3-blade impellers. In case of 1 × 3-blade and 1 × 6-blade impellers, the maximum stirrer speed may be reduced depending on the filling level. At maximum filling level, the stirrer speed is limited to 70 rpm at 1,000 L scale.

\*\* Alternative lower flow ranges are available upon request.

available      – not available      () optional, needs to be ordered separately

## Process Control Continued

	BIOSTAT STR® 50	BIOSTAT STR® 200	BIOSTAT STR® 500	BIOSTAT STR® 1000	BIOSTAT STR® 2000
<b>Sensors &amp; Measurement</b>					
Temperature probe Pt 100	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
■ Measurement range	0 - 150°C 32 - 302°F	0 - 150°C 32 - 302°F	0 - 150°C 32 - 302°F	0 - 150°C 32 - 302°F	0 - 150°C 32 - 302°F
pH, single-use	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
■ Measurement range	6.0 - 8.0	6.0 - 8.0	6.0 - 8.0	6.0 - 8.0	6.0 - 8.0
■ Recalibration function	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
pH, electro-chemical, multi-use and single-use	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
■ Measurement range	4 - 10	4 - 10	4 - 10	4 - 10	4 - 10
DO, single-use	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
■ Measurement range	0 - 110%	0 - 110%	0 - 110%	0 - 110%	0 - 110%
■ Recalibration function	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
DO optical or polarographic, reusable	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
■ Measurement range	0 - 100%	0 - 100%	0 - 100%	0 - 100%	0 - 100%
Load cells	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Balance substrate	(up to 4)	(up to 4)	(up to 4)	(up to 4)	(up to 4)
BioPAT® ViaMass	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
BioPAT® Trace (Glucose   Lactate sensor)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
■ Accuracy	0.1 g/L	0.1 g/L	0.1 g/L	0.1 g/L	0.1 g/L
<b>Accessories</b>					
Ladder   Platform	-	-	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )
Filter line IN	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Filter line OUT	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SU exhaust cooler	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	-
Holder for conventional probes	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )

available

- not available

() optional, needs to be ordered separately

## Communication | Interface of Control Tower

	BIOSTAT STR® 50	BIOSTAT STR® 200	BIOSTAT STR® 500	BIOSTAT STR® 1000	BIOSTAT STR® 2000
OPC UA	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
OPC DA (convertor)	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )
RS232	up to 4	up to 4	4	4	5
Industrial Ethernet	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Analog IN	4	4	4	4	4
Analog OUT	4	4	4	4	4
DeltaV™ connectivity	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )
RemoteIO connectivity (on request)	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )
native DeltaV™ solution (on request)	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )	( <input type="checkbox"/> )

available      - not available      () optional, needs to be ordered separately

## Flexsafe STR® Bags

### Geometrical Data of Flexsafe STR® Bags

Flexsafe STR®	50 L	200 L	500 L	1,000 L	2,000 L
Total volume [L]	68	280	700	1,300	2,800
Max. working volume [L]	50	200	500	1,000	2,000
Min. working volume [L]	12.5	50	125	250	500
Turndown ratio	1: 4	1: 4	1: 4	1: 4	1: 4
Bag diameter, d <sub>1</sub> [mm]	370	585	815	997	1,295
Bag height, h <sub>1</sub> [mm]	666	1,055	1,467	1,800	2,330
Ratio, h <sub>1</sub> /d <sub>1</sub>	1.8	1.8	1.8	1.8	1.8
Liquid height, h <sub>2</sub> [mm]	480	783	1,005	1,360	1,670
Ratio, h <sub>2</sub> /d <sub>1</sub>	1.3	1.34	1.23	1.36	1.29
Impeller diameter, d <sub>2</sub> [mm]	143	225	310	379	492
Ratio, d <sub>2</sub> /d <sub>1</sub>	0.38	0.38	0.38	0.38	0.38
Distance between impellers [mm]	186	300	403	493	640
Volume bottom impeller, fully immersed V <sub>1</sub> [L]	12.5	50	125	250	395
Volume top impeller does not touches surface V <sub>2</sub> [L]	23	95	247	450	953
Volume top impeller, fully immersed V <sub>3</sub> [L]*	31	125	327	597	1,266
Bag packaging dimensions W × D × H (mm   in.)	395 × 395 × 1170   15.5 × 15.5 × 46	395 × 595 × 1440   15.5 × 23.4 × 56.7	764 × 1153 × 1320   30.1 × 45.4 × 52	764 × 1153 × 1440   30.1 × 45.4 × 56.7	995 × 1195 × 1715   39.2 × 47 × 67.5

\* Volume varied depending on the bag configuration (e.g. impeller and sparger configuration).

### Example of Flexsafe STR® 200 L Basic Configuration

Pos.	Designation	Tubing Material	Tubing Termination	Remarks
A	Overlay aeration	Si(Pt)	Opta connector	Without spare port
B	Substrate line 1	Si - TPE	MPC quick coupling	
C	Substrate line 2	Si - TPE	MPC quick coupling	
D	Sparger aeration	Si(Pt)	Opta connector	Without spare port
E	Substrate 3-4	Si - TPE, C-Flex®*	MPC quick coupling	2 lines via Y, dip tube
F	Base addition	Si - TPE, C-Flex®*	Clave connector	
G	Antifoam addition	Si - TPE, C-Flex®*	Clave connector	
H	Gas out (exhaust)	ID, Si(Pt)	Opta connector	2 lines via Y
I	Temperature sensor	Si(Pt)	N/A	Reusable sensor (Pt100)
J   K	DO   pH Sensor	N/A	N/A	Optical sensor
L	Small-volume sampling	Si - TPE, Si(Pt)	Clave connector	
M	Bottom-drain harvest	Si - TPE, C-Flex®*	MPC quick coupling	

\* C-Flex® is a registered trademark of Saint-Gobain Performance Plastics Corporation.

# Ordering Information

Description	Order No.
<b>BIOSTAT STR® Bioreactor System</b>	
System will be configured according to process requirements based on broad choice of configuration options. Please contact your sales representative.	
<b>Flexsafe STR® Single-use Bioreactor Bag</b>	
Standard Design	
Flexsafe STR® 50 L	FRS132668
Flexsafe STR® 200 L	FRS132670
Flexsafe STR® 500 L	FRS124057
Flexsafe STR® 1,000 L	FRS124058
Flexsafe STR® 2,000 L	FRS132367
Customer-specific designs: Please contact your sales representative.	
<b>Bag Tester Fleece</b>	
For 50 L STR®	DZ050L-S2SIT
For 200 L STR®	DZ200L-S2SIT
For 500 L STR®	DZ500L-S2SIT
For 1,000 L STR®	DZ001K-S2SIT
For 2,000 L STR®	DZ002K-S2SIT

Description	Order No.
<b>Accessories</b>	
Filter line IN 50   200 L	DS200L-SBFLI
Filter line OUT 50   200 L	DS200L-SBFLO
Filter line IN 500  1,000 L	DS001K-SBFLI
Filter line OUT 500  1,000 L	DS001K-SBFLO
Filter line IN 2,000 L	DS002K-SBFLI
Filter line OUT 2,000 L	DS002K-SBFLO
Lab Cart	1ZG-0032
Lifting device & adapter (STR® 2,000)	2ZG-0009
BIOSTAT STR® ViaMass Support 60mm	DZ-----AD1
BIOSTAT STR® ViaMass Support 100mm	DZ-----AD2
Autoclave tray for re-usable probes	1ZG-0034
<b>Single-Use Exhaust Cooler</b>	
SU exhaust cooler STR® 50   200	DS200L-EC
SU exhaust cooler STR® 500  1,000	DS001K-EC

## Sales and Service Contacts

For further contacts, visit  
[sartorius.com](http://sartorius.com)

### Germany

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49 551 308 0

### USA

Sartorius Stedim North America Inc.  
565 Johnson Avenue  
Bohemia, NY 11716  
Toll-Free +1 800 368 7178