



**Sartobind<sup>®</sup>**  
Membrane Adsorbers  
for Rapid Purification  
of Proteins



# Sartobind – the pace maker in Membrane Adsorber technology

## Unique microporous structure

Sartobind Membrane Adsorbers display a macroporous structure with a pore size of 0.45 or  $> 3 \mu\text{m}$ . That's orders of magnitudes larger than conventional chromatographic gel matrices. Molecules are transported by convective flow to ligands.

## Characteristics of Membrane Adsorbers (MA)

- : Ready-to-use units
- : Simple handling
- : Pore sizes  $> 3$  and  $0.45 \mu\text{m}$
- : Negligible diffusion limitation
- : Low bed heights with 0.3 up to 16 mm
- : Scalable to process dimension with Sartobind SingleSep disposable capsules and Sartobind System reusable modules
- : Robust high performance separations
- : No bed cracking, channeling, air entrapment
- : Flow rate of ion exchange membranes  $> 80 \text{ ml/min}$  100 kPa (linear flow rate:  $> 4,800 \text{ cm/h}$ )
- : Chemistries: strong and weak ion exchange, coupling, affinity and metal chelate ligands

## Low unspecific adsorption

The basis for all Sartobind membranes is a stabilized reinforced cellulose. It is made from regenerated cellulose and during the production to Sartobind it runs through a number of stabilization and grafting steps until a chromatographic matrix is formed on the cellulose backbone. In principle any ligands known from conventional chromatography can be bound covalently on the matrix.

## Speed up 100 times

In a simple experiment using a Sartobind ion exchange unit with  $5 \text{ cm}^2$  membrane area and a Luer Lock syringe you can achieve a flow rate of about 10 ml per 0.5 second by hand which corresponds to a linear flow rate of more than 14,000 cm/h. Even under these conditions you'll attain capture of your protein. Just try it (sample: cytochrome c, buffer: 10 mM sodium phosphate pH 7.0, unit: Sartobind S 5). In production scale (more than 10,000 liters), a typical speed up factor is about 25 measured in direct comparison to conventional column technology (reference: Walter, J. K. in: *Bioseparation and Bioprocessing, Strategies and Considerations for Advanced Economy in Downstream Processing of Biopharmaceutical Proteins*. G. Subramanian (ed.) Wiley VCH, Vol. II p. 447–460 (1998). Flow rate does not affect break through performance (see fig. 2).




Fig. 2: Break through performance of a Sartobind MultiSep 35 ml cation exchange module at monoclonal antibody purification with various flow rates. For other published data see also: H.L. Knudsen et al., Genentech Inc., *J. Chromatogr. A* 907 (2001) 145–154.

## Sartobind membrane types

- : Sartobind S, Q, C and D ion exchange
- : Sartobind IDA (iminodiacetic acid) metal chelate
- : Sartobind Aldehyde
- : Sartobind Epoxy
- : Sartobind Protein A (recombinant)
- : Other ligands on request

## Sartobind applications

- Purification and Concentration
- : Proteins, viruses, viral particles, monoclonal antibodies, oligonucleotides
- Contaminant removal
- : DNA, endotoxins, viruses, host cell proteins

◀ Fig. 1: Cytochrome c is eluted from a Sartobind S 75 Membrane Adsorber

# ... for robust separations

## Chemical compatibility

The housing of Sartobind SingleSep capsules is polypropylene. The housing of Sartobind MA 15, 75 and 100 is polysulfone which is stable to many standard solvents applied in chromatography. Sartobind ion exchange membranes are compatible with alcohols such as ethanol, isopropanol, glycerol, and denaturing solvents such as 8 M urea and 8 M guanidine HCl and can be cleaned with 1 N sodium hydroxide.

## Constant capacity

The robustness of Sartobind Membrane Adsorbers in ion exchange and affinity chromatography has been tested in consecutive runs of 1000 chromatographic cycles, (Fig. 3). Sartobind can be reused many times, and affinity membranes maintain their binding capacity over many cycles reflecting the chemical stability of the ligands.




Fig. 3: Purification cycles of bovine serum albumin on Sartobind Q 75 strong anion exchanger repeated 1000 times. Flow rate: 120 cm/h, cycle time: 10 min, equilibration buffer: 20 mM phosphate buffer pH 7,0, sample: 5 ml bovine serum diluted 1:20 with equilibration buffer, elution buffer: 20 mM phosphate buffer pH 7,0 + 1 N NaCl, regeneration after each 100 cycles with 1 N NaOH

# Formats from lab to process



Fig. 4: The Demo kit for the demonstration of the high speed and selectivity of Sartobind contains everything for an ion exchange chromatography purification

## Demo kit

A demonstration kit for educational purposes is available to demonstrate a protein purification by hand in two minutes. Everything needed is included in the box: Buffer, colored protein | impurity solution, buffers, syringe and recipe for reuse.



Fig. 5: Sartobind MAs may be used by hand or with a chromatographic system via Luer Lock adapters.

## Laboratory: Sartobind MA units

You may use Sartobind Membrane Adsorbers for any ion exchange or affinity chromatography which require high speed and simple operation. Four sizes of laboratory units can be chosen. Sartobind MA 5 (1 membrane layer) are for principle tests and should be disposed after use. If you have to work with larger mg quantities, please choose Sartobind MA 15 (3 layers), MA 75 (15 layers) or MA 100 (5 layers). These may be reused hundreds of times.



Fig. 6: Sartobind membrane



Fig. 7: Vivapure centrifugal units

## Laboratory: Sartobind A4 sheets

The sheet (624 cm<sup>2</sup>) format allows you to design your own Membrane Adsorber formats. An application for Sartobind A4 sheet is e.g. DNA fragment isolation from agarose gels.

## Centrifugal format

If you don't need flow control but you have to screen a large number of proteins in parallel, use the Vivapure<sup>®</sup> centrifugal units.



# SingleSep™ the tool for contaminant removal




Fig. 12: Sartobind SingleSep Scale up from ml to 10000 liter scale

## High speed polishing

Sartobind SingleSep capsules are designed to remove charged contaminants from therapeutic proteins at accelerated flow rates by ion exchange membrane chromatography. The high throughput is a direct result of negligible mass transfer effects and is made possible by the  $>3 \mu\text{m}$  macroporous membrane with 4 mm (15 layer) bed height.

## Capsule design

The Sartobind SingleSep looks like a standard filter capsule except that the adsorbing membrane is reeled on a core to form a cylinder. The flow is from top into the outside channel and then perpendicular through the membrane layers to the center of the cylinder and leaves the capsule at the outlet.

## Efficient

- Higher throughput (g/h) for trace impurity removal

## Economical

- Saves capital
- No hardware investment & maintenance
- No column packing, testing, regeneration
- No re-use validation
- Less unspecific binding – higher yield
- Less labor
- Buffer consumption may be decreased 95%

## Easy to use

- Disposable
- Simple and fast set up
- Handling like a filter capsule

# A standard for virus purification

## Unique for large proteins, DNA and viruses

Sartobind membranes are unique for purification of large proteins and virus particles. Conventional gel beads possess 95% of the binding capacity within the bead and large particles cannot enter their small pores. In contrast Sartobind membranes display huge capacities for large molecules. E.g. for DNA 10 times more dynamic capacity than conventional beads (reference: Sartobind Application Note DNA removal) and for virus purifications gels cannot compete at all. Sartobind membranes have been successfully applied for the purification of alphaherpesvirus, adenovirus and adenoassociated virus. The Vivapure AdenoPACK adenovirus purification and concentration kits offer researchers a convenient way to clean up within 1–2 hours. The system is scalable to production format.




Fig. 13: Viruses can bind to the internal membrane structure



Fig. 14: Adeno PACK for convenient adenovirus purification and concentration

# Technical Data

For liquid chromatography systems  
or Luer Lock syringe



	MA 5	MA 15	MA 75	MA 100
Application	Screening, disposable	Screening, application optimization, reusable	Down-scale trial, screening, application optimization, reusable	Screening, application optimization, reusable
Membrane material	Stabilized reinforced cellulose			
Adsorption area [cm <sup>2</sup> ]	5	15	75	100
Number of layers	1	3	15	5
Bed height [mm]	0.275	0.8	4.0	1.4
Bed volume* [ml]	0.14	0.41	2.1	2.8
Membrane diameter [mm]	25	25	25	50
Housing material	MBS copolymer	Polysulfone	Polysulfone	Polysulfone
Inlet connector	Female Luer Lock	Female Luer Lock	Female Luer Lock	Female Luer Lock
Outlet connector	Male Luer Lock	Male Luer Lock	Male Luer Lock	Male Luer Lock
Minimum static binding capacity* [mg/unit]	4 / S 5 4 / Q 5 3 / D 5	12 / S 15 12 / Q 15 9 / D 15	60 / S 75 60 / Q 75 45 / C 75 45 / D 75 7.5 / IDA 75 2.2 / Epoxy 75 6 / Protein A 75	80 / S 100 80 / Q 100 60 / D 100
Flow rate** at 0.1 MPa (1 bar   14.5 psi) [ml/min]	> 150	> 50	> 20	> 75
Dead volume [ml]	0.8	1.0	1.3	4.2
Maximum pressure [MPa]	0.4	0.6	0.6	0.6
Storage before use at	Room temperature	Room temperature	Room temperature Protein A 75 at +4°C	Room temperature

1 ml membrane volume is equal to 36.4 cm<sup>2</sup> membrane area

\* Reference proteins ion exchange: S, C: lysozyme, Q, D: BSA, IDA: cytochrome C, epoxy. IgG, Protein A: IgG

\*\* Ion exchange

## Overview Membrane Types and Formats

Membrane Description Type	Pore size (µm)	Capacity* (mg/ml)	Sartobind Lab Scale		Sartobind Process Scale			Sartobind Direct
			A4 sheet	MA 5, 15, 100	MA 75	SingleSep	System	
S Sulfonic acid	>3	29	•	•	•	•	•	(•)
Q Quaternary ammonium	>3	29	•	•	•	•	•	(•)
C Carboxylic acid	>3	22	•		•		•	
D Diethylamine	>3	22	•	•	•		•	
IDA Iminodiacetic acid	>3	3.6	•		•		(•)	(•)
Epoxy Epoxy	0.45	1.1	•		•		(•)	
Aldehyde Aldehyde	0.45	1	•					
Protein A Recombinant Protein A	0.45	3			•		(•)	(•)

\* Minimum static binding capacity for reference protein, • available, (•) Sartobind Special

**For centrifuge**



**96 IEX**

**8 Strip**

**Mini**

**Maxi**

	<b>96 IEX</b>	<b>8 Strip</b>	<b>Mini</b>	<b>Maxi</b>
Application	Screening, disposable	Screening, disposable	Screening, application optimization, disposable	Screening, application optimization, disposable
Membrane material	Stabilized reinforced cellulose			
Membrane diameter [mm]	5	5.5	6.4	22.5
Housing material	Polypropylene			
Minimum static binding capacity* [mg/unit]	0.2	IEX: 1 IDA: 130	M: 1 H: 4	M: 12–16 H: 48–64
Maximum volume [ml] (device inner volume)	0.8	0.3	M: 0.5 H: 0.4	M: 20 H: 19
Spin speed [g]	1000	IEX: 1000 IDA: 500	M: 500 H: 2000	M: 500 H: 500
Storage before use at	Room temperature	Room temperature	Room temperature Protein A at +4°C	Room temperature

IEX: ion exchange (S, Q, C, D), IDA: iminodiacetic acid

M: medium, H: high

\* Reference proteins ion exchange: S, C: lysozyme, Q, D: BSA, IDA: cytochrome C

# Ordering Information

## Luer Lock units | sheets | kits

Order No.	Sartobind MA	Quantity
S5F	Sartobind S 5	15
Q5F	Sartobind Q 5	15
D5F	Sartobind D 5	15
S15X	Sartobind S 15	2
Q15X	Sartobind Q 15	2
D15X	Sartobind D 15	2
S75X	Sartobind S 75,	2
Q75X	Sartobind Q 75,	2
C75X	Sartobind C 75,	2
D75X	Sartobind D 75,	2
S100X	Sartobind S 100	1
Q100X	Sartobind Q 100	1
D100X	Sartobind S 100	1
93IDA426DB-12--V	Sartobind IDA 75	2
93EPOX06DB-12--V	Sartobind Epoxy 75	2
93PR-A06DB-12--V	Sartobind Protein A 75	2

Order No.	Sartobind A4 Sheet	Quantity
94IEXS42-001	Sartobind S A4 Sheet	1
94IEXQ42-001	Sartobind Q A4 Sheet	1
94IEXC42-001	Sartobind C A4 Sheet	1
94IEXD42-001	Sartobind D A4 Sheet	1
94IDA-42-001	Sartobind IDA A4 Sheet	1
94EPOX06-001	Sartobind Epoxy A4 Sheet	1
94ALD-06-001	Sartobind Aldehyde A4 Sheet	1

Order No.	Sartobind SingleSep	Quantity
92IEXQ42DN-11	Sartobind Q SingleSep nano 2.5 ml	1
92IEXQ42DN-11--A	Sartobind Q SingleSep nano 2.5 ml	4
92IEXQ42D4-SS--A	Sartobind Q SingleSep mini 7 ml	4
92IEXS42D4-SS--A	Sartobind S SingleSep mini 7 ml	4

Order No.	Accessories	Quantity
VS-AVPQ102	Vivapure Adeno PACK 100 RT*	
VS-AVPQ501	Vivapure Adeno PACK 500	
90-KIT-01	Sartobind Demo Kit	
17002---140	Pair of Luer Lock adapters, black Tefzel to connect Sartobind MA units to a liquid chromatography system with M6 thread	2
16517-----E	Syringe filter holders, 25 mm diameter, polycarbonate for filter or adsorber membrane	12
16214	Syringe filter holder 25 mm diameter, stainless steel for filter or adsorber membrane	1

### Related products for production

Sartobind SingleSep capsules  
Sartobind System modules  
Sartobind Direct

\* Vivapure Adeno PACK 100 RT does not contain Benzonase®

**96 well plates | 8 strips | spin columns**

<b>Order No.</b>	<b>Vivawell</b>	<b>Quantity</b>
VW96IS02	Vivawell-96 S-IEX purification plate	2 plates
VW96IQ02	Vivawell-96 Q-IEX purification plate	2 plates
VW96IC02	Vivawell-96 C-IEX purification plate	2 plates
VW96ID02	Vivawell-96 D-IEX purification plate	2 plates
VW08IS02	Vivawell 8 Strip S-IEX centrifugal purification strips	24 strips (2 plates)
VW08IQ02	Vivawell 8 Strip Q-IEX centrifugal purification strips	24 strips (2 plates)
VW08IC02	Vivawell 8 Strip C-IEX centrifugal purification strips	24 strips (2 plates)
VW08ID02	Vivawell 8 Strip D-IEX centrifugal purification strips	24 strips (2 plates)

<b>Order No.</b>	<b>Vivapure Mini Spin Columns</b>	<b>Spin columns/ Centrifuge tubes</b>
VS-IX01SM24	Vivapure S Mini M	24/48
VS-IX01SH24	Vivapure S Mini H	24/48
VS-IX01QM24	Vivapure Q Mini M	24/48
VS-IX01QH24	Vivapure Q Mini H	24/48
VS-IX01CM24	Vivapure C Mini M	24/48
VS-IX01CH24	Vivapure C Mini H	24/48
VS-IX01DM24	Vivapure D Mini M	24/48
VS-IX01DH24	Vivapure D Mini H	24/48

<b>Order No.</b>	<b>Vivapure Maxi Spin Columns</b>	<b>Spin columns/ Centrifuge tubes</b>
VS-IX20SM08	Vivapure S Maxi M	8/16
VS-IX20SH08	Vivapure S Maxi H	8/16
VS-IX20QM08	Vivapure Q Maxi M	8/16
VS-IX20QH08	Vivapure Q Maxi H	8/16
VS-IX20CM08	Vivapure C Maxi M	8/16
VS-IX20CH08	Vivapure C Maxi H	8/16
VS-IX20DM08	Vivapure D Maxi M	8/16
VS-IX20DH08	Vivapure D Maxi H	8/16

# Sales and Service Contacts

For further contacts, visit [www.sartorius.com](http://www.sartorius.com)

## Europe

### Germany

Sartorius AG  
Weender Landstrasse 94-108  
37075 Goettingen

Phone +49.551.308.0  
Fax +49.551.308.3289

[www.sartorius.com](http://www.sartorius.com)

Sartorius BBI Systems GmbH  
Schwarzenberger Weg 73-79  
34212 Melsungen

Phone +49.5661.71.3400  
Fax +49.5661.71.3702

[www.sartorius.com](http://www.sartorius.com)

### U.K.

Sartorius Ltd.  
Longmead Business Park  
Blenheim Road, Epsom  
Surrey KT19 9 QQ

Phone +44.1372.737100  
Fax +44.1372.720799

### Austria

Sartorius Ges.m.b.H. Wien  
Franzosengraben 12  
A-1030 Vienna

Phone +43.1.7965763.18  
Fax +43.1.796576344

### Belgium

Sartorius Technologies N.V.  
Leuvensesteenweg, 248/B  
1800 Vilvoorde

Phone +32.2.756.06.80  
Fax +32.2.756.06.81

### Denmark

Sartorius Nordic A/S  
Himmelev Bygade 49  
4000 Roskilde

Phone +45.70.23.4400  
Fax +45.46.30.4030

### France

Sartorius S.A.S.  
4, rue Emile Baudot  
91127 Palaiseau Cedex

Phone +33.1.6919.2100  
Fax +33.1.6920.0922

### Italy

Sartorius S.p.A.  
Via dell'Antella, 76/A  
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41  
Fax +39.055.63.40.526

### Ireland

Sartorius Limited  
Unit 41, The Business Centre  
Stadium Business Park  
Ballycoolin Road  
Dublin 11

Phone +353.1.8089050  
Fax +353.1.8089388

### Netherlands

Sartorius Technologies B.V.  
Edisonbaan 24  
3439 MN Nieuwegein

Phone +31.30.6025080  
Fax +31.30.6025099

### Spain

Sartorius, S.A.  
C/Isabel Colbrand 10 -12,  
Planta 4, Oficina 121  
Polígono Industrial de Fuencarral  
28050 Madrid

Phone +34.91.3586102  
Fax +34.91.3588804

### Switzerland

Sartorius Schweiz AG  
Lerzenstrasse 21  
8953 Dietikon

Phone +41.1.746.50.00  
Fax +41.1.746.50.50

## America

### USA

Sartorius North America, Inc.  
131 Heartland Blvd.  
Edgewood, New York 11717

Phone +1.631.254.4249  
Toll-Free +1.800.3687178  
Fax +1.631.254.4253

Sartorius BBI Systems, Inc.  
2800 Baglyos Circle  
Bethlehem, PA 18020

Phone +1.610.866.4800  
Fax +1.610.866.4890

### Argentina

Sartorius Argentina S.A.  
Calle Avalos 4251 (B1605ECS) Munro  
Buenos Aires

Phone +54.11.4721.0506  
Fax +54.11.4762.2333

### Brazil

Sartorius do Brasil Ltda.  
Rua Santo André, 331  
09020-230 Santo André  
São Paulo

Phone +55.11.4438.3833  
Fax +55.11.4438.2355

### Mexico

Sartorius de México, S.A. de C.V.  
Circuito Circunvalación Poniente No. 149  
Ciudad Satélite  
53100 Naucalpan, Estado de México

Phone +52.55.62.1102  
Fax +52.55.62.2942

## Asia | Pacific

### China

Beijing Sartorius Instrument & System  
Engineering Co., Ltd.

- Beijing Rep. Office -  
Dong Hu Qu, Wang Jing  
Industrial Zone  
Chao Yang District  
100102 Beijing, P.R.C.  
P.O. Box 8516

Phone +86.10.6439.2552  
Fax +86.10.6439.2726

Sartorius Ltd.

Unit 1110-12, Lu Plaza,  
2 Wing Yip Street  
Kwun Tong, Kowloon, Hong Kong

Phone +852.2774.2678  
Fax +852.2766.3526

### India

Sartorius India Private Ltd.  
10, 6th Main, 3rd Phase Peenya  
KIADB Industrial Area  
Bangalore - 560 058

Phone +91.80.2839.1963 | 0461  
Fax +91.80.2839.8262

### Japan

Sartorius K.K.  
KY Building, 8-11  
Kita Shinagawa 1-chome  
Shinagawa-ku  
Tokyo 140-0001

Phone +81.3.3740.5407  
Fax +81.3.3740.5406

### South Korea

Sartorius Korea Biotech  
B-1023, Paragon  
17-2, Jungja-Dong, Bundang-Gu  
Sungnam, Gyunggi-Do  
463-811, South Korea

Phone +82.31.782.7011  
Fax +82.31.782.7666

### Malaysia

Sartorius (Malaysia) Sdn. Bhd.  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur

Phone +60.3.8996.0622  
Fax +60.3.8996.0755

### Singapore

Sartorius Singapore Pte. Ltd.  
10, Science Park Road, The Alpha  
#02-25, Singapore Science Park 2  
Singapore 117684

Phone +65.6872.3966  
Fax +65.6778.2494

### Australia

Sartorius Australia Pty. Ltd.  
Unit 17/104 Ferntree Gully Road  
Waverley Business Park  
East Oakleigh, Victoria 3166

Phone +61.3.9590.8800  
Fax +61.3.9590.8828