

## Sartoflow Alpha

### Don't restrict your choice

Designed for use in laboratory and small scale pilot plants, this system is ideal for scale-up trials. Reliable operation is ensured by use of a rotary lobe pump together with zero dead leg valves and inline pressure transducers.

- Compact (60x50x60 cm)
- Low hold up volume (< 300 ml)
- Ready to use suitable for volumes from 1 to 50l
- Different tank sizes available
- Additional Instrumentation optional
- Data Recording available (21 CFR Part 11)
- Material certificates (3.1b, FDA)
- IQ/OQ documentation available

