

Flexboy® Bags: C-Flex

Standard Flexboys® 5L-50L are available with C-Flex tubing (see fig 4). Aseptic connections with C-Flex tubing are made using a Tube Fuser. Design description is as follows:



Item	Component
1	Bar Handle (red)
2	S71 Film
3	EVA Tubing 4.8 x 6.5 L100mm (3/16" x 1/4" L3.94")
4	Female Luer Lock connector
5	EVA Tubing 9.3 x 12 L100mm (11/30" X 9/19" L3.94")
6	Slide Clamp
7	Connector 9.53mm Luer (3/8") (PC [†])
8	C-Flex Tubing 9.3 x 12 L100mm (11/30" x 9/19" L3.94")
9	Male Luer Lock Cap
10	Male Luer Lock cap with built-in latex free septum

If a totally different design is required, please visit customisation.

Flexboy® Bag Dimensions ***

Volume	5L	10L	20L	50L
Length (L) mm	376 (14.80")	621 (24.45")	654 (25.75")	790 (31.10")
Width (W) mm	332 (13.07")	300 (11.81")	431 (16.97")	580 (22.83")
Length inc. Tubing (T) mm	473 (18.62")	718 (28.27")	749 (29.57")	887 (34.92")
Film Surface Area cm ²	1929 (298,99 sq in)	3528 (546,84 sq in)	4826 (748,03 sq in)	8106 (1256,43 sq in)

[†]PC=Polycarbonate, ***Due to flexible nature of bags, dimensions are approximate values and should not be considered as precise specifications.