

### Stackable Palletank® Range

The Stackable Palletank® System (STS) was developed for customers who wish to leave the containers in close proximity to process equipment. STS Palletanks® therefore allow bag positioning and manipulation whilst stacked in a process zone using securely latched, removable, panels. The panels facilitate in-situ bag positioning and manipulation enabling continuous processing without movement of the containers. As a result a forklift is not routinely required in the manufacturing area. Sterile connections can be made as required using the aseptic connectors, portable LAF cabinets or tube fusers.

STS Palletanks® are available to order in 200L, 500L, and 1000L sizes. The non-stackable and non-foldable Palletank® range is available for process volumes of 1500L and above. STS Palletanks® can be customized, and are delivered ready assembled with the following features:

- Stacking corners
- Integrated Pallet-base and trays
- Bottom gate(s) – one gate on 200-500L versions, and two gates for 1000L
- Specially designed Dolly
- Front panel with safety latches for bag positioning and manipulation



Stackable Palletank® with Removable Front Panels