

S71 Multilayer Mono-Chamber Bag

- The Sartorius Stedim Biotech multi-layer bag is made of S71 dry-laminated multi-layer film.
- The fluid contact layer remains EVAM[®] (EVA).
- The S71 bag offers extremely low oxygen permeability, due to its high gas barrier properties.
- Oxygen permeability of the S71 bag is approximately 70 times lower than standard EVA bags.
- These high gas barrier properties are important for hospital pharmacies and compounding centers who produce bags for intermediate storage and where longer shelf-life is necessary.

