



### Gotove puferne otopine

#### Otopine pufera sa tri različita pH:

PY-Y21	4.00 ± 0.01 pri 25°C, 500 ml
PY-Y21-6	4.00 ± 0.01 pri 25°C, 6 boca, svaka od 90 ml
PY-Y22	7.00 ± 0.01 pri 25°C, 500 ml
PY-Y22-6	7.00 ± 0.01 pri 25°C, 6 boca, svaka od 90 ml
PY-Y23	10.00 ± 0.01 pri 25°C, 500 ml

Svaka otopina je druge boje, kako bi se spriječile zamjene.

Poseban, dvostruki vrat boce nudi prednosti prilikom rada u laboratoriju i izvan njega.

#### pH puferne kapsule

50 kapsula u pakiranju; svaku kapsulu potrebno je otopiti u 100 ml destilirane vode

pH = 4.01± 0.02 pri 25°C	PY-Y01
pH = 7.00± 0.02 pri 25°C	PY-Y02
pH = 9.00± 0.02 pri 25°C	PY-Y03
pH = 10.00± 0.02 pri 25°C	PY-Y04

#### Standardi konduktiviteta, slijedivi do NIST Standarda

0.084 mS/cm ± 1.0% pri 25°C, 500 ml	PY-Y10
0.147 mS/cm ± 1.0% pri 25°C, 500 ml	PY-Y11
1.413 mS/cm ± 1.0% pri 25°C, 500 ml	PY-Y12
12.88 mS/cm ± 1.0% pri 25°C, 500 ml	PY-Y13
Otopina elektrolita, KCl (3 mol/l), 500 ml	PY-Y07
Otopina za čuvanje, za pH kombinirane elektrode, 500 ml	PY-Y05
Otopina za čišćenje, pepsin/kloridna kiselina, 500 ml	PY-Y06

Sartorius Croatia - Libra elektronik d.o.o.

Savska 45A, 10 290 Zaprešić

Ovlašteno zastupstvo i servis Sartorius RH

[Sartorius.Croatia@email.t-com.hr](mailto:Sartorius.Croatia@email.t-com.hr), [service@sartorius.hr](mailto:service@sartorius.hr)

[www.Sartorius.hr](http://www.Sartorius.hr)

Tel: 01 3340-290/291/293/295/296, fax: 01 3340-299

P.J. Solin, Zgon 80, 21210 Solin, 091 563 2266



sartorius