

Contents of the Weight Sets

- Weights without marking
- Weights with marking

Range	Contents	mg	g	kg
1 mg – 5 g Total contents: 11.11 g 16 weights	1	●	●	
	2	● ○	● ○	
	5	●	●	
	10	●		
	20	● ○		
	50	●		
	100	●		
	200	● ○		
500	●			
1 mg – 100 g Total contents: 211.11 g 21 weights	1	●	●	
	2	● ○	● ○	
	5	●	●	
	10	●	●	
	20	● ○	● ○	
	50	●	●	
	100	●	●	
	200	● ○		
500	●			
1 mg – 200 g Total contents: 611.11 g 23 weights	1	●	●	
	2	● ○	● ○	
	5	●	●	
	10	●	●	
	20	● ○	● ○	
	50	●	●	
	100	●	●	
	200	● ○	● ○	
500	●			
1 mg – 1 kg Total contents: 2,111.11 g 25 weights	1	●	●	●
	2	● ○	● ○	
	5	●	●	
	10	●	●	
	20	● ○	● ○	
	50	●	●	
	100	●	●	
	200	● ○	● ○	
500	●	●		
1 mg – 5 kg Total contents: 11,111.11 g 28 weights	1	●	●	●
	2	● ○	● ○	● ○
	5	●	●	●
	10	●	●	
	20	● ○	● ○	
	50	●	●	
	100	●	●	
	200	● ○	● ○	
500	●	●		
1 g – 1 kg Total contents: 2,110 g 13 weights	1		●	●
	2		● ○	
	5		●	
	10		●	
	20		● ○	
	50		●	
	100		●	
	200		● ○	
500		●		
1 g – 5 kg Total contents: 11,110 g 16 weights	1		●	●
	2		● ○	● ○
	5		●	●
	10		●	
	20		● ○	
	50		●	
	100		●	
	200		● ○	
500		●		
1 g – 10 kg Total contents: 21,110 g 17 weights	1		●	●
	2		● ○	● ○
	5		●	●
	10		●	●
	20		● ○	
	50		●	
	100		●	
	200		● ○	
500		●		

*From 1 kg and up with gloves